

 **Bioaktif**[®]
Support For Motion

**20
24**
CATALOG



Biofix[®]

Orthopedic Soft Products



Hello everyone,

Established in 1995, Bioaktif, while operating as an orthopedic prosthesis orthosis centre, started producing orthopedic soft products in 1999 in order to reduce the dependence on imports of ready-made orthopedic products in our country and to produce new products. Meanwhile, considering that the issue of raw materials is the most important factor in producing quality products and product development, he founded Bioteks narrow weaving and textile limited company in 2005 and started the production of orthopedic elastic fabrics.

In 2007, it established its machinery for the production of flat knitted bandages and pioneered the sector.

Bioaktif, which produces orthopedic soft products under the Biofix brand, orthopedic slippers under the Footflex brand, and clothes suitable for use after plastic surgery under the Bodyflex Liposuction Garments brand, has taken its place in Türkiye and the world markets today. It provides services by adhering to innovation, quality and commitments in order to constantly improve and produce products at world standards.

Bioaktif, which has a wide production range, exports its products to 20 countries with its export network in its sector. Bioaktif, which has been developing and growing in Türkiye and the world market since 1995, continues to develop and produce world-class products.

As Bioaktif, our priority is to always exceed our customers' expectations, take new initiatives and find the most appropriate solutions to patient problems.

Yours faithfully,
Mehmet BÜYÜKÇANGA



bioaktif



Contents

Cervical Zone Orthosis

3-6

Shoulder and Elbow Orthosis

7-17

Wrist Orthosis

18-26

Waist and Upper Body Orthosis

27-41

Hip Orthosis

42-43

Knee Orthosis

44-53

Ankle Orthosis

54-60

Foot Health

61-65

Orthopedic Pillows

66-68

Walking Aids

69-70

Biofix[®]
Orthopedic Soft Products



CERVICAL ZONE ORTHOSIS



BA 506

PHILADELPHIA COLLAR

Functions and Features:

- Ensures cervical zone support and stabilization.
- Anatomically appropriate for thorax, mandibula and occiput.
- Stabilizes cervical spines in accordance with the anatomical position.
- Consists of two parts (front and back).
- Has plastic enhancers in front and back.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcros on the sides.
- Reduced perspiration thanks to the air holes.

Indication:

- Cervical trauma
- Disk hernia
- Cervical stabilization
- Support and stabilization of cervical spine in pre and post-operative periods.



Size

| XS/0 | S/1 | M/2 | L/3 | XL/4 |
|-------|-------|-------|-------|-------|
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm |

BA 507

PHILADELPHIA COLLAR (Open Tracheostomy)

Functions and Features:

- Ensures cervical zone support and stabilization.
- Anatomically appropriate for thorax, mandibula and occiput.
- Stabilizes cervical spines in accordance with the anatomical position.
- Consists of two parts (front and back).
- Has plastic enhancers in front and back.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcros on the sides.
- Reduced perspiration thanks to the air holes.

Indication:

- Cervical trauma
- Disk hernia
- Cervical stabilization
- Support and stabilization of cervical spine in pre and post-operative periods.



Size

| XS/0 | S/1 | M/2 | L/3 | XL/4 |
|-------|-------|-------|-------|-------|
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm |



BA 508

NELSON PLASTOZOTE COLLAR

Functions and Features:

- Supports cervical zone.
- Limits the neck in the direction of flexion.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcro

Indication:

- Disk hernia
- Cervical muscle spasm.
- Cervical spine support.
- For cervical disc hernias, injuries, and protection in post-trauma period.



Size

| XS/0 | S/1 | M/2 | L/3 | XL/4 |
|-------|-------|-------|-------|-------|
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm |



BA 509

CERVICAL COLLAR WITH CHIN

Functions and Features:

- Ensures cervical zone support and immobilization.
- Covers the chin to limit flexion.
- Made of plastic.
- Covered with artificial leather.
- Has a foam pad appropriate for the anatomical structure of the chin.

Indication:

- In cervical disc hernia,
- Disk hernia
- Cervical muscle spasm, Cervical spine support
- Injuries, for recovery and protection in post-surgery period.



| Size | | | | | |
|-------|-------|-------|-------|-------|-------|
| XS/0 | S/1 | M/2 | L/3 | XL/4 | XXL/5 |
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm | 50 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm | 54 cm |



BA 510

CERVICAL COLLAR WITHOUT CHIN

Functions and Features:

- Ensures cervical zone support and immobilization.
- Made of plastic, Covered with artificial leather.
- Adjustable neck height.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm,
- Cervical spine support.
- Adjusted according to the size of neck circle.



| Size | | | | | |
|-------|-------|-------|-------|-------|-------|
| XS/0 | S/1 | M/2 | L/3 | XL/4 | XXL/5 |
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm | 50 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm | 54 cm |



BA 511

FOAM COLLAR

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



| Size | | | | |
|-------|-------|-------|-------|-------|
| XS/0 | S/1 | M/2 | L/3 | XL/4 |
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm |



BA 511 U

FOAM CERVICAL COLLAR ORTHOSIS

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



| Size | |
|-------|--|
| 6 cm | |
| 8 cm | |
| 10 cm | |
| 12 cm | |





BA 515

NELSON COLLAR (FABRIC COATED)

Functions and Features:

- Ensures cervical zone support and stabilization. Made of hypoallergenic plastazote material.
- Covered with cotton, non-sweating fabric.

Indication:

- Traumatic cases, Osteoarthritis, Cervical tendinitis.
- Muscle spasm, disk hernia, for protection in post-trauma period.



Size

| XS/0 | S/1 | M/2 | L/3 | XL/4 |
|-------|-------|-------|-------|-------|
| 28 cm | 33 cm | 38 cm | 42 cm | 46 cm |
| 32 cm | 37 cm | 42 cm | 46 cm | 50 cm |

FOAM COLLAR UNI.

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size (Universal)

BA 516 U



BA 511 P



FOAM COLLAR PEDIATRIC

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size (Universal)

BA 513

MINERVA CTO

Functions and Features:

- Ensures cervical zone support and stabilization.
- Cervical supports that support chin and occipital zone are made of plastic, and the product has body supports that support chest and back. A connection was ensured between body and cervical supports using metal pieces.
- The product gets support from the shoulder through chin and shoulder parts, and ensures extraction and fraction to the cervical zone, if required.

Indication:

- Used for protective purposes after Osteochondrosis, Spondylose, Spondylarthrose, Tortikolis and other operations.



Size

| S | M | L |
|-------|-------|-------|
| 34 cm | 38 cm | 42 cm |
| 38 cm | 42 cm | 46 cm |





SHOULDER & ELBOW ORTHOSIS

Biofix[®]
Orthopedic Soft Products



BA 413 U

EPIROM COMFORT ROM ELBOW ORTHOSIS UNI.

Functions and Features:

- An elbow orthosis with adjustable angle, designed in accordance with the anatomical elbow shape.
- Inner surface that touches the skin is made of foam-backed fabric that reduces sweating.
- The mechanism used for angle adjustments is an easy-to-use pushing system.
- Orthosis controlled elbow ensures flexion and extraction.
- Flexion (0-15-30-45-60-75-90-105) and extraction (-15-0-15-30-45-60-75-90)
- Angle can be adjusted. Orthosis length can be adjusted according to the length of arm.
- Aluminum support that avoids the wrist to bend is available and length can be adjusted.
- Shoulder belt offers comfort in use.

Indication:

- Support for rehabilitation in elbow contracture.
- Protection and support for rehabilitation in elbow fracture.
- Limited and controlled movement in the elbows.
- Conservative treatment after elbow dislocation, in post-op period for ligament instability.
- Post-op rehabilitation after elbow endoprosthesis.



Size (Universal Right and Left)

BA 412

HUMERUS SERMIANTO ORTHOSIS UNI.

Functions and Features:

- Made of 2 mm polyethylene, and covered with foam fabric on the inside.
- This fabric can be removed and washed.
- Ensures full harmony by surrounding the humerus with integrated band system.

Indication:

- With the aim of the protection of Humerus proximal fractures after operation, plaster treatment and stabilization.



Size (Universal)

| S/1 | M/2 |
|-------|-------|
| 19 cm | 26 cm |
| 25 cm | 30 cm |



BA 20505-1

SHOULDERFIX PLUS

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Elastic tapes on the shoulder.
- Strengthens the muscles by compressing the tissues with its elastic structure, avoids swelling.
- Limits rotation and spreading movements of the shoulders.
- Avoids humerus head to distance from shoulder joint, and limits the injuring movements.

Indication:

- Shoulder bone fractions • Shoulder arthritics.
- Rehabilitation after shoulder implant prosthesis treatment
- Post-operative and post-traumatic irritations.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 22 cm | 24 cm | 27 cm | 30 cm | 33 cm |
| 24 cm | 27 cm | 30 cm | 33 cm | 37 cm |



BA-20505-2



BA 110

DELUX ARM SLING (With Waist Belt)

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.



| Size | XS | S | M | L | XL |
|------|-------|-------|-------|-------|-------|
| | 25 cm | 30 cm | 35 cm | 40 cm | 45 cm |
| | 30 cm | 35 cm | 40 cm | 45 cm | 50 cm |



BD 110 U

ARM SLING (With Waist Belt) UNI.

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.



Size (Universal)

BA 110 P

DELUX ARM SLING (With Waist Belt) PEDIATRIC

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.



Size (Universal)





BA 110 A

DELUX ARM SLING (With Waist Belt)

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.

| Size | XS | S | M | L | XL |
|------|-------|-------|-------|-------|-------|
| | 25 cm | 30 cm | 35 cm | 40 cm | 45 cm |
| | 30 cm | 35 cm | 40 cm | 45 cm | 50 cm |

LUX ARM SLING (With Waist Belt)

Function and Features:

- Made of air-permeable airtex spacer fabric.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Easy to wear and remove.
- Shoulder and waist belts do not stretch, and do not cause discomfort for the patient.

Indication:

- Used when immobilization of arm and shoulder is required.

| Size | XS | S | M | L | XL |
|------|-------|-------|-------|-------|-------|
| | 25 cm | 30 cm | 35 cm | 40 cm | 45 cm |
| | 30 cm | 35 cm | 40 cm | 45 cm | 50 cm |

BA 115



BD 115 U



ARM SLING ADJUSTABLE

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.

| | |
|--|------------------|
| | Size (Universal) |
|--|------------------|

BA 116


ECO ARM SLING (With Waist Belt)

Function and Features:

- Made of air-permeable airtec spacer fabric.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Designed to be used in both arms.

Indication:

- Used when immobilization of arm and shoulder is required.

| Size | XS | S | M | L | XL |
|---|----------------|----------------|----------------|----------------|----------------|
|  | 25 cm 30 cm | 30 cm 35 cm | 35 cm 40 cm | 40 cm 45 cm | 45 cm 50 cm |



BA 140


ACTIVE SLING 30° (Arm Sling with Pad)

Function and Features:

- Used in immobilization and positioning of shoulder.
- Made of 6 mm, spongy, air-permeable, foam-backed fabric.
- The pad that abducts the shoulder at 30° is made of 25x12 foam.
- Designed to ensure comfort in long-term treatment of arm, and to offer armpit air circulation.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

Indication:

- Rotator cuff, capsular rupture, and bankart repair. Used in shoulder joint dislocation and soft tissue injury.

| Size | S | M | L |
|---|----------------|----------------|----------------|
|  | 30 cm 35 cm | 35 cm 40 cm | 40 cm 45 cm |



BA 140 A


ACTIVE SLING 30° (Arm Sling with Pad)

Function and Features:

- Used in immobilization and positioning of shoulder.
- Made of 6 mm, spongy, air-permeable, foam-backed fabric.
- The pad that abducts the shoulder at 30° is made of 25x12 foam.
- Designed to ensure comfort in long-term treatment of arm, and to offer armpit air circulation.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

Indication:

- Rotator cuff, capsular rupture, and bankart repair. Used in shoulder joint dislocation and soft tissue injury.

| Size | S | M | L |
|---|----------------|----------------|----------------|
|  | 30 cm 35 cm | 35 cm 40 cm | 40 cm 45 cm |





BA 150

ACTIVE SLING 45°-60° (Arm Sling with Pad)

Functions and Features:

- Used in immobilization and positioning of shoulder.
- Six mm inner surface is made of foam-backed fabric that reduces sweating.
- Ensured 45° or 60° abduction angle on shoulder.
- Designed to ensure comfort in long-term treatment of arm.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

Indication:

- In shoulder prosthesis operations, rotator cuff, capsular rupture, and bankart repairs. Used in shoulder joint dislocation and soft tissue injury.

| Size | S | M | L |
|------|-------|-------|-------|
| | 30 cm | 35 cm | 40 cm |
| | 35 cm | 40 cm | 45 cm |

VELPOFIX (Velpeau Bandage)

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.

| | Size (Right-Left) | | | | |
|--|-------------------|-------|-------|--------|--------|
| | XS | S | M | L | XL |
| | 60 cm | 70 cm | 80 cm | 90 cm | 100 cm |
| | 70 cm | 80 cm | 90 cm | 100 cm | 110 cm |

BA 160



BD 160 U

VELPOFIX (Velpeau Bandage) UNI.

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.

| | Size (Universal) |
|--|------------------|
| | |

BD 160 E

VELPOFIX (Velpeau Bandage) ECO

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.



Size (Universal)



BA 160 P

VELPOFIX (Velpeau Bandage) PEDIATRIC

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.



Size (Universal)



BA 170

CLAVICULAFIX (Clavicle Bandage)

Function and Features:

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

| XS | S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|--------|
| 60 cm | 70 cm | 80 cm | 90 cm | 100 cm | 110 cm |
| 70 cm | 80 cm | 90 cm | 100 cm | 110 cm | 120 cm |





BD 170 U

CLAVICULAFIX (Clavicle Bandage) UNI.

Function and Features:

- Made of 8 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size (Universal)

CLAVICULAFIX (CLAVICLE BANDAGE) PEDIATRIC

Function and Features:

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

| XS | XXS | XXXS |
|-------|-------|---------|
| 50 cm | 40 cm | < 40 cm |
| 70 cm | 50 cm | |

BA 170 P



BA 171

CLAVICULAFIX PLUS (Clavicle Bandage)

Functions and Features:

- Designed in accordance with anatomic structure.
- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

| XS | S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|--------|
| 60 cm | 70 cm | 80 cm | 90 cm | 100 cm | 110 cm |
| 70 cm | 80 cm | 90 cm | 100 cm | 110 cm | 120 cm |

EPIFIX ACTIVE (Epicondylitis Elbow Band)

BA 185

Functions and Features:

- Produced as a standard size elbow band from velcro-compatible fabric, and offers epicondylitis support.
- Epicondylitis support is anatomically ensured with a silicon pad.
- The support reduces the stress on tendon, and ensures relief of pain.

Indication:

- Medial and lateral epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



BA 186-1

EPITEC COMFORT (Epicondylitis Elbow Band)

Functions and Features:

- Produced as a standart size elbow band from velcro-compatible fabric, and offers epicondylitis support.
- Epicondylitis support is anatomically ensured with a silicon pad.
- The support reduces the stress on tendon, and ensures relief of pain.

Indication:

- Medial and lateral epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



BA 186-2

EPICONDYLITIS ELBOW BAND UNI.

BD 180 U

Function and Features:

- Produced as a standard size elbow band from soft fabric, and offers epicondylitis support.
- Epicondylitis support is ensured with a pad.
- This support
- Reduces the stress.

Indication:

- Medial and lateral epicondylitis, epicondylitis humeri (tennis elbow)



Size (Universal)



BA 20301-1



EPIFIX PLUS (Silicon Backed Elbow Support With Extra Band)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that apply pressure on condyls relieve the pain.

Indication:

- Lateral and Medial Epicondylitis (Tennis Elbow) Elbow dislocation. Chronic post-operative and post-traumatic inflammation.



| Size | | | | | |
|-------|-------|-------|-------|-------|-------|
| S | M | L | XL | XXL | |
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm | |



BA 20301-2

BA 20302-3

EPIFIX BASIC (Elastic Knitted Elbow Support)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Appropriate for use in sport.

Indication:

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



| Size | | | | | |
|-------|-------|-------|-------|-------|-------|
| S | M | L | XL | XXL | |
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm | |



BA 20302-4



EPIFIX BASIC (Elastic Knitted Elbow Support)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Appropriate for use in sport.

Indication:

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



| Size | | | | | |
|-------|-------|-------|-------|-------|-------|
| S | M | L | XL | XXL | |
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm | |



BA 20302-5

EPIFIX BASIC (Elastic Knitted Elbow Support)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Appropriate for use in sport.

Indication:

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |

BA 905

EPITEX LUX ELBOW SLEEVE

Functions and Features:

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

Indications:

- Elbow Dislocation
- Post-operative or Posttraumatic inflammation
- Osteoarthritis, Arthritis
- Tenosynovitis



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |



BD 30603 U

EPICONDYLITIS ELBOW BANDAGE WITH PAD UNI.

Functions and Features:

- The bandage that ensures epicondylitis support is made of air-permeable, three dimensional airtec spacer fabric.
- Reduces edema with compression.
- Reduces sweating, and avoids allergy.
- Epicondylitis support is ensured with an extra band and pad. This support reduces the stress on tendon, and ensures relief of pain.

Indication:

- Medial and lateral Epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



WRIST ORTHOSIS

Biofix[®]
Orthopedic Soft Products



BA 20601-1

MANU COMFORT (Knitted Wrist Splint)

Functions and Features:

- Produced with an advanced technology using smart knitting-waving comfort polyamide threads.
- Air-permeable, breathable. Does not hold moisture, transfers the moisture to outside.
- Hygienic, non-allergic, anti-bacterial.
- Shaped in accordance with the anatomical structure, fits the hand and wrist.
- Offers high comfort and harmony. A formable thermoplastic splint is available.
- The silicon pad on the volar side of the hand ensures support, and relieves the pain.

Indication:

- Tenosynovitis
- Carpal instability
- CTS
- After plaster treatment of wrist fractures
- Rheumatoid arthritis
- Ganglion pain
- Soft tissue traumas on the wrist



Size (Right-Left)

| S | M | L | XL | XXL |
|----------------|----------------|----------------|----------------|----------------|
| 14 cm 16 cm | 16 cm 18 cm | 18 cm 20 cm | 20 cm 22 cm | 22 cm 24 cm |



BA 20601-2

WRIST SLEEVE LUX

Functions and Features:

- It is produced with circular knitting technique providing balanced compression. The inner surface is made of cotton outer surface polyamide thread.
- Air permeable and moist regulator.
- High wearing comfort.

Indication:

- Wrist soft tissue injuries
- Edema and pain around the wrist
- Wrist support



Size

| S | M | L | XL | XXL |
|----------------|----------------|----------------|----------------|----------------|
| 14 cm 16 cm | 16 cm 18 cm | 18 cm 20 cm | 20 cm 22 cm | 22 cm 24 cm |

BA 907



BA 460-1



MANU LACE (Wrist Splint With Lace)

Functions and Features:

- Supports sensorimotor movements, and ensures muscle stability of the wrist.
- Made of air-permeable, good-for-skin, elastic textile fabric.
- Ensures comfortable stabilization thanks to the adjustable laces that wrap the wrist.
- Limits the hand-wrist movement. Fully adjusts with the anatomical structure of the wrist.
- Keeps your wrist at neutral position, relieves you pain, helps through rehabilitation process

Indication:

- Acute and chronic irritation (e.g. tendovaginitis, tendomyopathy, minor joint effusion, arthralgia).
- Post-operative and post-trauma irritation.
- Wrist Injuries.
- Rheumatoid Arthritis.
- Carpal Tunnel Syndrome
- Degenerative Wrist Disorders (e.g. arthritis)
- Sense of weakness.



Size (Right-Left)

| S/1 | L/2 |
|----------------|----------------|
| 14 cm 18 cm | 18 cm 22 cm |

BD 460 U

MANU LACE (Wrist Splint With Lace) UNI.

Functions and Features:

- Supports sensorimotor movements, and ensures muscle stability of the wrist.
- Made of air-permeable, good-for-skin, elastic textile fabric.
- Ensures comfortable stabilization thanks to the adjustable laces that wrap the wrist.
- Limits the hand-wrist movement. Fully adjusts with the anatomical structure of the wrist.
- Keeps your wrist at neutral position, relieves you pain, helps through rehabilitation process

Indication:

- Acute and chronic irritation (e.g. tendovaginitis, tendomyopathy, minor joint effusion, arthralgia).
- Post-operative and post-trauma irritation.
- Wrist Injuries.
- Rheumatoid Arthritis.
- Carpal Tunnel Syndrome
- Degenerative Wrist Disorders (e.g. arthritis)
- Sense of weakness.



Size (Universal)

| S/1 | L/2 |
|-------|-------|
| 14 cm | 18 cm |
| 18 cm | 22 cm |



BA 440

MANUFIX (Elastic Wrist Splint)

Function and Features:

- Ensures wrist support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Has anatomically shaped aluminum support on the volar surface of the wrist.
- Easy to wear and remove thanks to the velcros.

Indication:

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 14 cm | 16 cm | 18 cm | 20 cm | 22 cm |
| 16 cm | 18 cm | 20 cm | 22 cm | 24 cm |



BA 30601

MANUAI RTEC (Wrist Splint)

Functions and Features:

- Ensures wrist support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Aluminum support that keeps the wrist at neutral position is available. Ensures wrist support and stabilization.
- Stability is increased with the extra tape on the wrist.

Indication:

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 14 cm | 16 cm | 18 cm | 20 cm | 22 cm |
| 16 cm | 18 cm | 20 cm | 22 cm | 24 cm |



BA 30601 P



MANUAI RTEC (Wrist Splint) PEDIATRIC

Functions and Features:

- Ensures wrist support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Aluminum support that keeps the wrist at neutral position is available. Ensures wrist support and stabilization.
- Stability is increased with the extra tape on the wrist.

Indication:

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

| | | |
|-------|-------|--------|
| XS | XXS | XXXS |
| 13 cm | 11 cm | <11 cm |
| 15 cm | 13 cm | |

BD 40601 U

UNIFIX MANU (Universal Wrist Splint) UNI.

Function and Features:

- Ensures wrist support and stabilization.
- Made of 3 mm airtec spacer fabric. Reduces the pain in the arthritis thanks to the heating feature of neoprene.
- An adjustable aluminum support on the volar surface of the wrist, and an extra bandage around the wrist are available.

Indication:

- For immobilization in wrist pains, for carpal tunnel syndrome, tendinitis, tenosynovitis, tenovaginitis, injuries, for protection in post-operative period, and after plaster treatment in fractures.



Size (Universal Right-Left)



BA 932



MANUFIX LUX WRIST SUPPORT

Function and Features:

- Air permeable and moist regulator.
- It covers the wrist, stabilizes and provides thumb movement.
- There is immobilization with ligament supports on the inner and outer sides of the palm.

Indication:

- Carpal Tunnel Syndrome
- Painful wrist conditions
- Wrist strain



Size (Universal)

| | |
|-------|-------|
| S/1 | L/2 |
| 14 cm | 18 cm |
| 18 cm | 22 cm |

BA 450

MANUFIX PLUS (Wrist Splint With Thumb)

Function and Features:

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

Indication:

- MCP joint disorders, CMC joint degenerations, De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 14 cm | 16 cm | 18 cm | 20 cm | 22 cm |
| 16 cm | 18 cm | 20 cm | 22 cm | 24 cm |



BA 30602

MANUAI RTEC PLUS (Wrist Splint With Thumb)

Function and Features:

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Ensures wrist and thumb support and stabilization.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

Indication:

- MCP joint disorders, CMC joint degenerations
- De Quervain, and other cases where carpal tunnel and De Quervain are observed together.

Size (Right-Left)

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 14 cm | 16 cm | 18 cm | 20 cm | 22 cm |
| 16 cm | 18 cm | 20 cm | 22 cm | 24 cm |



BA 30602 P

MANUAI RTEC PLUS (Wrist Splint With Thumb)

Function and Features:

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Ensures wrist and thumb support and stabilization.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

Indication:

- MCP joint disorders, CMC joint degenerations
- De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

| XS | XXS | XXXS |
|-------|-------|--------|
| 13 cm | 11 cm | <11 cm |
| 15 cm | 13 cm | |





BD 40602 U

UNIFIXMANU PLUS (Wrist Splint With Thumb Support) UNI.

Functions and Features:

- Ensures wrist and thumb support and stabilization.
 - Made of neoprene fabric - Reduces the pain in the arthritis thanks to the heating feature of neoprene.
 - An adjustable aluminum baleen is available for volar support.
 - Easy to apply on the patient thanks to its velcro tapes.
- Aluminum baleen is used in stabilization at neutral position. Produced as a single body.

Indication:

- Used in injuries, twists, rheumatic and traumatic arthritis, tendinitis, carpal tunnel syndrome pains, and for support in all plaster treatment cases requiring wrist rest.



Size (Universal Right-Left)

DEQUFIX LUX THUMB SUPPORT

Function and Features:

- Air permeable and moist regulator.
- Designed to cover wrist and thumb.
- Support and stabilization of thumb joints

Indication:

- Tendinitis
- Skier finger
- De Quervain Syndrome



Size (Universal)

| S/1 | L/2 |
|-------|-------|
| 14 cm | 18 cm |
| 18 cm | 22 cm |

BA 931



BA 451 U



DEQUFIX (Thumb Fixing Splint) UNI.

Function and Features:

- Ensures support and stabilization of thumb joints.
- Made of elastic textile fabric.
- Has a formable aluminum support.
- Fixes CMC and MC joints.

Indication:

- In treatment of chronic, post-operative, and post-traumatic tissue sensitivity.
- For stabilization in ulnar cruciate ligaments.
- Carpal Tunnel Syndrome.
- Rheumatoid arthritis.
- Root arthrosis.



Size (Universal Right-Left)

THUMB SPLINT

Functions and Features:

- Limits IBP / CMC joint movements.
- Made of soft and breathable fabric.
- Easy-to-use and adjustable thanks to elastic tapes.
- De Quervain tendinitis
- Gamekeeper thumb
- Collateral ligament immobilization



Size (Universal Right-Left)

BA 453



BD 40600 U

UNIFIX

Function and Features:

- Produced as a single body from neoprene material.
- Keeps the thumb at resting position.

Indications:

- Thumb collateral ligament injuries.
- Root arthrosis.
- Rheumatoid arthritis.



Size (Universal)



BA 455

BOXER FRACTURE SPLINT

Functions and Features:

- The Biofix Boxer Splint is an easy to adjust splint that positions and protects the 4th and 5th metacarpals in the treatment of "boxer fractures", and offers superior comfort and support to enhance compliance.
- The malleable palmar and dorsal aluminum stays are easily bent to immobilize the patient's hand.
- The products design protects and supports the fourth and fifth metacarpals without restricting flexion.

Indications:

- Boxer fracture



Size (Universal Right-Left)



BA 183

PATIENTS HANDS AND FEET STRAP BELT

Functions and Features:

- Made of soft fabric with padding, breathable and comfortable to wear
- Size: adjustable, suitable for wrist or ankle circumference: 8.6" ~ 11.4" / 22 ~ 29cm
- Velcro close, easy and convenient to use
- Use this medical restraint strap to prevent patients from self-injury to ensure safety and smooth care or treatment
- Unconscious bedridden patients, mental disorder people, cognitive impairment people etc. applicable to.

Size (Universal)



BA 430

MANUSTABIL SPASTIC WRIST SPLINT

Function and Features:

- Ensures neutral position in hands and wrists.
- Made of 3 mm polyethylene using vacuum method.
- Internal surface is covered with 3 mm plastozote.
- Durable, can be used in long-term treatments, easy-to-wear and remove.

Indication:

- For protections after fractures. For paralyzed patients, in wrist problems, and in cases where hand is required to be rested.



Size (Right-Left)

| S | M | L |
|-------|-------|-------|
| 14 cm | 16 cm | 18 cm |
| 16 cm | 18 cm | 20 cm |





BA 431

**MANU ANTI-SPASTISITE
(Anti-Spasticity Wrist Splint Separator)**

Functions and Features:

- Keeps the hand, wrist, and fingers at neutral position.
- Made of 3 mm polyethylene using vacuum method.
- Internal surface is covered with 3 mm plastazote.
- Durable, can be used in long-term treatments, easy-to-wear and remove.
- Ideal for patients with spasticity due to its channeled finger design and resistance to stress.
- Ensures full harmonization and control with the special-designed tape system.

Indication:

- For patients with spasticity, protective after fractures. For paralyzed patients, in wrist problems, and in cases where hand is required to be rested.



Size (Right-Left)

| S | M | L |
|----------------|----------------|----------------|
| 14 cm 16 cm | 16 cm 18 cm | 18 cm 20 cm |

FROG FINGER SPLINT

Function and Features:

- Used for fixing the fingers at required position.
- Internal surface is covered with soft plastazote.
- Can be formed if required.



Size

| S | M | L |
|---|---|---|
|---|---|---|

BA 101



BA 102



BASEBALL FINGER SPLINT

Function and Features:

- Used to keep the finger at extension position.
- Used in finger injuries, tendon splitting cases, and low joint cases.
- Size and numbers varies.



Size

| S | M | L |
|---|---|---|
|---|---|---|

BA 104

MALLET FINGER SPLINT

Function and Features:

- Mallet finger splint is a finger orthosis made of special material and is used to keep the DIP joint at extraction position.
- Finger is free from the point of PIP joint.
- Extensor tendon injuries.
- Finger Splint is used in fractures.
- Used in injuries and twists.

Indication:

- Extensor tendon treatment, mallet finger deformity. Finger bone fractures. Finger injuries. Finger nerve damages.



Size

| | | | | | | |
|---|---|---|---|---|-----|---|
| 1 | 2 | 3 | 4 | 5 | 5.5 | 6 |
|---|---|---|---|---|-----|---|



BA 118

DYNAMIC FINGER SPLINT

Function and Features:

- Made of stainless steel spring covered with eloxal.
- Points that touch the finger are made of special plastic with rounded edges.
- Produced in two sizes.
- The center point of the spring is attached using a screw, while the other edges are assembled with stainless rivet.
- Used in tendon cuts and nerve injuries.

Indication:

- Improvement of muscle immobility,
- Improvement of joint immobility,
- Ulnar nerve damages,
- Rehabilitation in post-operative period



Size

| | | |
|---|---|---|
| S | M | L |
|---|---|---|





W A I S T & U P P E R B O D Y O R T H O S I S

Biofix[®]
Orthopedic Soft Products



BA 20501-2

LUMBOSOFT (Knitted Lumbosacral Corset)

Functions and Features:

- Supports the lumbar region.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Compatible with the anatomy, and comfortable. Ensures anatomically appropriate support.
- Easy to wear and remove thanks to the velcros on the front side.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome.



| Size | | S | M | L | XL | XXL |
|-------|-------|-------|-------|--------|--------|--------|
| 75 cm | 83 cm | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |
| 83 cm | 91 cm | 91 cm | 99 cm | 107 cm | 115 cm | |



BA 20501-1

BA 240

LUMBOLIGHT PLUS (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



| Size | | S | M | L | XL | XXL |
|-------|-------|-------|-------|--------|--------|--------|
| 75 cm | 83 cm | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |
| 83 cm | 91 cm | 91 cm | 99 cm | 107 cm | 115 cm | |



BA 245

LUMBOLIGHT EVO

Functions and Features:

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbar.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 26 cm, Front height: 16 cm.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome,



| Size | | S | M | L | XL | XXL |
|-------|-------|-------|-------|--------|--------|--------|
| 75 cm | 83 cm | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |
| 83 cm | 91 cm | 91 cm | 99 cm | 107 cm | 115 cm | |



BD 240 U

LUMBOLIGHT UNI. (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



| Size | |
|---------|----------|
| S-M / 1 | L-XL / 2 |
| 75 cm | 99 cm |
| 99 cm | 115 cm |



BD 240 US

LUMBOLIGHT UNI. (Lumbosacral Corset With An Extra Belt)

Function and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic



| Size | |
|---------|----------|
| S-M / 1 | L-XL / 2 |
| 75 cm | 99 cm |
| 99 cm | 115 cm |



BA 20503-1

LUMBOSACRAL CORSET

Function and Features:

- Supports the lumbar region.
- A triangle silicon pad-support appropriate for the anatomy of the body is placed on the back side of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Compatible with the anatomy, and comfortable. Silicon support appropriate for the anatomical structure is available.
- Easy to wear and remove thanks to the velcros on the front side.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-dissectomy syndrome.



| Size | | | | | |
|-------|-------|-------|--------|--------|--|
| S | M | L | XL | XXL | |
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm | |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm | |



BA 20503-2



BD 243 U

WAIST SUPPORT BLACK (17 cm)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, nylon fabric that avoids sweating.
- Four plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove.

Indication:

- Acute waist pain, degenerative disc, paravertebral muscle spasm.



| Size | |
|---------|----------|
| S-M / 1 | L-XL / 2 |
| 75 cm | 99 cm |
| 99 cm | 115 cm |

**LUMBOLIGHT UNI ECO.
(Lumbosacral Corset With An Extra Belt)**

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, polyester fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove.

Indication:

- Acute waist pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



| Size | | | | | |
|-------|-------|-------|--------|--------|--|
| S | M | L | XL | XXL | |
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm | |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm | |

BD 240 E



BA 20502-2

LUMBOSOFT (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Easy-to-wear thanks to hand bands on the front right and left sides of the corset.
- Back height: 30 cm, Front height: 16 cm.
- Offers more rigid support and stabilization with the extra bands.

Indication:

- Lumbar spine instability, lumbar stress syndromes, post-operative treatment after disc surgery.



| Size | | | | | |
|-------|-------|-------|--------|--------|--|
| S | M | L | XL | XXL | |
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm | |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm | |



BA 20502-1

BA 246

LUMBOLIGHT EVO PLUS

Functions and Features:

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbar.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 26 cm, Front height: 16 cm.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, disorders, post-discectomy syndrome, ligamentous



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 241

LUMBOLIGHT PLUS (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 242

LUMBOLIGHT CORSET BLACK (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Four plastic baleens are available at the back side.
- Easy to wear and remove.

Indication:

- Acute waist pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |





BA 230

LUMBOLIGHT (Lumbosacral Corset 26 cm)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two adjustable steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset offers full support.

Indication:

- Can be used in lumbago, sciatic, and herniated disk problems, for post-operative support, and in physiotherapy as a stable support.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |

LUMBOFIX PLUS (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support, and ensure full harmonization with the body.

Indication:

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago - sciatic, spondylolisthesis, post-operative immobilization.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |

BA 260



BD 260 U

LUMBOFIX UNI. (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Produced in accordance with the anatomy of the body.
- Back height: 30 cm, Front height: 18 cm.

Indication:

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago sciatic, spondylolisthesis, post-operative immobilization.



Size

| S-M / 1 | L-XL / 2 |
|---------|----------|
| 75 cm | 99 cm |
| 99 cm | 115 cm |



BD-260 US

LUMBOFIX UNI. (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Produced in accordance with the anatomy of the body.
- Back height: 30 cm, Front height: 18 cm.

Indication:

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago sciatic, spondylolisthesis, post-operative immobilization.



Size

| S-M / 1 | L-XL / 2 |
|---------|----------|
| 75 cm | 99 cm |
| 99 cm | 115 cm |



BA 265

LUMBOFIX EVO

Functions and Features:

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbosacral.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 32 cm, Front height: 20 cm.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome, post-operative support and stabilization.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 933

STRONG LUMBAR CORSET

Functions and Features:

- Ensures stabilization and support of the lower back.

Indications:

- During the handling of loads, labor protection in continuous overexertion or lumbar



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |





BA 280

DORSOFIX PLUS (Dorsolumbar Corset)

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.

Indication:

- Can be used in compression fractures and herniated disks, for post-operative support, and in physiotherapy as a stable support.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |

DORSOFIX UNI. (Dorsolumbar Corset)

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.

Indication:

- Can be used in compression fractures and herniated disks, for post-operative support, and in physiotherapy as a stable support



Size

| S-M / 1 | L-XL / 2 |
|---------|----------|
| 75 cm | 99 cm |
| 83 cm | 115 cm |

BD 280 U



BA 270

DORSOFIX (Dorsolumbar Corset)

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support.

Indication:

- Compression fractures and herniated disks, post-operative support, and in physiotherapy, as a stable support.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |

BA 287

FABRIC TILSO

Functions and Features:

- Made of rigid and air-permeable fabric.
- Limits flexion-extraction of the body.
- Avoids lateral bending and rotation.
- Polyethylene rigid supports are available at the back and front sides.
- Limits extraction through thorax pad.

Indication:

- To mobilize the patient after surgery,
- In compression fractures, not requiring surgery.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 286



JEWEET / HYPEREXTENSION CORSET

Functions and Features:

- Limits the flexion of the body based on 3 point principle, and ensures hyperextension.
- Stabilizes the spine against lateral flexion and rotation movements.
- A mobile hip bar that ensure greater comfort while seated.
- Four different sizes at the side bars, which can be adjusted.
- Easy use with the B-fast lock.
- Adjustable chest bars.

Indication:

- Stable compression fractures in lumbar and inferior thoracic vertebra.
- Immobilization after intervertebral disc surgery.
- Axial stabilization in thoracic and osteoporotic compression fractures of lumbar spine.
- Temporary orthotic care after open surgeries on non-stable vertebral body fracture.
- Permanent orthotic care for spinal tumors and metastases which are not appropriate for surgery.
- Post-operative orthotic care after vertebral body tumor and metastases surgeries.
- Care after spinal decompression surgery with or without internal fixation.



Size

| S | M | L |
|-------|--------|--------|
| 75 cm | 90 cm | 105 cm |
| 90 cm | 105 cm | 120 cm |

BA 310

POSTURFIX

Functions and Features:

- Ensures support to fix the posture.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Two plastic baleens are available at the back side.
- Bands on the shoulder ensure postural support, and avoid kyphoses.

Indication:

- Used for weak and lazy thoracic zone muscles, and ad stimulant and straightener in posture disorders.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |





BA 225

PREGNANT CORSET

Functions and Features:

- Helps to avoid waist pain during pregnancy.
- Designed to support the abdominal area, and to lift it upwards.
- Ensures maximum comforts thanks to the ergonomic design of the corset.
- Support can be increased or decreased with the elastic bands on the sides.

Indication:

- Lumbar and abdominal support during pregnancy.



| Size | | S | M | L | XL | XXL |
|------|--|-------|-------|-------|--------|--------|
| | | 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| | | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |

ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S | M | L | XL | XXL |
|------|--|-------|-------|-------|--------|--------|
| | | 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| | | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |

BA 300



BD 300 U



ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S-M / 1 | L-XL / 2 |
|------|--|---------|----------|
| | | 75 cm | 99 cm |
| | | 99 cm | 115 cm |

BD 300 U S

ABDOMINAFIX BLACK UNI.

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S-M / 1 | L-XL / 2 |
|------|--|---------|----------|
| | | 75 cm | 99 cm |
| | | 99 cm | 115 cm |



BA 300 E

ABDOMINAFIX ECO

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic polyester fabric that avoids sweating.

Indications:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S | M | L | XL | XXL |
|-------|-------|-------|-------|--------|--------|--------|
| 75 cm | 83 cm | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |
| 83 cm | 91 cm | 91 cm | 99 cm | 107 cm | 115 cm | |



BA 3002

ABDOMINAFIX

Functions and Features:


- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S-M / 1 | L-XL / 2 |
|-------|-------|---------|----------|
| 75 cm | 99 cm | 99 cm | 115 cm |
| 99 cm | | | |



BA 3002-T

BA 3002-S



BA 3003

ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S-M / 1 | L-XL / 2 |
|-------|-------|---------|----------|
| 75 cm | 99 cm | 99 cm | 115 cm |
| 99 cm | | | |



BA 3003-T

BA 3003-S



BA 3004

ABDOMINAFIX

Functions and Features:


- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



| Size | | S-M / 1 | L-XL / 2 |
|-------|-------|---------|----------|
| 75 cm | 99 cm | 99 cm | 115 cm |
| 99 cm | | | |



BA 3004-S

BA 3004-T





BA 302


STERNAFIX (Sternum Corset)

Functions and Features:

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

Indication:

- To support sternum zone in post-operative period.



| Size | | | | |
|-------|-------|--------|--------|--------|
| S | M | L | XL | XXL |
| 75 cm | 85 cm | 95 cm | 105 cm | 115 cm |
| 85 cm | 95 cm | 105 cm | 115 cm | 125 cm |



BD 302 U


STERNAFIX (Sternum Corset)

Functions and Features:

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

Indication:

- To support sternum zone in post-operative period.



| Size | |
|---------|----------|
| S-M / 1 | L-XL / 2 |
| 75 cm | 99 cm |
| 99 cm | 115 cm |



BD 302 ES


STERNAFIX ECO (Sternum Corset)

Functions and Features:

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.

Indication:

- To support sternum zone in post-operative period.



| Size | |
|---------|----------|
| S-M / 1 | L-XL / 2 |
| 75 cm | 99 cm |
| 99 cm | 115 cm |

BD 302 ET



BA 306


STERNUM CORSET WITH STRAP

Functions and Features:

- Ensures immobilization after bypass operation.
- Made of highly air-permeable elastic fabric. Avoids allergy, transfers the moisture to outside.
- Can be used to increase the stabilization by tightening the corset with hand, provided that the patient is to cough after operation.
- Fixed with the shoulder bands to avoid sliding.

Indication:

- After bypass operation.



Size

| S | M | L | XL | XXL |
|-------|-------|--------|--------|--------|
| 75 cm | 85 cm | 95 cm | 105 cm | 115 cm |
| 85 cm | 95 cm | 105 cm | 115 cm | 125 cm |



BA 303

WOOL CORSET

Functions and Features:

- Ensures abdominal and waist support.
- Made of wool, elastic fabric.
- Closed on the front part.
- Keeps warm due to heating feature.

Indication:

- Used to keep the affected areas hot in rheumatic and sciatic pains.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 303 V

WOOL ELASTIC CORSET WITH VELCRO BAND

Functions and Features:

- Made of wool, elastic fabric.
- Used to keep waist and abdominal area, and kidneys warm, in muscle spasms, and for rheumatic and sciatic pains.
- Can be opened and closed with the velcro on the front side.

Indication:

- Used to keep the affected areas hot in rheumatic and sciatic pains.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 305

COLOSTOMY CORSET

Functions and Features:

- Manufactured of flexible and cotton fabric.
- 26 cm long corset.
- With the aid of the flexible fabric and 2 removeable plastic underwires at each sides of the stoma opening.
- Prevents any unwanted movements of the stoma bag and potential leakages which may occur during mundane and sporting activities.
- Offers perception of comfort and assurance to the patient. Suitable for wear underneath daily clothing.

Indication

- Protecting the ostomy collection bag
- Post-ostomy rehabilitation and support purposes
- Decrease risks of parastomal hernia



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |





BA 290

UMBILICAL HERNIA CORSET

Function and Features:

- A support that compresses the hernia area is available.
- Made of air-permeable, elastic, cotton fabric.


Indication:

- Used in umbilical hernias.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |



BA 500

SINGLE TRUSS

Functions and Features:

- Supports the hernia in inguinal area.
- The pad that supports the hernia is made of rubber. Covered with fabric.

Indication:

- Used for hernia, and in post-operative period.



Size Universal -Right / Left



BA 501

DOUBLE TRUSS

Functions and Features:


- Supports the hernia in inguinal area.
- The pad that supports the hernia is made of rubber. Covered with fabric.

Indication:

- Used for hernia, and in post-operative period.



Size Universal



BA 504

SUSPENSORY

Functions and Features:

- Avoids vesicas to lie over, and supports the ones that lied over.
- Made of cotton fabric. Has waist belt, and elastic bands that avoid sliding.

Indication:

- Used to support testicles in sports activities.
- Used in vesica varicoses and inflammation
- After varicocele operation.



Size

| 5 | 6 | 7 | 8 | 9 | 10 |
|-------|-------|-------|--------|--------|--------|
| 75 cm | 83 cm | 91 cm | 99 cm | 107 cm | 115 cm |
| 83 cm | 91 cm | 99 cm | 107 cm | 115 cm | 123 cm |

BA 419 P

FREJKA PILLOW

Functions and Features:

- Made of 12 mm plastozote material.
- Allows compulsory hip and leg movements of the baby, as it is made of elastic material.

Indication:

- Coksit, hip dislocation.



| Size | | S | M | L | XL | XXL |
|------|--|-------|-------|-------|-------|-------|
| | | 27 cm | 32 cm | 37 cm | 42 cm | 46 cm |
| | | 32 cm | 37 cm | 42 cm | 46 cm | 53 cm |



BA 416 P


PAVLIC BANDAGE

Functions and Features:

- Positions both hip joints at flexion and abduction.
- Avoids adduction and extension.
- This position and movement ensure spontaneous reduction of the hip, and acetabular development in displaced hip.
- Ensures femur head to stay in acetabulum during free movements of the child.
- Easy-to-wear and remove; made of velcro bands and cotton column.

Indication:

- Hip displacement, subluxation.



| Size | | S | M | L | XL |
|------|--|-------|-------|-------|-------|
| | | 24 cm | 34 cm | 43 cm | 52 cm |
| | | 34 cm | 43 cm | 52 cm | 61 cm |



BA 411 P

HIP-X CRUSER HIP ABDUCTION BRACE

Function and Features:

- Positions the hip joint at flexion and abduction.
- Ensures femur head to stay in acetabulum.
- External part is made of polyethylene while the internal part is made of air-permeable spongy fabric.
- All materials can be washed. Easy-to-apply, light.
- Single orthoses would be enough during treatment, as the adjustment room is wide.

Indication:

- Coxitis, hipdislocation.



| Size | | S | M | L |
|------|--|-------|-------|-------|
| | | 32 cm | 37 cm | 42 cm |
| | | 37 cm | 42 cm | 46 cm |





BA 409 P

PAVLIK HARNESS TUBINGER TYPE

Functions and Features:

- It has straps that are fastened around the baby's legs and held up by shoulder and chest straps.
- This holds the hips and knees up with the legs apart. This is the best position for the hip joint to be in. It allows contact between the thigh and hip bones and helps strengthen the muscles and ligaments of the hip while it is developing.
- Soft TPR material to make baby feel comfortable.
- Adjustable pearl strap fits most of the babies.



| Size | 1 | 2 | 3 |
|------|-------|-------|-------|
| | 24 cm | 43 cm | 61 cm |
| | 43 cm | 61 cm | 80 cm |

HIPFIX (Post-Operative Hip Abduction Orthosis)

Functions and Features:

- It is an anatomically designed hip abduction orthosis made of thermoplastic material, which stabilizes the hip joint.
- Internal surface that touches the skin is covered with towel cloth.
- Flexion can be adjusted up to 120 degree, using 15 degree steps
- Abduction angle can be adjusted.

Indication:

- Ensures stabilization in total hip prosthesis or revision.
- Used in femur head resection,
- After femoral or acetabular osteotomy,
- Recurrent hip dislocation cases.



Size Universal-Right / Left

BA 415





KNEE ORTHOSIS

Biofix[®]

Orthopedic Soft Products



BA 20103-1

GENUFIX STABLE (Knitted Kneepad With Hinge)

Function and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect.
- Supports and stabilizes the knee.
- Medial and lateral joint supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.

Indication:

- Collateral ligament injuries. Knee instability.
- ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



| Size | | S | M | L | XL | XXL |
|------|--------|-------|-------|-------|-------|-------|
| | Top | 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| | Bottom | 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |



BA 20103-2

BA 30301

GENUAIRTEC STABLE (Kneepad with Hinge)

Functions and Features:

- Ensures medial and lateral stabilization and support.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.

Indication:

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



| Size | | S | M | L | XL | XXL |
|------|--------|-------|-------|-------|-------|-------|
| | Top | 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| | Bottom | 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |



BA 30301 P

GENUAIRTEC STABLE (Kneepad with Hinge) PEDIATRIC

Functions and Features:

- Ensures medial and lateral stabilization and support.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.

Indication:

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



| Size | | XS | XXS | XXXS |
|------|--------|-------|-------|-------|
| | Top | 24 cm | 21 cm | 19 cm |
| | Bottom | 26 cm | 23 cm | 21 cm |



BD 40301 U

UNIFIX STABIL (Universal Knee Pad with Hinge)

Functions and Features:

- Supports and stabilizes the knee.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of 3.5 mm spongy Neoprene fabric as one size.

Indication:

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



Size Universal



BA 400

IMMOFIX (Knee Mobilizer)

Functions and Features:

- Ensures immobilization at 0 degree extension of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of three different panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

Indication:

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



Size

| S | M | L |
|-------|-------|-------|
| 28 cm | 31 cm | 34 cm |
| 31 cm | 34 cm | 37 cm |



S - 48 cm
M - 56 cm
L - 62 cm



BA 400 P

IMMOFIX (Knee Mobilizer) PEDIATRIC

Functions and Features:

- Ensures immobilization at 0 degree extension of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of three different panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

Indication:

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



Size

| XS | XXS | XXXS |
|-------|-------|-------|
| 24 cm | 21 cm | 19 cm |
| 26 cm | 23 cm | 21 cm |



XS - 40 cm
XXS - 30 cm
XXXS - 20 cm





BA 401

ECO KNEE IMMOBILIZER

Functions and Features:

- Ensures immobilization at 0 degree extension of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of one panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

Indication:

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.

| Size | S | M | L |
|------|----------------|----------------|----------------|
| | 28 cm 31 cm | 31 cm 34 cm | 34 cm 37 cm |

| | |
|--|-----------|
| | S - 48 cm |
| | M - 56 cm |
| | L - 62 cm |

GENUROM ECO LONG (Rom Ortois with Adjustable Angle)

Functions and Features:

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

Indication:

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move within a certain range. In cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture



Size Universal

BA 524



BA 524 K

GENUROM ECO SHORT (Rom Ortois with Adjustable Angle)

Functions and Features:

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg, side bars can be adjusted according to the height of the patient.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

Indication:

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move within a certain range, in cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture.



Size Universal



BA 525

GENUROM COMFORT (Rom Knee Ortois with Adjustable Angle)

Functions and Features:

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg. Side bars can be adjusted according to the height of the patient.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

Indication:

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move with in a certain range, in cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture.



Size Universal



BA 20101-1

**GENUFIX PLUS
(Ligament and Patella Supported Knitted Knee Pad)**

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect.
- One spiral baleen is available in medial and lateral

Indication:

- Advanced twists on the knee, swelling, post-operative and post-traumatic irritations, patella subluxation Chondromalacia, Patellar tendinitis



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |

BA 20101-2



BA 30201

GENUAIARTEC PLUS (Ligament and Patella Supported Knee Pad)

Function and Features:

- Ensures patella and ligament stabilization and support.
- Made of air-permeable, skin-friendly airtec (spacer) fabric that avoids sweating.
- Silicon support reduces pain and edema thanks to its massage effect.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization.

Indication:

- Advanced twists on the knee, swelling,
- Post-operative and post-traumatic irritations,
- Patella subluxation, Chondromalacia, Patellar tendinitis



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |





BA 30201 P

GENUAIARTEC PLUS (Ligament and Patella Supported Knee Pad)

Function and Features:

- Ensures patella and ligament stabilization and support.
- Made of air-permeable, skin-friendly airtec (spacer) fabric that avoids sweating.
- Silicon support reduces pain and edema thanks to its massage effect.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization.

Indication:

- Advanced twists on the knee, swelling,
- Post-operative and post-traumatic irritations,
- Patella subluxation, Chondromalacia, Patellar tendinitis



| Size | XS | XXS | XXXS |
|------|-------|-------|-------|
| | 24 cm | 21 cm | 19 cm |
| | 26 cm | 23 cm | 21 cm |

BA 30203

GENUAIARTEC LIGA (Lateral J Supported Knee Pad)

Functions and Features:

- Fixes the patella from lateral to medial using the "J" shaped silicon support.
- Spiral baleens in medial and lateral increase stability.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

Indication:

- Lateral patellar subluxation, dislocation or patella femoral disorder.



| Size (Right / Left) | S | M | L | XL | XXL |
|---------------------|-------|-------|-------|-------|-------|
| | 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| | 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |



BD 40201 U

UNIFIX LIGA (Ligament and Patella Supported Knee Pad)

Function and Features:

- Made of 3.5 mm, air-permeable, spongy Neoprene fabric.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization, and ensure additional support.

Indication:

- Advanced twists on the knee, swelling, post-operative and post-traumatic irritations, patella subluxation, Chondromalacia, Patellar tendinitis, sports activities.



Size Universal



BA 20102-1

**GENUFIX
(Patella Silicon Supported Weaved Knee Pad)**

Functions and Features:

- Supports and stabilizes the patella.
- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect

Indication:

- Advanced twists on the knee, swelling, patella subluxation, Chondromalacia, Patellar tendinitis.



| Size | | | | | |
|-------|-------|-------|-------|-------|--|
| S | M | L | XL | XXL | |
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm | |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm | |

BA 20102-2



BA 30101

GENUAIRTEC (Patella Silicon Supported Knee Pad)

Function and Features:

- Controls patilla with silicon support.
- Avoids patella compression thanks to open patella.
- Made of air-permeable, airX® fabric that avoids sweating.
- One spiral baleen is available in medial and lateral.

Indication:

- Patella subluxation, Chondromalacia, Patellar tendinitis



| Size | | | | | |
|-------|-------|-------|-------|-------|--|
| S | M | L | XL | XXL | |
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm | |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm | |



BD 40101 U

UNIFIX PATELLA (Patella Supported Knee Pad)

Functions and Features:

- Controls patilla with silicon support.
- Avoids patella compression thanks to open patella.
- Made of 3.5 mm, air-permeable, spongy Neoprene fabric.
- Standard size. Appropriate for all patients.

Indication:

- Patella subluxation, Chondromalacia, Patellar tendinitis



Size Universal





BA 30104

GENUAIIRTEC TENDON

Functions and Features:

- Patella tendon is supported with anatomically appropriate pad.
- Reduces the pain thanks to the compression.
- Ensures extra stabilization with the upper band.
- Made of air-permeable airtec fabric that avoids sweating.

Indication:

- Osgood-schlatter disease, Patellar tendinitis, Chondromalacia



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |

BD 30105 U

GENUAIIRTEC BAND UNI.

Functions and Features:

- Patella tendon is supported with an anatomically appropriate silicon pad.
- Reduces the pain thanks to the compression.
- Made of air-permeable airtec fabric that avoids sweating.

Indication:

- Osgood-schlatter disease, Patellar tendinitis, Chondromalacia



Size Universal



BA 940-1



GENU X PATELLAR BAND

Functions and Features:

- Support is provided on the patella tendon with a silicone pad suitable for the anatomic shape.
- It is aimed to reduce pain with compression.
- Air permeable and moist regulator.

Indication:

- Patellar tendinitis
- Chondromalacia



Size

| S | L |
|-------|-------|
| 28 cm | 34 cm |
| 31 cm | 37 cm |

BA 20104-3

GENUFIX BASIC (Elastic Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

Indication:

- To relief the pain on knees, and to reduce edema.



| Size | | | | | |
|-------|-------|-------|-------|-------|--|
| S | M | L | XL | XXL | |
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm | |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm | |



BA 20104-4

GENUFIX BASIC (Elastic Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

Indication:

- To relief the pain on knees, and to reduce edema.



| Size | | | | | |
|-------|-------|-------|-------|-------|--|
| S | M | L | XL | XXL | |
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm | |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm | |



BA 20104-5

GENUFIX BASIC (Elastic Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

Indication:

- To relief the pain on knees, and to reduce edema.



| Size | | | | | |
|-------|-------|-------|-------|-------|--|
| S | M | L | XL | XXL | |
| 28 cm | 31 cm | 34 cm | 37 cm | 40 cm | |
| 31 cm | 34 cm | 37 cm | 40 cm | 43 cm | |





BA 900

KNEE SLEEVE LUX

Functions and Features:

- It is produced with circular knitting technique that provides balanced compression and can stretch in both directions.
- Thanks to equal compression, it reduces oedema and oedema related pain.
- Air permeable and moist regulator.
- Supports the circulatory system and stimulates muscle regeneration.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

Indication:

- Oedema and pain in and around Patella



| Size | S | M | L | XL | XXL |
|------|-------|-------|-------|-------|-------|
| | 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| | 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |

GENUWOOL (Wool Knee Pad)

Functions and Features:

- Keeps the knees warm.
- Made of wool, elastic fabric.

Indication:

- Used in rheumatic, sciatic, and knee pains.



| Size | S | M | L | XL | XXL |
|------|-------|-------|-------|-------|-------|
| | 28 cm | 31 cm | 34 cm | 37 cm | 40 cm |
| | 31 cm | 34 cm | 37 cm | 40 cm | 43 cm |

BA 304



BA 40402 U



UNIFIX NEOPRENE CALF SUPPORT

Functions and Features:

- Made of air-permeable, spongy neoprene fabric to support the calf.
- Ensures support with compression feature, and reduces the pain.

Indication:

- In myotonia,
- For protective purposes during and after sportive activities,
- To warm and compress the area.



Size Universal

BD 40401 U

UNIFIX NEOPRENE THIGH SUPPORT

Functions and Features:

- Made of air-permeable, spongy neoprene fabric to support the calf.
- Ensures support with compression feature, and reduces the pain.

Indication:

- In myotonia,
- For protective purposes during and after sportive activities,
- To warm and compress the area.



Size Universal





ANKLE ORTHOSIS



BA 20201-1

MALLEOFIX (Malleolar Supported Elastic Weaved Ankle Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that press around the malleolar reduce the pain through massage effect, and avoid edema.

Indication:

- Acute ankle twists, in case of any pain, Soft tissue traumas, Arthritis cases.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |

BA-20201-2



BD 40501 U

MALLEOLAR SUPPORTED ANKLE SUPPORT UNI.

Function and Features:

- Made of non-allergic, air-permeable airtec fabric that avoids sweating.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.

Indication:

- Acute ankle twists, in case of any pain,



Size Universal

BD 589 U

ANKLE SUPPORT LIGAMENT SUPPORTED WITH SPIRAL BALEENS

Functions and Features:

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate spiral baleens
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

Indication:

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size Universal



BD 40503 U

ANKLE SUPPORT LIGAMENT SUPPORTED WITH PAD UNI.

Functions and Features:

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate rigid plastic pieces.
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

Indication:

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size Universal



BA 30503

MALLEOARTEC (Ligament Supported ankle Pad)

Functions and Features:

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate rigid plastic pieces.
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

Indication:

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |



BA 20203-1

ACHILLFIX (Achilles Tendon Supported Ankle Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- The silicons that are designed in accordance with the anatomical structure of Achilles tendon
- Offers support with massage and compression effect.
- The silicon heel pads that can be placed in the shoes reduces the Achilles tension, and relieves the pain.

Indication:

- Used in acute and chronic inflammation of Achilles tendon, post-operative and post-traumatic irritation, Achillodynia, and Achillobursitis.



| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |

BA-20203-2





BD 40502 U

UNIVERSAL (Neoprene Ankle Pad)

Functions and Features:

- Made of air-permeable, Neoprene fabric.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.
- (standard size)

Indication:

- Used in acute ankle twist,
- In case of pain, and in arthritis cases.



Size Universal

BA 20204-5

MALLEOFIX BASIC (Elastic Ankle Support)

Function and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Reduces the pain and removes the edema thanks to the compression feature and massage effect.

Indication:

- Used in acute ankle twists, in case of any pain, soft tissue traumas, and arthritis cases.



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |

BA 20204-4

BA 20204-3

BA 903

ANKLE SLEEVE LUX

Function and Features:

- It is produced with circular knitting technique providing balanced compression.
- The inner surface is made of cotton outer surface polyamide thread.
- Air permeable and moist regulator.
- It reduces oedema and oedema related pain with compression.
- It supports healthy and correct walking.

Indications:

- Acute ankle sprains
- Soft tissue traumas
- Osteoarthritis, Arthritis



Size

| S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|
| 19 cm | 21 cm | 23 cm | 25 cm | 27 cm |
| 21 cm | 23 cm | 25 cm | 27 cm | 29 cm |

BA 538

MALLEOAIR (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External parts are made of plastics, and shaped anatomically.
- Internal air bladder is blown up as required, and ensures full support.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, and post-operative cases.



Size Universal



BA 589

FOAMCAST (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External plastics are shaped anatomically. Fits the ankle.
- Latex foam coveren with fabric takes the anatomical shape of the ankle.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, ad post-operative cases.

Size Universal



BA 589 P

FOAMCAST (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External plastics are shaped anatomically. Fits the ankle.
- Latex foam coveren with fabric takes the anatomical shape of the ankle.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, ad post-operative cases.

Size Universal





BA 588

FOAMCAST PLUS (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization in ankles. Easy-to-use as it allows movement in sagittal plane. External plastics are shaped anatomically.
- Fits the ankle. Latex foam covered with fabric takes the anatomical shape of the ankle.
- Ensures extra stabilization with the extra bandage that covers the ankle to metatarsus.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, and post-operative cases.

Size Universal

ECO ROM WALKER

Functions and Features:

Stabilizes the ankle with angle adjustable hinge, and steel bars.

- Dorsiflexion and plantar flexion can be adjusted through 10 degree ranges.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.

Indication:

- Used in Achilles tendon repairs, acute twist of ankles, Serious soft tissue injuries, Stable fractures of angles, and non-healing fractures in aistal tibia and fibula.

Size Universal



| | S | M | L |
|--|----|----|----|
| | 36 | 39 | 42 |
| | 39 | 42 | 45 |

BA 519



BA 522



AIR ACCILLES BOOT LONG

Functions and Features:

- Ankle plantar flexium height can be adjusted using the pads on the heel.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.
- Air Inflatable

Indication:

- Used in Achilles tendon repairs, Acute twist of ankles, Serious soft tissue injuries, Stable fractures in ankles, And non-healing fractures in aistal tibia and fibula.



| Size | S | M | L |
|------|----|----|----|
| | 36 | 39 | 42 |
| | 39 | 42 | 45 |

BA 522 K

AIR ACCHILLES BOOT SHORT

Functions and Features:

- Ankle plantar flexion height can be adjusted using the pads on the heel.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.
- Air Inflatable

Indication:

- Used in Achilles tendon repairs, Acute twist of ankles, Serious soft tissue injuries, Stable fractures in ankles, And non-healing fractures in aistal tibia and fibula.



| Size | S | M | L |
|------|----|----|----|
| → ← | 36 | 39 | 42 |
| → ← | 39 | 42 | 45 |



BA 420

FOOT GUARD (Foot Position Splint)

Function and Features:

- Made of polypropylene material.
- Keeps the ankle at 90°.
- Used in hemiplegia, paraplegia, peroneal nerve paralysis, and ankle instability.
- Avoids the foot to experience plantar flexion in patients in intensive-care unit.
- Furthermore, the portable T band under the foot avoids it to tend to right or left.
- The T band can be removed when required.
- Avoids eschar injuries.



| Size | S | M | L |
|------|----|----|----|
| → ← | 36 | 39 | 42 |
| → ← | 39 | 42 | 45 |



BA 536

AFO

Function and Features:

- An orthosis made of polypropylene that can be worn with shoes and is used to keep the ankle angle at 90 degrees.
- Used in mild foot deformities caused by nerve injuries such as hemiplegia and disc hernia, and by peroneal nerve disorders.
- Has three different sizes.
- Right and left options are available.
- State the size of your calf.



| Size (Right - Left) | S | M | L |
|---------------------|----|----|----|
| → ← | 36 | 39 | 42 |
| → ← | 39 | 42 | 45 |



BA 583

ROCKER PLASTER SLIPPER

Functions and Features:

- Makes it easy to walk on your plaster covered feet thanks to its special base (Rocker).
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.



| Size | | | |
|------|----|----|--|
| S | M | L | |
| 36 | 39 | 42 | |
| 39 | 42 | 45 | |



BA 583-2

PLASTER SLIPPER

Functions and Features:

- Makes it easy to walk on your plaster covered feet thanks to its special base (Rocker).
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.



| Size | | | |
|------|----|----|--|
| S | M | L | |
| 36 | 39 | 42 | |
| 39 | 42 | 45 | |



BA 108

HALLUX VALGUS SHOES

Functions and Features:

- Post-operative shoe designed in accordance with the anatomical form of the foot for the purpose of supporting mobilisation after Hallux Valgus surgery.
- It is manufactured from skin friendly, antibacterial materials. Its nonskid, platform sole speeds up the healing process preventing pressure on the operated region.
- The foot contact surface is soft. It does not put stress on metatarsus and phalanges, normalizes the extra stress on the heel as well. The upper surface is totally made of breathable fabric.
- Specially designed hook and loop fastener allows it to be used on both feet. It reduces possible post-operative edema on foot.



| Size | | | |
|------|----|----|--|
| S | M | L | |
| 36 | 39 | 42 | |
| 39 | 42 | 45 | |





BA 701

SILICONE INSOLE

Function and Features:

- Used for painful plantars and heels, hyperkeratosis, plantar fasciitis, and as pain absorber.
- Avoids the pains caused by certain disorders such as pescausus.
- Avoids partial pressure and tension on the feet due to remaining standing position for a long time, and posture disorders.
- As the density is adjusted according to the foot anatomy (blue parts are less dense), the product protects sensitive and boned areas by increasing the adaptation to plantar.
- Absorbs the impacts on knee and ankle joints in dynamic activities such as jogging, running, and jumping.
- Avoids the pain likely to occur on the same joint due to imbalanced pressure.
- Can be used during sports, daily life, and all activities.
- Made of 100% silicone, and is hypoallergenic.
- Physical life can be prolonged by washing the product with cold water every other day, and drying.



Size


| XS | S | M | L | XL | XXL |
|----|----|----|----|----|-----|
| 32 | 35 | 38 | 41 | 44 | 47 |
| 34 | 37 | 40 | 43 | 46 | 49 |

BA 701-1

SILICONE INSOLE WITH ARCH SUPPORTED

Function and Features:

- Used for painful plantars and heels, hyperkeratosis, plantar fasciitis, and as pain absorber.
- Avoids the pains caused by certain disorders such as pescausus.
- Avoids partial pressure and tension on the feet due to remaining standing position for a long time, and posture disorders.
- As the density is adjusted according to the foot anatomy (blue parts are less dense), the product protects sensitive and boned areas by increasing the adaptation to plantar.
- Absorbs the impacts on knee and ankle joints in dynamic activities such as jogging, running, and jumping.
- Avoids the pain likely to occur on the same joint due to imbalanced pressure.
- Can be used during sports, daily life, and all activities.
- Made of 100% silicone, and is hypoallergenic.
- Physical life can be prolonged by washing the product with cold water every other day, and drying.



Size

| XS | S | M | L | XL | XXL |
|----|----|----|----|----|-----|
| 32 | 35 | 38 | 41 | 44 | 47 |
| 34 | 37 | 40 | 43 | 46 | 49 |

BA 546

FABRIC COVERED SILICONE SUPPORT

Function and Features:

- In addition to the features of the product coded BA 701,
- Ensures comfort as it is fabric covered, avoids odor, and reduces sweating.

Indication:

- Flat feet
- Plantar fasciitis
- Pain on the heels
- Pain on metatarsus



Size

| S | M | L | XL | XXL |
|----|----|----|----|-----|
| 35 | 38 | 41 | 44 | 47 |
| 37 | 40 | 43 | 46 | 49 |

BA 682

SILICONE HEEL PAD

Function and Features:

- Used for heel support.
- Raises the heel in achilles tendon operations, and reduces the stress on achilles ligament.
- Can be used in jogging as well.
- Can be easily used in training and daily shoes.
- The product is portable, and should be used in all shoes.
- It is appropriate for the patients receiving heel treatment to use the products in the house.
- Made of 100% silicone. The product is orthopedic and hypoallergenic.



Size

| S | M | L |
|----|----|----|
| 35 | 38 | 41 |
| 37 | 40 | 43 |



BA 683

SILICONE INTERDIGITAL APPARATUS

Functions and Features:

- Preserves the health gained (treatment process) using Hallux Valgus Night Splint in posture disorders of the toe.
- Avoids the irritation caused by overlapping toes.
- Made of silicon material.
- Anti-allergenic.
- Easily washed with water (cold water is recommended).
- Furthermore, ease of use is ensured if powdered with talcum powder.



Size

| S | M | L |
|----------|----------|----------|
| 35 37 | 38 40 | 41 43 |



BA 537

SILICONE METATARSUS PILLOW

Functions and Features:

- Made of 10% gel silicone.
- Maximum elastic.
- The product is orthopedic and hypoallergenic.
- Relieves the pain in metatarsus areas, and avoid rash.
- Passed through the second toe and placed on the area where pain and rash is felt.
- Avoids callus on the same area.



Size Universal



BA 105

TRANSVERS ARCH SUPPLEMENT

Functions and Features:

- Used to support longitudinal arch in patients with pes planus or fallen arches



Size (1 - 2)



BA 106

METATARTUS PAD

Functions and Features:

- Used to support metatarsal arch in patients with pes planus or fallen arches



Size (1 - 2)





BA 685

BUNION SHIELD

Function and Features:

- Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



Size **Universal**

BA 541

EPIN HEEL SUPPORT

Functions and Features:

- Serve as a pad on the heel, and reduces the stress on the Achilles tendon;
- Reduces the irritation and pain (Epin Calcaneus) on the heel.
- Made of polyurethane, and has three different sizes.
- Bottom part is 100% rubber, not sponge.



BA 703



CORK INSOLE WITH EPIN CALCANEUS

Function and Features:

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



Size

33-34 | 35-36 | 37-38 | 39-40 | 41-42 | 43-44

BA 702

CORK SUPPORT

Function and Features:

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



Size

33-34 | 35-36 | 37-38 | 39-40 | 40-41 | 42-43 | 44-45



BA 800

INSOLE WITH STEEL SUPPORT

Functions and Features:

- Used in the treatments of Pes Planus, Pes Cavus, Metatarsalgia, and Epin Calcainen deformities.
- Supported with stainless steel, and has apparent supports.
- Preserves elasticity and dynamism for a long time.
- The outer surface is covered with artificial leather.



Size

33 - 34 | 35 - 36 | 37 - 38 | 39 - 40 | 41 - 42 | 43 - 44



BA 544 U

ELASTIC HALLUX VALGUS NIGHT SPLINT

Functions and Features:

- Avoids Hallux Valgus deformity to progress with the pull applied on the toe.
- Can be used in shoes.

Indication:

- Used for Bunion, hallux valgus, din post-operative period.



Size Universal (Right - Left)



BA 542

HALLUX VALGUS NIGHT SPLINT

Functions and Features:

- Internal surface is plastazote, while the outer surface is rigid plastic.
- Avoids Hallux Valgus deformity to progress with three point principle.
- Used to bend toe outwards in post-operative process in the Halluks Valgus problems.
- Used in bed, not with shoes.



Size Universal



BD 542 U

HALLUX VALGUS NIGHT SPLINT UNI.

Functions and Features:

- Internal surface is plastazote, while the outer surface is rigid plastic.
- Avoids Hallux Valgus deformity to progress with three point principle.
- Used to bend toe outwards in post-operative process in the Halluks Valgus problems.
- Used in bed, not with shoes.



Size Universal




Biofix[®]

ORTHOPEDIC VISCO PILLOWS

Bamboo
Pillowcase

Memory foam
Visco



MEMORY



ERGONOMIC



ECOLOGICAL

Product Info:

- Pillow case is available. Case can be removed and washed

BA-609
VISCO NECK PILLOW



BA-610
VISCO CONTOUAR NECK PILLOW



BA-615
VISCO ANATOMIC NECK PILLOW



BA-612
VISCO SEATING WHEEL



BA-611
VISCO WAIST PILLOW



BA-612 - P
POLIURETHANE SEATING WHEEL



BA-616
VISCO ANATOMIC WAIST PILLOW



BA-613
VISCO SEATING WHEEL U FORM



BA-619
VISCO SEATING CUSHION COCCYX



BA-617
VISCO CYLINDER PILLOW



BA-618
VISCO CYLINDER PILLOW HALF



BA-620
VISCO SEATING CUSHION



BA-621
VISCO SUPPORT PILLOW



BA-614
VISCO TRAVEL NECK PILLOW



BA-622
VISCO VEHICLE NECK PILLOW



BA 548

MOBILE WALKER

Product Info:

- An appropriate walker for elderly, patients having troubles with walking, patients experienced foot paralysis, and it can be used to rehabilitate the patients receiving physiotherapy treatment.
- It can be folded thanks to the mobile joints of the angled pipes, and ensures easy and safe walking by lifting two legs and pushing forward.
- Easy-to-use.
- Offers an easy experience thanks to the high quality, soft sponge knobs appropriate for the anatomical structure of the hand.
- Height can be adjusted according to the demands of the patients, thanks to the click-fit legs.
- Does not slide on wet surface thanks to the rubber bottom supports.
- Does not scratch or damage the gentle surfaces such as parquet.
- Metal parts can be wiped with a moist duster.
- Ensures easy and safe walking by lifting two legs and pushing forward.



BA 549

ARM PIT CRUTCH

Product Info:

- Arm pit support is available. No burr likely to harm your arm pit is found on the plastic support.
- Handle support is available.
- Handle height can be adjusted.
- Height can be adjusted.



BA 550

FOREARM CRUTCH

Product Info:

- Arm support is made of rigid plastic.
- Hence, no burr likely to harm your arm is found on the plastic support.
- Made of light aluminum alloyed profile.
- Height can be adjusted.





BA 551

TRIPOD WALKING STICK

Product Info:

- Balance walking stick, tripod aluminum.
- Made of aluminum.
- Height can be adjusted.
- Handle is appropriate for the shape of the hand.
- Has three legs.
- Can be used in both sides (right and left).
- Ensures stable walking support.



BA 552

SINGLE BALANCE STICK

Product Info:

- Single Aluminum Walking Stick
- Made of aluminum.
- Height can be adjusted.
- A handle appropriate for the shape of the hand is available.

BA 553



DENNIS BROWN ORTHESIS

Functions and Features:

- Controls dorsiflexion and abduction angle of the foot by assembling special PEV boot on an aluminum bar.
- Abduction angle is stable, while the external rotation is adjusted.
- Ankle dorsiflexion angle should be adjusted by the technician according to the instructions of the doctor.

Indication:

- Pes ekinovarus (PEV), Hip Dysplasia

Size:

- The product is combined according to shoe size and shoulder width of the patient.
To order the product separately: BA 553 - Dennis Brown Splint *BA 553-1 Dennis Brown Boot.

BA 852

TOILET SEAT BOOSTER FOR CLOSET

Product Info:

- BA-853 Toilet Seat Booster (With Legs)
- Toiler Seat Booster with Cover
- Toiler seat booster is first class, white, and plastic.
- Used for the patients having trouble in standing up and sitting down after hip, knee, and waist surgery.
- Can be used in all toilets.
- Boosts the sitting height of the toilet for 11 cm.
- The toilet seat booster has key apparatus, anti-bacterial, and does not hold any bacteria.
- The product has a cover, and this cover can open and close the seat.



Privately produced for you...



Our company has been producing for orthopaedic market since 1999.

We continue to insert for secure product quality and variety.

We use appropriate material to produce healthy products.

Starting with supply of gain, we are completing all phases of production line, such as narrow fabric weaving. Flat knitting and sewing.

So maintaining full products quality, we are capable of producing semi-finished and finished good which are privately required by our costumers.



We have advantages of ;

- Height quality,
- Rich variety of products,
- Being closed to the market,
- Quick services capability,
- Reasonable prices.

NOTES

A large, light gray rectangular area with a torn paper edge effect, containing horizontal dashed lines for writing notes. The area is intended for the user to take notes on the product.





BİOAKTİF Tıbbi Ürünler ve Teknik Tekstil Sanayi Ticaret Limited Şirketi
Haraççı Mah. Bayrak Cad. No: 26 34381 Arnavutköy - İstanbul - Turkey
Phone: +90 212 567 92 34 / +90 533 215 41 84
Fax: +90 212 674 88 92 E-mail: info@bioaktif.com
www.bioaktif.com