

# Bioaktif®

Medical Compression Products



[www.bioaktif.com](http://www.bioaktif.com)



Hello everyone,

Established in 1995, Bioaktif, while operating as an orthopedic prosthesis orthosis centre, started producing orthopedic soft products in 1999 in order to reduce the dependence on imports of ready-made orthopedic products in our country and to produce new products. Meanwhile, considering that the issue of raw materials is the most important factor in producing quality products and product development, he founded Bioteks narrow weaving and textile limited company in 2005 and started the production of orthopedic elastic fabrics.

In 2007, it established its machinery for the production of flat knitted bandages and pioneered the sector.

Bioaktif, which produces orthopedic soft products under the Biofix brand, orthopedic slippers under the Footflex brand, and clothes suitable for use after plastic surgery under the Bodyflex Liposuction Garments brand, has taken its place in Türkiye and the world markets today. It provides services by adhering to innovation, quality and commitments in order to constantly improve and produce products at world standards.

Bioaktif, which has a wide production range, exports its products to 20 countries with its export network in its sector. Bioaktif, which has been developing and growing in Türkiye and the world market since 1995, continues to develop and produce world-class products.

As Bioaktif, our priority is to always exceed our customers' expectations, take new initiatives and find the most appropriate solutions to patient problems.

Yours faithfully,  
Mehmet BÜYÜKÇANGA



Varifix Medical Compression Stockings are made from high quality polyamid and elastan yarns. Has a soft texture that breathes on the skin and offer perfect fit and excellent wearing comfort. Due to the firm knitted texture is optimal for patients with weak connective tissue and for the serious treatment.



### Features:



- CCL 1,2,3
- Microfilament yarns give the stocking breathability
- High moisture transportation for pleasant wearing comfort
- Perfect fit because of soft stretch characteristics
- Easy put on and take off
- Soft, elastic toe and heel
- Ensures high wearing comfort and high efficacy
- It is produced in seven sizes and two lengths
- Opaque appearance
- Skin friendly - Tested to Oeko-Tex Standard 100



## VARIFIX BELOW KNEE VARICOSE STOCKINGS

- BA 716 Varicose stocking, below knee, open toe, beige, CCL1
- BA 717 Varicose stocking, below knee, closed toe, beige, CCL1
- BA 718 Varicose stocking, below knee, open toe, beige, CCL2
- BA 719 Varicose stocking, below knee, closed toe, beige, CCL2
- BA 740 Varicose stocking, below knee, open toe, black, CCL2
- BA 742 Varicose stocking, below knee, closed toe, black, CCL2


COLOR



Closed toe



Open toe



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50

\*Please contact us to customize for products.

## VARIFIX ABOVE KNEE VARICOSE STOCKINGS

- BA 720 Varicose stocking, above knee, open toe with silicone, beige, CCL1
- BA 721 Varicose stocking, above knee, closed toe with silicone, beige, CCL1
- BA 722 Varicose stocking, above knee, open toe with silicone, beige, CCL2
- BA 723 Varicose stocking, above knee, closed toe with silicone, beige, CCL2
- BA 741 Varicose stocking, above knee, closed toe with silicone, black, CCL2
- BA 743 Varicose stocking, above knee, open toe with silicone, black, CCL2


COLOR



Closed toe



Open toe



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	46-71	60-84	64-79

\*Please contact us to customize for products.



## VARIFIX VARICOSE PANTYHOSE

- BA 724 Varicose stocking, pantyhose, open toe, CCL1
- BA 725 Varicose stocking, pantyhose, closed toe, CCL1
- BA 726 Varicose stocking, pantyhose, open toe, CCL2
- BA 727 Varicose stocking, pantyhose, closed toe, CCL2

COLOR


## VARIFIX MATERNITY VARICOSE PANTYHOSE

- BA 726 H Varicose stocking, maternity, open toe, CCL1
- BA 727 H Varicose stocking, pantyhose, closed toe, CCL1

COLOR



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	46-71	60-84	64-79




\*Please contact us to customize for products.




**VARIFIX TRAVEL VARICOSE STOCKINGS**

BA 748 Varicose stocking, travel stocking, beige, CCL1  
BA 750 Varicose stocking, travel stocking, black, CCL1

COLOR



Closed toe



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50

\*Please contact us to customize for products.

## **VARIFIX TED ANTI EMBOLISM STOCKINGS**

Varifix Anti Embolism TED Stockings are economical solution for patient and hospitals.

Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients.

Anti Embolism stocking which are produced as knee length, thigh length and thigh length with waist belt types, has a wide measurement range. Stocking has an inspection hole which controls the blood circulation under the toe.

Different colours are given to the inspection hole of each size of stockings provide easy separation.

Thigh length anti embolism stocking has dot silicone topband which prevents the stocking from slipping and dot not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti embolism stocking.

### **Features:**

- Compression Level: 18-22 mmHg
- Colour coded inspection hole
- Ensures high wearing comfort and high efficacy
- Easy put on and take off
- Soft, elastic toe and heel
- It is produced in six sizes and three lengths
- Soft and smooth touch
- With upper inspection hole or open toe are produced upon request
- Skin friendly - Tested to Oeko-Tex Standard 100








## VARIFIX TED ANTI EMBOLISM STOCKINGS KNEE LENGTH

BA 728 Anti embolism stocking, below knee




Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50



## VARIFIX TED ANTI EMBOLISM STOCKINGS THIGH LENGTH

BA 729 Anti embolism stocking, above knee




Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77

\*Please contact us to customize for products.



## VARIFIX TED ANTI EMBOLISM STOCKINGS THIGH LENGTH WITH SILICONE BAND


BA 729 S Anti embolism stocking with silicone band, above knee



Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77

## VARIFIX TED ANTI EMBOLISM STOCKINGS WITH WAIST BELT

BA 730 Anti embolism stocking with waist belt, inguinal length, above knee



Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77

\*Please contact us to customize for products.

## COMPREFIX® LYMPHEDEMA GARMENTS

Lymphedema compression garments are a special type of stocking used to treat a condition called lymphedema, caused by a malfunction in the lymphatic system. Lymphedema is a condition in which the normal drainage of lymph fluid is disrupted and typically manifests as swelling, pain and discomfort in the legs or arms.

The recommended duration of use for this product is 6 months, after which it is recommended to renew it under the supervision of a doctor.

- Ensures soft and comfortable experience
- It offers the highest level of comfort and freedom of movement while providing the necessary working pressure for effective treatment.
- Made from a breathable material , non - irritating to the skin and requires minimal maintenance.
- Multiple options are available to customize according to the person's specific needs and provides treatment for the whole body.






**COMPREFIX BELOW KNEE LYMPHEDEMA STOCKINGS**

- CF201AD Lymphedema stocking, below knee, CCL1
- CF202AD Lymphedema stocking, below knee, CCL2
- CF203AD Lymphedema stocking, below knee, CCL3


Leg Circumference (cm)	SIZES				Custom Size
	1 (S)	2 (M)	3 (L)	4 (XL)	
cB	20-22	22-24	24-26	26-28	Accept
cC	31-34	34-37	37-40	40-43	
cD	29-32	32-35	35-38	38-41	



**COMPREFIX ABOVE KNEE LYMPHEDEMA STOCKINGS**

- CF301AG Lymphedema stocking, above knee, CCL1
- CF302AG Lymphedema stocking, above knee, CCL2
- CF303AG Lymphedema stocking, above knee, CCL3

Leg Circumference (cm)	SIZES				Custom Size
	1 (S)	2 (M)	3 (L)	4 (XL)	
cB	20-22	22-24	24-26	26-28	Accept
cC	31-34	34-37	37-40	40-43	
cD	29-32	32-35	35-38	38-41	
cG	47-51	51-55	55-59	59-63	







## COMPREFIX LYMPHEDEMA GLOVES

CF501HFA Lymphedema gloves, CCL1  
CF502HFA Lymphedema gloves, CCL2  
CF503HFA Lymphedema gloves, CCL3


Hand Circumference (cm)	SIZES				Custom Size
	1 (S)	2 (M)	3 (L)	4 (XL)	
cA / cB	19	20.5	22	23.5	Accept
cC	16	17.5	19	20.5	



## COMPREFIX LYMPHEDEMA SLEEVE

CF601CG Lymphedema sleeve, CCL1  
CF602CG Lymphedema sleeve, CCL2  
CF603CG Lymphedema sleeve, CCL3

Arm Circumference (cm)	SIZES				Custom Size
	1 (S)	2 (M)	3 (L)	4 (XL)	
cG	27-31	31-35	35-39	39-43	Accept
cE	23-25	25-28	28-31	31-34	
cD	22-24	24-27	27-30	30-33	
cC	16-18	18-20	20-22	22-24	



## FOAM COMPRESS<sup>®</sup> ADJUSTABLE COMPRESSION PRODUCTS



Foam Compress Adjustable Compression Products are compression products with short stretching properties, designed to reduce edema.

Foam Compress products can be adjusted according to the user thanks to their velcro.

The built - in primer ensures that the application is swift and uncomplicated, which enhances compliance with treatment.

### Lower Extremity



### Upper Extremity



## FOAM COMPRESS® ADJUSTABLE COMPRESSION PRODUCTS

### Lower Extremity



**BFC103**  
Thigh Wrap



**BFC102**  
Knee Wrap



**BFC101**  
Calf Wrap



**BFC100**  
Ankle Wrap

Foam Compress Adjustable Lower Extremity Compression Products can be used individually or in combination to alleviate edema in the foot, calf, knee, and thigh.


Provides a precise and custom fit for optimal effectiveness.

It is effortless to put on and take off which makes it a useful and practical product.

**FOAM COMPRESS<sup>®</sup> ADJUSTABLE COMPRESSION PRODUCTS**

**Lower Extremity Size Chart**

Circumference Measurements (cm)	SIZES			
	1 (S)	2 (M)	3 (L)	4 (XL)
cA	25-27	28-30	31-35	35-39
cB	25-32	30-38	37-44	45-52
cC	40-52	48-60	55-65	60-70
cF	50-65	60-75	70-85	85-100
cG	60-75	70-80	75-90	85-100





## FOAM COMPRESS® ADJUSTABLE COMPRESSION PRODUCTS

### Upper Extremity



**BFC105**  
Arm Wrap



**BFC104**  
Hand Wrap

Foam Compress Adjustable Upper Extremity Compression Products can be used individually or in combination to alleviate edema in the foot, calf, knee, and thigh.


Provides a precise and custom fit for optimal effectiveness.


It is effortless to put on and take off which makes it a useful and practical product.

**FOAM COMPRESS® ADJUSTABLE COMPRESSION PRODUCTS**

**Upper Extremity Size Chart**

Circumference Measurements (cm)	SIZES			
	1 (S)	2 (M)	3 (L)	4 (XL)
cA	<18	<18-20	<20-22	<22
cC	17-19	20-22	23-25	24-27
cE	26-31	30-35	34-39	39-43
cG	28-33	33-38	43-48	48-52





**Bioaktif**<sup>®</sup>  
Medical Compression Products

# Varifix<sup>®</sup>

Medical compression stockings

## Comprefix<sup>®</sup>

Lymphedema Products

## Foam Compress<sup>®</sup>

Adjustable Compression Products

 **Bioaktif**  
Support For Motion

**Bioaktif Tıbbi Ürünler ve Teknik Tekstil San. Tic. Ltd. Şti.**  
Haraççı Mah. Bayrak Cad. No: 26 Arnavutköy - İstanbul  
Tel: +90 212 567 92 34 Pbx [www.bioaktif.com](http://www.bioaktif.com)