

Medical Compression Products

Güçlü destek, Özgür hareket
Strong support, freedom of movement

 **Bioaktif**[®]
Support For Motion

Varifix®

Medical compression stockings



Destek Çorapları
Support Stockings

Uyguladığı dereceli basınç sayesinde kanın tek bir noktada birikmesini engelleyerek kan dolaşımının hızlanmasına yardımcı olur.

- Hafif varis ve ödem durumlarında
 - Uzun süreli seyahatlerde varis ve venöz tromboz oluşumunu engelleme amaçlı
 - Ortopedik ameliyatlardan sonra
- Ürünün terapötik etkisini en uygun şekilde göstermesi için doğru beden seçilmesi önemlidir.

Helps to accelerate blood circulation by preventing accumulation of blood at one point thanks to its graduated compression.

Indications:


- In presence of light varicose veins and oedemas
 - Prevention of varicose veins and venous thrombosis during long distance trips
 - After ortopaedic surgeries
- It is important to select the correct size in order for the product to exert its therapeutic effect properly.

Ürünleri özelleştirmek için lütfen bizimle iletişime geçin.


Please contact us to customize the products.



- Renk / Color
- BA-707** 70 Den. Diz altı / Below knee, open toe ●
 - BA-708** 70 Den. Diz üstü / Below knee, closed toe ●
 - BA-709** 70 Den. Külotlu / Below knee, open toe ●
 - BA-710** 140 Den. Diz altı / Below knee, closed toe ●
 - BA-711** 140 Den. Diz üstü / Below knee, closed toe ●
 - BA-712** 140 Den. Külotlu / Below knee, open toe ●



- Renk / Color
- BA-713-BA** 280 Den. Diz altı / Below knee, open toe ●
 - BA-713-BK** 280 Den. Diz altı / Below knee, closed toe ●
 - BA-714** 280 Den. Diz üstü / Below knee, closed toe ●
 - BA-715** 280 Den. Külotlu / Below knee, closed toe ●



Bacak Çevresi Leg Circumference	Bedenler / Sizes						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	22-25	23-28	26-32	28-35	30-37	32-38	34-39
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cE	32-38	33-41	37-44	39-47	41-50	43-52	45-54
cF	36-46	37-51	42-57	45-60	49-65	52-68	55-70
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79

Varifix®

Medical compression stockings



Varis Çorapları Compression Stockings

**Varifix Medikal Varis Çorapları,
yüksek kaliteli poliamid ve elastan ipliklerden üretilmiştir.
Cilt üzerinde nefes alan, mükemmel uyum ve mükemmel taşıma konforu sunan yumuşak bir dokuya sahiptir.
Sıkı örgü dokusu nedeniyle bağ dokusu zayıf olan hastalar ve ciddi tedaviler için idealdir.**

Varifix Medical Compression Stockings are manufactured from high-quality polyamide and elastane yarns. They have a soft texture that is breathable on the skin, offering an excellent fit and superior wearing comfort. Due to their tight knit structure, they are ideal for patients with weak connective tissue and those undergoing serious treatments.

Özellikler

- CCL 1,2,3
- Mikrofilament iplikler çorabın hava almasını sağlar.
- Mükemmel kullanım konforu için yüksek nem aktarımı sağlar.
- Yumuşak esneme özellikleri nedeniyle mükemmel uyum sağlar.
 - Kolayca giyilip çıkarılabilir.
 - Soft elastik burun ve topuğa sahiptir.
- Yüksek giyim konforu ve yüksek etki sağlar.
 - Yedi beden ve iki uzunlukta üretilir.
 - Opak görünüm.
- Cilt dostudur - Oeko-Tex Standard 100'e göre test edilmiştir.



Features

- CCL 1,2,3
- Microfilament yarns allow the sock to breathe.
- Provides high moisture transfer for excellent wearing comfort.
- Provides a perfect fit due to its soft stretch properties.
- Easy to put on and take off.
- Features a soft elastic toe and heel.
- Provides high wearing comfort and high impact.
- Produced in seven sizes and two lengths.
- Opaque appearance.
- Skin-friendly - tested according to Oeko-Tex Standard 100.

Ürünleri özelleştirmek için lütfen bizimle iletişime geçin.
Please contact us to customize the products.

Diz Altı Varis Çorapları / Below Knee Varicose Stockings

Varifix®
Medical compression stockings



Açık burun
Open toe



- BA-716** Diz altı, burnu açık CCL1 / Below knee, open toe CCL1
- BA-717** Diz altı, burnu kapalı CCL1 / Below knee, closed toe CCL1
- BA-718** Diz altı, burnu açık CCL2 / Below knee, open toe CCL2
- BA-719** Diz altı, burnu kapalı CCL2 / Below knee, closed toe CCL2
- BA-740** Diz altı, burnu kapalı CCL2 / Below knee, closed toe CCL2
- BA-742** Diz altı, burnu açık CCL2 / Below knee, open toe CCL2
- BA-764** Diz altı, burnu açık CCL3 / Below knee, open toe CCL3
- BA-767** Diz altı, burnu kapalı CCL3 / Below knee, closed toe, CCL3

Renk / Color



Kapalı burun
Closed toe

Ölçüler Dimensions cm	Bedenler / Sizes						
	1(XS)	2(S)	3(M)	4(L)	5(XL)	6(2XL)	7(3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50





Kapalı burun
Closed toe


- BA-720** Diz üstü, burnu açık, silikonlu CCL1 / Above knee, open toe, with silicone CCL1
BA-721 Diz üstü, burnu kapalı, silikonlu CCL1 / Above knee, closed toe, with silicone CCL1
BA-722 Diz üstü, burnu açık, silikonlu CCL2 / Above knee, open toe, with silicone CCL2
BA-723 Diz üstü, burnu kapalı, silikonlu CCL2 / Above knee, closed toe, with silicone CCL2
BA-741 Diz üstü, burnu kapalı, silikonlu CCL2 / Above knee, closed toe, with silicone CCL2
BA-743 Diz üstü, burnu açık, silikonlu CCL2 / Above knee, open toe, with silicone CCL2
BA-765 Diz üstü, burnu açık, silikonlu CCL3 / Above knee, open toe, with silicone CCL3
BA-768 Diz üstü, burnu kapalı, silikonlu CCL3 / Above knee, closed toe, with silicone CCL3

Renk / Color



Açık burun
Open toe

Ölçüler Dimensions cm	Bedenler / Sizes						
	1(XS)	2(S)	3(M)	4(L)	5(XL)	6(2XL)	7(3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79



Varis K lotlu orapları / Varicose Pantyhose

Varifix®
Medical compression stockings



- BA-724** K lotlu orap, burnu aık, CCL1 / Pantyhose, open toe, CCL1
BA-725 K lotlu orap, burnu kapalı, CCL1 / Pantyhose, closed toe, CCL1
BA-726 K lotlu orap, burnu aık, CCL2 / Pantyhose, open toe, CCL2
BA-727 K lotlu orap, burnu kapalı, CCL2 / Pantyhose, closed toe, CCL2
BA-766 K lotlu orap, burnu aık, CCL3 / Pantyhose, open toe, CCL3
BA-769 K lotlu orap, burnu kapalı, CCL3 / Pantyhose, closed toe, CCL3

Renk / Color



Hamile Varis orapları / Maternity Varicose Pantyhose




- BA-726-H** Hamile orabı, burnu aık, CCL1 / Pregnant, open toe, CCL1
BA-727-H Hamile orabı, burnu kapalı, CCL1 / Pregnant, closed toe, CCL1

Renk / Color



�c�ler Dimensions cm	Bedenler / Sizes						
	1(XS)	2(S)	3(M)	4(L)	5(XL)	6(2XL)	7(3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79



Renk / Color

BA-748 Seyehat çorabı, burnu kapalı, CCL1 / Travel socks, closed toe. CCL1



BA-750 Seyehat çorabı, burnu kapalı, CCL1 / Travel socks, closed toe. CCL1



Bedenler / Sizes	1	2	3	4	5
	37-38	39-40	41-42	43-44	45-46



Anti-Embolism
HOSPITAL

Bioaktif Anti-Emboli Çorapları hasta ve hastaneler için ekonomik bir çözümdür. Hastalara maksimum konfor ve fayda sağlamak amacıyla tasarlanmıştır. Diz altı, diz üstü ve belden kemerli olarak üretilen Anti Emboli çorapları geniş bir ölçü aralığına sahiptir.

Çorapta ayak parmaklarının altındaki kan dolaşımını kontrol eden kontrol deliği bulunmaktadır.

Her beden çorabın kontrol deliğine farklı renkler verilmiş olup, çorapların kolay ayırt edilmesi sağlanır.

Diz üstü silikonlu anti emboli çorabı, çorabın kaymasını önleyen ve nokta kompresyon profilini bozmayan noktalı silikon üst banda sahiptir.

Ayrıca belden kemerli anti emboli çorabı ile diz üstü çorap dikiminde mükemmel elastikiyet ve yumuşaklığa sahip, bacağı rahatsız etmeyen elastik bant kullanılır.

Varifix Anti-Embolism Stockings are economical solution for patient and hospitals. Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients. Anti Embolism stocking which are produced as knee length, thigh length and thigh length with waist belt types, has a wide measurement range. Stocking has an inspection hole which controls the blood circulation under the toe.

Different colours are given to the inspection hole of each size of stockings provide easy separation.

Thigh length anti embolism stocking has dot silicone topband which prevents the stocking from slipping and dot not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti embolism stocking.

Özellikler

- Kompresyon Seviyesi: 18-22 mmHg
 - Renk kodlu muayene deliği
- Yüksek giyim konforu ve yüksek etki sağlar
 - Kolayca giyilip çıkarılabilir
 - Yumuşak, elastik burun ve topuk
- Altı boyut ve üç uzunluk olarak üretilmektedir
 - Yumuşak ve pürüzsüz dokunuş
- İsteğe göre üst kontrol delikli veya burnu açık olarak da üretilmektedir
- Cilt dostudur - Oeko-Tex Standard 100'e göre test edilmiştir



Features

- Compression Level: 18-22 mmHg
- Color-coded inspection hole
- Ensures high wearing comfort and high efficacy
- Easy to put on and take off
- Soft, elastic toe and heel
- It is produced in six sizes
- Soft and smooth touch
- With upper inspection hole or open toe are produced upon request
- Skin-friendly - tested according to Oeko-Tex Standard 100.


Ürünleri özelleştirmek için lütfen bizimle iletişime geçin.

Please contact us to customize the products.

BA-728 Anti-emboli çorabı, diz altı
Anti-embolism stocking, below knee




Ölçüler Dimensions cm	Bedenler / Sizes					
	1(S)	2(M)	3(L)	4(XL)	5(2XL)	6(3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-51
cD	28-35	32-38	34-41	36-44	38-47	40-50



BA-729 Anti-emboli çorabı, diz üstü - Silikonsuz
Anti-embolism socks, knee-high - Silicone-free


Ölçüler Dimensions cm	Bedenler / Sizes					
	1(S)	2(M)	3(L)	4(XL)	5(2XL)	6(3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77






BA-729-S Anti-emboli çorabı, diz üstü - Silikonlu
 Anti-embolism socks, knee-high - with silicone

Ölçüler Dimensions cm	Bedenler / Sizes					
	1(S)	2(M)	3(L)	4(XL)	5(2XL)	6(3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77



BA-730 Anti-emboli çorabı, kasık boyu - Bel kemerli
 Anti-embolism stockings, crotch length - with waist belt.

Ölçüler Dimensions cm	Bedenler / Sizes					
	1(S)	2(M)	3(L)	4(XL)	5(2XL)	6(3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77



Comprefix®

Lymphedema Products



Lenfödem Giysileri
Lymphedema Garments

Lenfödem kompresyon giysileri, lenfatik sistemdeki bir işlev bozukluğunun neden olduğu lenfödem adı verilen hastalığı tedavi etmek için kullanılan özel giysilerdir. Lenfödem, lenf sıvısının normal drenajının bozulduğu, tipik olarak bacaklarda veya kollarda şişlik, ağrı ve rahatsızlık olarak kendini gösteren bir durumdur. Ürün için önerilen kullanım süresi 6 ay olup, sonrasında doktor kontrolünde yenilenmesi tavsiye edilir.

The Comprefix Lymphedema Garments are a special type of stocking used to treat a condition called lymphedema, caused by a malfunction in the lymphatic system. The Comprefix Lymphedema Garments are produced from polyamide and elastan yarns with using flat knitting or circular knitting techniques. Lymphedema is a condition in which the normal drainage of lymph fluid is disrupted and typically manifests as swelling, pain and discomfort in the legs or arms. The recommended duration of use for this product is 6 months, after which it is recommended to renew it under the supervision of a doctor.

Özellikler

- Soft ve rahat bir deneyim sağlar.
- Etkin tedavi için gerekli basıncı sağlarken, en üst düzeyde konfor ve hareket özgürlüğü sunar.
- Nefes alabilen bir malzemeden yapılmıştır, cildi tahriş etmez ve minimum düzeyde bakım gerektirir.
- Kişinin özel ihtiyaçlarına göre özelleştirilebilen birden fazla seçenek mevcuttur ve tüm vücuda yönelik tedavi sağlar.



Features

- Ensures soft and comfortable experience.
- It offers the highest level of comfort and freedom of movement while providing the necessary working pressure for effective treatment.
- Made from a breathable material, non-irritating to the skin and requires minimal maintenance.
- Multiple options are available to customize according to the person's specific needs and provides treatment for the whole body.

Ürünleri özelleştirmek için lütfen bizimle iletişime geçin.
Please contact us to customize the products.


Diz Altı Lenfödem Çorapları / elow Knee Lymphedema Stockings

Comprefix®
Lymphedema Products



- CF201AD** Lenfödem çorabı, diz altı, CCL1 / Lymphedema stocking, below knee, CCL1
CF202AD Lenfödem çorabı, diz altı, CCL2 / Lymphedema stocking, below knee, CCL2
CF203AD Lenfödem çorabı, diz altı, CCL3 / Lymphedema stocking, below knee, CCL3

Ölçüler Dimensions cm	Bedenler / Sizes				Özel Ölçü Custom Size
	1(S)	2(M)	3(L)	4(XL)	
cB	20-22	22-24	24-26	26-28	✓
cC	31-34	34-37	37-40	40-43	
cD	29-32	32-35	35-38	38-41	




Diz Üstü Lenfödem Çorapları / Above Knee Lymphedema Stockings

Comprefix®
Lymphedema Products



- CF301AG** Lenfödem çorabı, diz üstü, CCL1 / Lymphedema stocking, above knee, CCL1
CF302AG Lenfödem çorabı, diz üstü, CCL2 / Lymphedema stocking, above knee, CCL2
CF303AG Lenfödem çorabı, diz üstü, CCL3 / Lymphedema stocking, above knee, CCL3

Ölçüler Dimensions cm	Bedenler / Sizes				Özel Ölçü Custom Size
	1(S)	2(M)	3(L)	4(XL)	
cB	20-22	22-24	24-26	26-28	✓
cC	31-34	34-37	37-40	40-43	
cD	29-32	32-35	35-38	38-41	
cG	47-51	51-55	55-59	59-63	



Lenfödem Eldivenleri / Seamless Lymphedema Gloves

Comprefix®
Lymphedema Products



CF501HFA Lenfödem eldiveni, CCL1 / Lymphoedema glove, CCL1

CF502HFA Lenfödem eldiveni, CCL2 / Lymphedema gloves, CCL2


CF503HFA Lenfödem eldiveni, CCL3 / Lymphedema gloves, CCL3

SCF502HFA Lenfödem eldiveni, dikişsiz, parmakları açık, CCL2 / Lymphedema glove, seamless, open fingers, CCL2

SCF501HDA Lenfödem eldiveni, dikişsiz, parmaksız, CCL1 / Lymphoedema glove, seamless, fingerless, CCL1

SCF502HDA Lenfödem eldiveni, dikişsiz, parmaksız, CCL2 / Lymphoedema glove, seamless, fingerless, CCL2

Ölçüler Dimensions cm	Bedenler / Sizes				Özel Ölçü Custom Size
	1(S)	2(M)	3(L)	4(XL)	
cA / cB	19	20,5	22	23,5	✓
cC	16	17,5	19	20,5	



Lenfödem Kollukları / Lymphoedema Sleeves

Comprefix®
Lymphedema Products




CF601CG Lenfödem kolluğu, CCL1 / Lymphoedema sleeve, CCL1

CF602CG Lenfödem kolluğu, CCL2 / Lymphoedema sleeve, CCL2

CF603CG Lenfödem kolluğu, CCL3 / Lymphoedema sleeve, CCL3

Ölçüler Dimensions cm	Bedenler / Sizes				Özel Ölçü Custom Size
	1(S)	2(M)	3(L)	4(XL)	
cC	16-18	18-20	20-22	22-24	✓
cD	22-24	24-27	27-30	30-33	
cE	23-25	25-28	28-31	31-34	
cG	27-31	31-35	35-39	39-43	



BA-731 Lenfödem kolluğu, dikişsiz, başparmaklı, CCL1 / Lymphoedema sleeve, seamless, thumb cuff, CCL1

BA-732 Lenfödem kolluğu, dikişsiz, başparmaklı, CCL2 / Lymphoedema sleeve, seamless, thumb cuff, CCL2

BA-734 Lenfödem kolluğu, dikişsiz, bilek boy, CCL1 / Lymphedema sleeve, seamless, wrist length, CCL1

BA-735 Lenfödem kolluğu, dikişsiz, bilek boy, CCL2 / Lymphedema sleeve, seamless, wrist length, CCL2

BA-738 Dikişsiz kol, başparmak ve omuz askılı, CCL2 / Seamless sleeve, with thumb & shoulder strap, CCL2

BA-780 Dikişsiz kol, bileğe ve omuz askısına kadar, CCL2 CCL2 / Seamless sleeve, to wrist & shoulder strap, CCL2

EWC Easy Wrap Compression

Adjustable Compression Products



Ayarlanabilir Kompresyon Ürünleri
Adjustable Compression Products

EWC wrap Compression Ayarlanabilir Kompresyon Ürünleri, ödemi azaltmak için tasarlanmış, kısa esneme özelliğine sahip kompresyon ürünleridir. Foam Compress ürünleri cırt cırtlı yapısı sayesinde kullanıcıya göre ayarlanabilir. Dahili astar, uygulamanın hızlı ve karmaşık olmamasını sağlar, bu da tedaviye uyumu artırır.

EWC wrap Compression Adjustable Compression Products are compression products with a short stretch feature designed to reduce edema. Foam Compress products can be adjusted to the user thanks to their Velcro structure. The internal liner ensures that the application is quick and uncomplicated, which increases compliance with the treatment.

Özellikler Alt Ekstremitte

- Alt Ekstremitte Kompresyon Ürünleri ayak, baldır, diz ve uyluktaki ödemi hafifletmek için tek tek veya kombinasyon halinde kullanılabilir.
- Optimum etkinlik için hassas ve özel bir uyum sağlar.
- Giyip çıkarmanın zahmetsiz olması sebebiyle kullanışlı ve pratik bir üründür.

Üst Ekstremitte

- Üst Ekstremitte Kompresyon Ürünleri kol ve eldeki ödemi hafifletmek için tek tek veya kombinasyon halinde kullanılabilir.
- Optimum etkinlik için hassas ve özel bir uyum sağlar.
- Giyip çıkarmanın zahmetsiz olması sebebiyle kullanışlı ve pratik bir üründür.



Features Lower Extremity

- Lower Extremity Compression Products Can be used individually or in combination to relieve edema in the foot, calf, knee, and thigh.
- Provides a precise and customized fit for optimum effectiveness.
- Easy to put on and take off, making it a convenient and practical product.

Upper Extremity

- Upper Extremity Compression Products can be used individually or in combination to relieve swelling in the arms and hands.
- Provides a precise and customized fit for optimum effectiveness.
- Easy to put on and take off, making it a convenient and practical product.

Ürünleri özelleştirmek için lütfen bizimle iletişime geçin.
Please contact us to customize the products.



BFC100
Ayak Bileđi Sargısı / Ankle Wrap



BFC101
Baldır Sargısı / Calf Wrap




BFC102
Diz Sargısı / Knee Bandage



BFC103
Uyluk Sargısı / Thigh Bandage

Ölçüler Dimensions cm	Bedenler / Sizes				
	1(XS)	2(S)	3(M)	4(L)	5(XL)
cA	21-24	25-27	28-30	31-34	35-28
cB	18-25	25-32	30-38	37-44	45-54
cC	30-42	40-52	48-60	55-65	66-78
cF	40-50	50-65	60-75	70-85	80-95
cG	50-60	60-75	70-80	75-90	90-100






BFC105
Kol Sargısı / Arm Wrap



BFC104
El Sargısı / Hand Wrap

Ölçüler Dimensions cm	Bedenler / Sizes				
	1(XS)	2(S)	3(M)	4(L)	5(XL)
cA	<18	18-20	20-22	22>	24>
cC	15-17	17-19	20-22	23-25	26-28
cE	22-26	26-30	29-34	33-39	40-47
cG	25-29	29-33	34-39	40-45	46-51



Varifix®
Medical compression stockings

Comprefix®
Lymphedema Products

EWC
Easy Wrap Compression
Adjustable Compression Products

 **Bioaktif®**
Support For Motion

Bioaktif Tıbbi Ürünler ve Teknik Tekstil San. Tic. Ltd. Şti.
Haraççı Mah. Bayrak Cad. No: 26 Arnavutköy - İstanbul
Tel: +90 212 567 90 34 Pbx www.bioaktif.com