

Bioaktif®

Medical Compression Products

www.bioaktif.com



Hello everyone,

Bioaktif was established in 1995, originally operated as prosthesis orthosis workshop. Our company began to produce orthopedic soft products in 1999 in order to reduce import-dependency in ready-to-use orthopedic products and to produce new ones.

In the meantime, in 2005 Bioteks narrow woven and textile limited company was established and began to produce orthopedic elastic fabrics as raw material is the most important factor in production of quality products and product development. In 2007, machinery for flatted knitted bangades was established, and the company became a pioneer in the sector.

Orthopedic soft products and orthopedic memory pillows under the brand of Biofix, orthopedic slippers under the brand of Footflex, after-surgery garments and after surgery body supports under the Bodyflex, medical compression stockings under the Varifix brand hit stores not only in Turkish market but also in the markets all around the world.

Our service is innovative and of high quality. We stick to our comitments in order to develop our company constantly and to produce world class products. Our always prior objective is to go beyond the expectations of our customers, to maximize the customer satisfaction and find the most appropriate solutions to patient problems.

Yours faithfully,
Mehmet BÜYÜKÇANGA



Varifix Medical Compression Stockings are made from high quality polyamid and elastan yarns. Has a soft texture that breathes on the skin and offer perfect fit and excellent wearing comfort. Due to the firm knitted texture is optimal for patients with weak connective tissue and for the serious treatment.

Features:

- Class 1: 18 - 21 mmHg
- Class 2: 23 - 32 mmHg
- Class 3: 34 - 46 mmHg
- Microfilament yarns give the stocking breathability
- Wearing comfort
- Perfect fit because of soft stretch characteristics
- Easy put on and take off
- Soft, elastic toe and heel
- Ensures high wearing comfort and high efficacy
- It is produced in seven sizes
- Skin friendly - Tested to Oeko-Tex Standard 100





VARIFIX BELOW KNEE VARICOSE STOCKINGS

- BA 716 Varicose stocking, below knee, open toe, beige, CCL1
- BA 717 Varicose stocking, below knee, closed toe, beige, CCL1
- BA 718 Varicose stocking, below knee, open toe, beige, CCL2
- BA 719 Varicose stocking, below knee, closed toe, beige, CCL2
- BA 740 Varicose stocking, below knee, closed toe, black, CCL2
- BA 742 Varicose stocking, below knee, open toe, black, CCL2
- BA 764 Varicose stocking, below knee, open toe, beige, CCL3
- BA 767 Varicose stocking, below knee, closed toe, beige, CCL3


COLOR



Closed toe



Open toe



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50

*Product customization available upon request.

VARIFIX ABOVE KNEE VARICOSE STOCKINGS

- BA 720 Varicose stocking, above knee, open toe, with silicone, beige, CCL1
 BA 721 Varicose stocking, above knee, closed toe, with silicone, beige, CCL1
 BA 722 Varicose stocking, above knee, open toe, with silicone, beige, CCL2
 BA 723 Varicose stocking, above knee, closed toe, with silicone, beige, CCL2
 BA 741 Varicose stocking, above knee, closed toe, with silicone, black, CCL2
 BA 743 Varicose stocking, above knee, open toe, with silicone, black, CCL2
 BA 765 Varicose stocking, above knee, open toe, with silicone, beige, CCL3
 BA 768 Varicose stocking, above knee, closed toe, with silicone, beige, CCL3


COLOR



Closed toe



Open toe



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	46-71	60-84	64-79

*Product customization available upon request.



VARIFIX VARICOSE PANTYHOSE

- BA 724 Varicose stocking, pantyhose, open toe, CCL1
- BA 725 Varicose stocking, pantyhose, closed toe, CCL1
- BA 726 Varicose stocking, pantyhose, open toe, CCL2
- BA 727 Varicose stocking, pantyhose, closed toe, CCL2
- BA 766 Varicose stocking, pantyhose, open toe, CCL3
- BA 769 Varicose stocking, pantyhose, closed toe, CCL3

COLOR 


VARIFIX MATERNITY VARICOSE PANTYHOSE

- BA 726 H Varicose stocking, pantyhose, open toe, CCL1
- BA 727 H Varicose stocking, pantyhose, closed toe, CCL1

COLOR 



Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	46-71	60-84	64-79



*Product customization available upon request.




VARIFIX TRAVEL VARICOSE STOCKINGS


BA 748 Varicose stocking, travel stocking, beige, CCL1

BA 750 Varicose stocking, travel stocking, black, CCL1

COLOR  



Closed toe



cD

cC

cB

Leg Circumference (cm)	SIZES						
	1 (XS)	2 (S)	3 (M)	4 (L)	5 (XL)	6 (2XL)	7 (3XL)
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50

*Product customization available upon request.

VARIFIX ANTI EMBOLISM STOCKINGS

Varifix Anti Embolism Stockings are economical solution for patient and hospitals.

Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients.

Anti Embolism stocking which are produced as knee length, thigh length and thigh length with waist belt types, has a wide measurement range. Stocking has an inspection hole which controls the blood circulation under the toe.

Different colours are given to the inspection hole of each size of stockings provide easy separation.

Thigh length anti embolism stocking has dot silicone topband which prevents the stocking from slipping and dot not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti embolism stocking.

Features:


- Compression Level: 15-18 mmHg
- Colour coded inspection hole
- Ensures high wearing comfort and high efficacy
- Easy put on and take off
- Soft, elastic toe and heel
- It is produced in six sizes
- Soft and smooth touch
- With upper inspection hole or open toe are produced upon request
- Skin friendly - Tested to Oeko-Tex Standard 100





VARIFIX ANTI EMBOLISM STOCKINGS KNEE LENGTH

BA 728 Anti embolism stocking, below knee




Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50



VARIFIX ANTI EMBOLISM STOCKINGS THIGH LENGTH

BA 729 Anti embolism stocking, above knee




Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77

*Product customization available upon request.



VARIFIX ANTI EMBOLISM STOCKINGS THIGH LENGTH WITH SILICONE BAND


BA 729 S Anti embolism stocking with silicone band, above knee



Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77

VARIFIX ANTI EMBOLISM STOCKINGS WITH WAIST BELT

BA 730 Anti embolism stocking with waist belt, inguinal length



Leg Circumference (cm)	SIZES					
	1 (S)	2 (M)	3 (L)	4 (XL)	5 (2XL)	6 (3XL)
cB	20-22	22-24	24-26	26-28	28-30	30-32
cC	30-37	33-40	35-43	37-46	39-49	41-52
cD	28-35	32-38	34-41	36-44	38-47	40-50
cG	46-61	49-64	52-68	55-71	58-74	61-77

*Product customization available upon request.

COMPREFIX® LYMPHEDEMA GARMENTS

The Comprefix Lymphedema Garments are a special type of stocking used to treat a condition called lymphedema, caused by a malfunction in the lymphatic system. The Comprefix Lymphedema Garments are produced from polyamide and elastan yarns with using flat knitting or circular knitting techniques.

Lymphedema is a condition in which the normal drainage of lymph fluid is disrupted and typically manifests as swelling, pain and discomfort in the legs or arms.

The recommended duration of use for this product is 6 months, after which it is recommended to renew it under the supervision of a doctor.

- Ensures soft and comfortable experience.
- It offers the highest level of comfort and freedom of movement while providing the necessary working pressure for effective treatment.
- Made from a breathable material , non - irritating to the skin and requires minimal maintenance.
- Multiple options are available to customize according to the person's specific needs and provides treatment for the whole body.

COMPREFIX® LYMPHEDEMA GARMENTS

Flat Knitting



Circular Knitting






COMPREFIX BELOW KNEE LYMPHEDEMA STOCKINGS

- CF201AD Lymphedema stocking, below knee, CCL1
CF202AD Lymphedema stocking, below knee, CCL2
CF203AD Lymphedema stocking, below knee, CCL3


Leg Circumference (cm)	SIZES				
	1 (S)	2 (M)	3 (L)	4 (XL)	Custom Size
cB	20-22	22-24	24-26	26-28	✓
cC	31-34	34-37	37-40	40-43	
cD	29-32	32-35	35-38	38-41	



COMPREFIX ABOVE KNEE LYMPHEDEMA STOCKINGS

- CF301AG Lymphedema stocking, above knee, CCL1
CF302AG Lymphedema stocking, above knee, CCL2
CF303AG Lymphedema stocking, above knee, CCL3

Leg Circumference (cm)	SIZES				
	1 (S)	2 (M)	3 (L)	4 (XL)	Custom Size
cB	20-22	22-24	24-26	26-28	✓
cC	31-34	34-37	37-40	40-43	
cD	29-32	32-35	35-38	38-41	
cG	47-51	51-55	55-59	59-63	




*Product customization available upon request.



COMPREFIX LYMPHEDEMA GLOVES

CF501HFA Lymphedema gloves, CCL1
CF502HFA Lymphedema gloves, CCL2
CF503HFA Lymphedema gloves, CCL3


Hand Circumference (cm)	SIZES				
	1 (S)	2 (M)	3 (L)	4 (XL)	Custom Size
cA / cB	19	20.5	22	23.5	✓
cC	16	17.5	19	20.5	



COMPREFIX LYMPHEDEMA SLEEVE

CF601CG Lymphedema sleeve, CCL1
CF602CG Lymphedema sleeve, CCL2
CF603CG Lymphedema sleeve, CCL3

Arm Circumference (cm)	SIZES				
	1 (S)	2 (M)	3 (L)	4 (XL)	Custom Size
cG	27-31	31-35	35-39	39-43	✓
cE	23-25	25-28	28-31	31-34	
cD	22-24	24-27	27-30	30-33	
cC	16-18	18-20	20-22	22-24	



*Product customization available upon request.



COMPREFIX SEAMLESS LYMPHEDEMA GLOVES

SCF502HFA Seamless lymphedema gloves, open finger, CCL2



COMPREFIX SEAMLESS LYMPHEDEMA ARMSLEEVES

BA 732 Seamless lymphedema sleeve, with thumb, CCL2

BA 735 Seamless lymphedema sleeve, to wrist, CCL2



COMPREFIX SEAMLESS LYMPHEDEMA ARMSLEEVES WITH SHOULDER STRAP

BA 738 Seamless sleeve, with thumb & shoulder strap, CCL2

BA 780 Seamless sleeve, to wrist & shoulder strap, CCL2

*Product customization available upon request.

FOAM COMPRESS[®] ADJUSTABLE COMPRESSION PRODUCTS



Lower Extremity



BFC103
Thigh Wrap



BFC102
Knee Wrap



BFC101
Calf Wrap



BFC100
Ankle Wrap

Foam Compress Adjustable Lower Extremity Compression Products can be used individually or in combination to alleviate edema in the foot, calf, knee, and thigh.

Provides a precise and custom fit for optimal effectiveness.

It is effortless to put on and take off which makes it a useful and practical product.

FOAM COMPRESS® ADJUSTABLE COMPRESSION PRODUCTS

Upper Extremity



BFC105 / Arm Wrap



BFC104 / Hand Wrap

Foam Compress Adjustable Upper Extremity Compression Products can be used individually or in combination to alleviate edema in the foot, calf, knee, and thigh.


Provides a precise and custom fit for optimal effectiveness.

It is effortless to put on and take off which makes it a useful and practical product.

FOAM COMPRESS® ADJUSTABLE COMPRESSION PRODUCTS


Lower Extremity Size Chart

Circumference Measurements (cm)	SIZES			
	1 (S)	2 (M)	3 (L)	4 (XL)
cA	25-27	28-30	31-35	35-39
cB	25-32	30-38	37-44	45-52
cC	40-52	48-60	55-65	60-70
cF	50-65	60-75	70-85	85-100
cG	60-75	70-80	75-90	85-100



Upper Extremity Size Chart

Circumference Measurements (cm)	SIZES			
	1 (S)	2 (M)	3 (L)	4 (XL)
cA	<18	<18-20	<20-22	<22
cC	17-19	20-22	23-25	24-27
cE	26-31	30-35	34-39	39-43
cG	28-33	33-38	43-48	48-52





Bioaktif®
Medical Compression Products

Varifix®

Medical compression stockings

Comprefix®

Lymphedema Products

Foam Compress®

Adjustable Compression Products

 **Bioaktif**
Support For Motion

Bioaktif Tıbbi Ürünler ve Teknik Tekstil San. Tic. Ltd. Şti.
Haraççı Mah. Bayrak Cad. No: 26 Arnavutköy - İstanbul
Tel: +90 212 567 92 34 Pbx www.bioaktif.com