



Orthopedic  
Soft Products

Hello everyone,

Our company was established in 1995.

Originally operated as prosthesis orthosis workshop, the company began to produce orthopedic soft products in 1999 in order to reduce import-dependency in ready-to-use orthopedic products and to produce new ones.

In the meantime, in 2005 Bioteks narrow woven and textile limited company was established and began to produce orthopedic elastic fabrics as raw material is the most important factor in production of quality products and product development.

In 2007 , machinery for flatbed knitted bangades was established, and the company became a pioneer in the sector.

Orthopedic soft products under the brand of Biofix hit stores not only in Turkish markets but also in the markets all around the world.

Our service is innovative and of high quality, and we stick to our comitments in order to develop our company constantly, and to produce world class products.

Our prior objective is to go beyond the expectations of our customers, and to maximize the customer satisfaction.

Yours faithfully,

Mehmet BÜYÜKÇANGA





## *Contents*

*Cervical Zone Orthoses*

*Shoulder and Elbow Orthoses*

*Wrist Orthoses*

*Waist and Upper Body Orthoses*

*Hip Orthoses*

*Knee Orthoses*

*Ankle Orthoses*

*Foot Health*

*Varicose Stockings*

*Waist and Neck Pillows*

*Walking Aids*



**Biofix**<sup>®</sup>  
Orthopedic Soft Products

BA-506 PHILADELPHIA CERVICAL COLLAR

**Functions and Features:**

- Ensures cervical zone support and stabilization.
- Anatomically appropriate for thorax, mandibula and occiput.
- Stabilizes cervical spines in accordance with the anatomical position.
- Consists of two parts (front and back).
- Has plastic enhancers in front and back.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcros on the sides.
- Reduced perspiration thanks to the air holes.

**Indication:**

- Cervical trauma
- Disk hernia
- Cervical stabilization
- Support and stabilization of cervical spine in pre and post-operative periods.



Size					
XS/0	S/1	M/2	L/3	XL/4	
28 cm	33 cm	38 cm	42 cm	46 cm	
32 cm	37 cm	42 cm	46 cm	50 cm	



BA-507 PHILADELPHIA COLLAR  
(OPEN TRACHEOSTOMY)

**Functions and Features:**

- Ensures cervical zone support and stabilization.
- Anatomically appropriate for thorax, mandibula and occiput.
- Stabilizes cervical spines in accordance with the anatomical position.
- Consists of two parts (front and back).
- Has plastic enhancers in front and back.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcros on the sides.
- Reduced perspiration thanks to the air holes.

**Indication:**

- Cervical trauma
- Disk hernia
- Cervical stabilization
- Support and stabilization of cervical spine in pre and post-operative periods.



Size					
XS/0	S/1	M/2	L/3	XL/4	
28 cm	33 cm	38 cm	42 cm	46 cm	
32 cm	37 cm	42 cm	46 cm	50 cm	



BA-508 NELSON PLASTOZOTE COLLAR

**Functions and Features:**

- Supports cervical zone.
- Limits the neck in the direction of flexion.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcro

**Indication:**

- Disk hernia
- Cervical muscle spasm.
- Cervical spine support.
- For cervical disc hernias, injuries, and protection in post-trauma period.



Size					
XS/0	S/1	M/2	L/3	XL/4	
28 cm	33 cm	38 cm	42 cm	46 cm	
32 cm	37 cm	42 cm	46 cm	50 cm	



BA-509 CERVICAL COLLAR WITH CHIN

**Functions and Features:**

- Ensures cervical zone support and immobilization.
- Covers the chin to limit flexion.
- Made of plastic.
- Covered with artificial leather.
- Has a foam pad appropriate for the anatomical structure of the chin.

**Indication:**

- In cervical disc hernia,
- Disk hernia
- Cervical muscle spasm, Cervical spine support
- Injuries, for recovery and protection in post-surgery period.



Size						
XS/0	S/1	M/2	L/3	XL/4	XXL/5	
28 cm	33 cm	38 cm	42 cm	46 cm	50 cm	
32 cm	37 cm	42 cm	46 cm	50 cm	54 cm	



BA-510 CERVICAL COLLAR WITHOUT CHIN

**Functions and Features:**

- Ensures cervical zone support and immobilization.
- Made of plastic, Covered with artificial leather.
- Adjustable neck height.

**INDICATION:**

- Neck pains, disk hernia, cervical muscle spasm,
- Cervical spine support.
- Adjusted according to the size of neck circle.



Size					
S/1	M/2	L/3	XL/4	XXL/4	
32 cm	36 cm	40 cm	43 cm	46 cm	
35 cm	39 cm	43 cm	46 cm	45 cm	



BA-511 FOAM COLLAR CERVICAL ORTHOSIS

**Functions and Features:**

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

**INDICATION:**

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size				
S/1	M/2	L/3	XL/4	
32 cm	36 cm	40 cm	43 cm	
35 cm	39 cm	43 cm	46 cm	



BA-511-U FOAM COLLAR CERVICAL ORTHOSIS

**Functions and Features:**

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

**INDICATION:**

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size				
S/1	M/2	L/3	XL/4	
32 cm	36 cm	40 cm	43 cm	
35 cm	39 cm	43 cm	46 cm	



6 cm  
8 cm  
10 cm  
12 cm



BA-515 NELSON COLLAR  
(FABRIC COATED)

**Functions and Features:**

- Ensures cervical zone support and stabilization. Made of hypoallergenic plastazote material.
- Covered with cotton, non-sweating fabric.

**Indication:**

- Traumatic cases, Osteoarthritis, Cervical tendinitis.
- Muscle spasm, disk hernia, for protection in post-trauma period.



Size					
XS/0	S/1	M/2	L/3	XL/4	
28 cm	33 cm	38 cm	42 cm	46 cm	
32 cm	37 cm	42 cm	46 cm	50 cm	



BA-514 **CERVITEC ADJUSTABLE  
CERVICAL ORTHOSIS**

**Functions and Features:**

- Ensures cervical zone support and stabilization.
- Chin part appropriate for the anatomical structure is produced by covering steel wire with polyurethane foam, while adjustable occipital posterior support is made of polyethylene with aluminum baleen, and covered with anti-allergic fabric.
- Adjustable flexion and extension level.
- Ensures maximum comfort for the patient, and avoids sweating and allergy.

**INDICATION:**

- Cervical trauma, disk hernia, osteoarthritis, cervical tendinitis, muscle spasm, cervical spine support in pre and post-operative period.
- Adjustable according to the size of neck circle.



Size Small(1):33-40 Large(2):41-48



BA-513 **MINEVRA CERVICAL ORTHOSIS**

**Functions and Features:**

- Ensures cervical zone support and stabilization.
- Cervical supports that support chin and occipital zone are made of plastic, and the product has body supports that support chest and back. A connection was ensured between body and cervical supports using metal pieces.
- The product gets support from the shoulder through chin and shoulder parts, and ensures extraction and fraction to the cervical zone, if required.

**Indication:**

- Used for protective purposes after Osteochondrosis, Spondylose, Spondylarthrose, Tortikolis and other operations.



**Size**

M/2	L/3
38 cm	42 cm
42 cm	46 cm



BA-413 EPIROM COMFORT

**Functions and Features:**

- An elbow orthosis with adjustable angle, designed in accordance with the anatomical elbow shape.
- Inner surface that touches the skin is made of foam-backed fabric that reduces sweating.
- The mechanism used for angle adjustments is an easy-to-use pushing system.
- Orthosis controlled elbow ensures flexion and extraction.  
Flexion (0-15-30-45-60-75-90-105) and extraction (-15-0-15-30-45-60-75-90)
- Angle can be adjusted. Orthosis length can be adjusted according to the length of arm.
- Aluminum support that avoids the wrist to bend is available and length can be adjusted.
- Shoulder belt offers comfort in use.

**Indication:**

- Support for rehabilitation in elbow contracture.
- Protection and support for rehabilitation in elbow fracture.
- Limited and controlled movement in the elbows.
- Conservative treatment after elbow dislocation, in post-op period for ligament instability.
- Post-op rehabilitation after elbow endoprosthesis.



Size (Universal Right and Left)



BA-412 HUMERUS SERMIANTO ORTHOSIS

**Functions and Features:**

- Made of 2 mm polyethylene, and covered with foam fabric on the inside.
- This fabric can be removed and washed.
- Ensures full harmony by surrounding the humerus with integrated band system.

**Indication:**

- With the aim of the protection of Humerus proximal fractures after operation, plaster treatment and stabilization.



Size (Right-Left)

	M	L
	19 cm	26 cm
	25 cm	30 cm



BA-20505 SHOULDERFIX PLUS

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Elastic tapes on the shoulder.
- Strengthens the muscles by compressing the tissues with its elastic structure, avoids swelling.
- Limits rotation and spreading movements of the shoulders.
- Avoids humerus head to distance from shoulder joint, and limits the injuring movements.

**Indication:**

- Shoulder bone fractions • Shoulder arthritics.
- Rehabilitation after shoulder implant prosthesis treatment
- Post-operative and post-traumatic irritations.



Size

S	M	L	XL	XXL
22 cm	24 cm	27 cm	30 cm	33 cm
24 cm	27 cm	30 cm	33 cm	37 cm



BA-30604 SHOULDERFIX

**Functions and Features:**

- Made of three-dimensional knitted, four-direction elastic, 3 mm thick airtec fabric.
- Elastic tapes on the shoulder.
- Velcro fastener and other plastic materials used in the product should be produced using certain raw materials that do not affect radiography.
- Strengthens the muscles by compressing the tissues with its elastic structure, avoids swelling.
- Limits rotation and spreading movements of the shoulders.
- Avoids humerus head to distance from shoulder joint, and limits the injuring movements.

**Indication:**

- Shoulder bone fractions • Shoulder arthritics
- Rehabilitation after shoulder implant prosthesis treatment
- Post-operative and post-traumatic irritations.



Size

S	M	L	XL	XXL
22 cm	24 cm	27 cm	30 cm	33 cm
24 cm	27 cm	30 cm	33 cm	37 cm






BA-140 **ACTIVE SLING 30°**  
(ARM SLING WITH PAD)

**Function and Features:**

- Used in immobilization and positioning of shoulder.
- Made of 6 mm, spongy, air-permeable, foam-backed fabric.
- The pad that abducts the shoulder at 30° is made of 25x12 foam.
- Designed to ensure comfort in long-term treatment of arm, and to offer armpit air circulation.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

**Indication:**

- Rotator cuff, capsular rupture, and bankart repair. Used in shoulder joint dislocation and soft tissue injury.

Size	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm




BA-140-A **ACTIVE SLING PLUS 30°**  
(ARM SLING WITH PAD)

**Function and Features:**

- Used in immobilization and positioning of shoulder.
- Made of 6 mm, spongy, air-permeable, foam-backed fabric.
- The pad that abducts the shoulder at 30° is made of 25x12 foam.
- Designed to ensure comfort in long-term treatment of arm, and to offer armpit air circulation.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

**Indication:**

- Rotator cuff, capsular rupture, and bankart repair. Used in shoulder joint dislocation and soft tissue injury.

Size	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm



BA-145 **SEOFIX BRACE**  
(EXTERNAL ROTATION ORTHOSIS)

**Functions and Features:**

- Used in immobilization and positioning of shoulder.
- Foam pad that ensure harmony with the patient holds the arm at 15° - 30° external rotation.
- Rotation reducer shoulder sling is available to fix the shoulder at the desired angle.
- The pad is designed in accordance with the anatomy of the body.
- Covered with air-permeable special velcro fabric.

**Indication:**

- Anterior capsular ruptures. Positioning in post-operative period.
- After posterior shoulder repairs.
- Capsular and anterior gleno humeral dislocations.
- Post-operative rotator cuff, and bankart repairs

Size (Right-Left)	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm




BA-150 **ACTIVE SLING 45° - 60°**  
(ARM SLING WITH PAD)

**Functions and Features:**

- Used in immobilization and positioning of shoulder.
- Six mm inner surface is made of foam-backed fabric that reduces sweating.
- Ensured 45° or 60° abduction angle on shoulder.
- Designed to ensure comfort in long-term treatment of arm.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

**Indication:**

- In shoulder prosthesis operations, rotator cuff, capsular rupture, and bankart repairs. Used in shoulder joint dislocation and soft tissue injury.

Size	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm



BA-110 **DELUX ARM SLING**  
(WITH WAIST BELT)

**Function and Features:**

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

**Indication:**

- Used when immobilization of arm and shoulder is required.



	XS	S	M	L	XL
	25 cm	30 cm	35 cm	40 cm	45 cm
	30 cm	35 cm	40 cm	45 cm	50 cm



BA-110-A **DELUX ARM SLING PLUS**  
(WITH WAIST BELT)

**Function and Features:**

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

**Indication:**

- Used when immobilization of arm and shoulder is required.



	XS	S	M	L	XL
	25 cm	30 cm	35 cm	40 cm	45 cm
	30 cm	35 cm	40 cm	45 cm	50 cm



BA-115 **LUX ARM SLING**  
(WITH WAIST BELT)

**Function and Features:**

- Made of air-permeable airtec spacer fabric.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Easy to wear and remove.
- Shoulder and waist belts do not stretch, and do not cause discomfort for the patient.

**Indication:**

- Used when immobilization of arm and shoulder is required.



Size	XS	S	M	L	XL
	25 cm	30 cm	35 cm	40 cm	45 cm
	30 cm	35 cm	40 cm	45 cm	50 cm



BA-116 **ECO ARM SLING**  
(WITH WAIST BELT)

**Function and Features:**

- Made of air-permeable airtec spacer fabric.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Designed to be used in both arms.

**Indication:**

- Used when immobilization of arm and shoulder is required.



Size	XS	S	M	L	XL
	25 cm	30 cm	35 cm	40 cm	45 cm
	30 cm	35 cm	40 cm	45 cm	50 cm



**BA-160 VELPOFIX**  
(VELPEAU BANDAGE)

**Functions and Features:**

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

**Indication:**

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.



**Size (Right-Left)**

XS	S	M	L	XL
60 cm	70 cm	80 cm	90 cm	100 cm
70 cm	80 cm	90 cm	100 cm	110 cm



**BA-170 CLAVICULAFIX**  
(CLAVICLE BANDAGE)

**Function and Features:**

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

**Indication:**

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



**Size**

XS	S	M	L	XL	XXL
60 cm	70 cm	80 cm	90 cm	100 cm	110 cm
70 cm	80 cm	90 cm	100 cm	110 cm	120 cm



**BA-171 CLAVICULAFIX PLUS**  
(CLAVICLE BANDAGE)

**Functions and Features:**

- Designed in accordance with anatomic structure.
- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

**Indication:**

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



**Size**

XS	S	M	L	XL	XXL
60 cm	70 cm	80 cm	90 cm	100 cm	110 cm
70 cm	80 cm	90 cm	100 cm	110 cm	120 cm



**BA-185 EPITEC ACTIVE**  
(EPICONDYLITIS ELBOW BAND)

**Functions and Features:**

- Produced as a standard size elbow band from velcro-compatible fabric, and offers epicondylitis support.
- Epicondylitis support is anatomically ensured with a silicon pad.
- The support reduces the stress on tendon, and ensures relief of pain.

**Indication:**

- Medial and lateral epicondylitis (Tennis Elbow), Golfer's Elbow



**Size (Universal)**



**BA-186 EPITEC COMFORT**  
(EPICONDYLITIS ELBOW BAND)

**Functions and Features:**

- Produced as a standard size elbow band from velcro-compatible fabric, and offers epicondylitis support.
- Epicondylitis support is anatomically ensured with a silicon pad.
- The support reduces the stress on tendon, and ensures relief of pain.

**Indication:**

- Medial and lateral epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



**BA-20301 EPITEC**  
(SILICON BACKED ELASTIC ELBOW PROTECTOR)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that apply pressure on condyls relieve the pain.

**Indication:**

- Lateral and Medial Epicondylitis (Tennis Elbow) Elbow dislocation. Chronic post-operative and post-traumatic inflammation.



Size

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



**BA-20301-B EPITEC PLUS**  
(SILICON PAD ELASTIC ELBOW WITH EXTRA BAND)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that apply pressure on condyls relieve the pain.

**Indication:**

- Lateral and Medial Epicondylitis (Tennis Elbow) Elbow dislocation. Chronic post-operative and post-traumatic inflammation.



Size

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



**BA-30603 EPITEX**  
(EPICONDYLITIS ELBOW BANDAGE WITH PAD)

**Functions and Features:**

- The bandage that ensures epicondylitis support is made of air-permeable, three dimensional airtec spacer fabric.
- Reduces edema with compression.
- Reduces sweating, and avoids allergy.
- Epicondylitis support is ensured with an extra band and pad. This support reduces the stress on tendon, and ensures relief of pain.

**Indication:**

- Medial and lateral Epicondylitis (Tennis Elbow), Golfer's Elbow



Size

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm





BA-430 MANUSTABIL SPASTIC WRIST SPLINT

**Function and Features:**

- Ensures neutral position in hands and wrists.
- Made of 3 mm polyethylene using vacuum method.
- Internal surface is covered with 3 mm plastozote.
- Durable, can be used in long-term treatments, easy-to-wear and remove.

**Indication:**

- For protections after fractures. For paralyzed patients, in wrist problems, and in cases where hand is required to be rested.



**Size (Right-Left)**

M	L
16 cm	18 cm
18 cm	20 cm



BA-431 MANU ANTI-SPASTISITE  
(ANTI-SPASTISITE WRIST SPLINT SAPARATOR)

**Functions and Features:**

- Keeps the hand, wrist, and fingers at neutral position.
- Made of 3 mm polyethylene using vacuum method.
- Internal surface is covered with 3 mm plastazote.
- Durable, can be used in long-term treatments, easy-to-wear and remove.
- Ideal for patients with spasticity due to its channeled finger design and resistance to stress.
- Ensures full harmonization and control with the special-designed tape system.

**Indication:**

- For patients with spasticity, protective after fractures. For paralyzed patients, in wrist problems, and in cases where hand is required to be rested.



**Size (Right-Left)**

M	L
16 cm	18 cm
18 cm	20 cm



BA-435 POLIETILEN DYNAMIC WRIST SPLINT

**Function and Features:**

- Manufactured of thermoplastic material with an inner terry cloth coating.
- Keeps wrist and fingers in passive extension
- Allows active finger flexion.
- In extensor tendon injuries provides early controlled motion after tendon repair and prevents the formation of adhesions.

**Indication:**

- Extensor tendon injuries



**Size (Right-Left)**

M	L
14 cm	18 cm
18 cm	20 cm



BA-436 KLEINERT WRIST SPLINT

**Function and Features:**

- Manufactured of thermoplastic material with an inner terry cloth coating.
- Keeps wrist and fingers in passive extension
- Allows passive flexion and active extensions.
- In flexor tendon injuries provides early controlled motion after tendon repair and prevents the formation of adhesions.

**Indication:**

- Flexor tendon injuries



**Size (Right-Left)**

M	L
14 cm	18 cm
18 cm	20 cm



**BA-440 MANUFIX**  
(ELASTIC WRIST SPLINT)

**Function and Features:**

- Ensures wrist support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Has anatomically shaped aluminum support on the volar surface of the wrist.
- Easy to wear and remove thanks to the velcros.

**Indication:**

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



**Size (Right-Left)**

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



**BA-441 MANUFIX LONG**  
(ELASTIC WRIST SPLINT)

**Function and Features:**

- Ensures wrist support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Has anatomically shaped aluminum support on the volar surface of the wrist.
- Easy to wear and remove thanks to the velcros.

**Indication:**

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



**Size (Right-Left)**

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



**BA-20601 MANU COMFORT**  
(KNITTED WRIST SPLINT)

**Functions and Features:**

- Produced with an advanced technology using smart knitting-waving comfort polyamide threads.
- Air-permeable, breathable. Does not hold moisture, transfers the moisture to outside.
- Hygienic, non-allergic, anti-bacterial.
- Shaped in accordance with the anatomical structure, fits the hand and wrist.
- Offers high comfort and harmony. A formable thermoplastic splint is available.
- The silicon pad on the volar side of the hand ensures support, and relieves the pain.

**Indication:**

- Tenosynovitis
- Carpal instability
- CTS
- After plaster treatment of wrist fractures
- Rheumatoid arthritis
- Ganglion pain
- Soft tissue traumas on the wrist



**Size (Right-Left)**

S	M	L	XL	XXL
14 cm	16 cm	18 cm	20 cm	22 cm
16 cm	18 cm	20 cm	22 cm	24 cm



**BA-932 MANUFIX LUX WRIST SUPPORT**

**Function and Features:**

- Air permeable and moist regulator.
- It covers the wrist, stabilizes and provides thumb movement.
- There is immobilization with ligament supports on the inner and outer sides of the palm.

**Indication:**

- Carpal Tunnel Syndrome
- Painful wrist conditions
- Wrist strain



**Size (Right-Left)**

S	M	L	XL
14 cm	20 cm		
20 cm	24 cm		



BA-30601 **MANUAIARTEC**  
(WRIST SPLINT)

**Functions and Features:**

- Ensures wrist support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Aluminum support that keeps the wrist at neutral position is available. Ensures wrist support and stabilization.
- Stability is increased with the extra tape on the wrist.

**Indication:**

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-460-1 **MANU LACE**  
(WRIST SPLINT WITH LACE)

**Functions and Features:**

- Supports sensorimotor movements, and ensures muscle stability of the wrist.
- Made of air-permeable, good-for-skin, elastic textile fabric.
- Ensures comfortable stabilization thanks to the adjustable laces that wrap the wrist.
- Limits the hand-wrist movement. Fully adjusts with the anatomical structure of the wrist.
- Keeps your wrist at neutral position, relieves you pain, helps through rehabilitation process.

**Indication:**

- Acute and chronic irritation (e.g. tendovaginitis, tendomyopathy, minor joint effusion, arthralgia).
- Post-operative and post-trauma irritation.
- Wrist Injuries.
- Rheumatoid Arthritis.
- Carpal Tunnel Syndrome
- Degenerative Wrist Disorders (e.g. arthritis)
- Sense of weakness.



Size (Right-Left)

S	L
14 cm	18 cm
16 cm	20 cm



BA-450 **MANUFIX PLUS**  
(WRIST SPLINT WITH THUMP)

**Function and Features:**

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

**Indication:**

- MCP joint disorders, CMC joint degenerations, De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-30602 **MANUAIARTEC PLUS**  
(WRIST SPLINT WITH THUMP)

**Function and Features:**

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Ensures wrist and thumb support and stabilization.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

**Indication:**

- MCP joint disorders, CMC joint degenerations
- De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm





BA-451 DEQUFIX  
(THUMP FIXING SPLINT)

**Function and Features:**

- Ensures support and stabilization of thumb joints.
- Made of elastic textile fabric.
- Has a formable aluminum support.
- Fixes CMC and MC joints.

**Indication:**

- In treatment of chronic, post-operative, and post-traumatic tissue sensitivity.
- For stabilization in ulnar cruciate ligaments.
- Carpal Tunnel Syndrome.
- Rheumatoid arthritis.
- Root arthrosis.



Size (Right-Left)



BA-452 THUMB SUPPORT  
(WITH A SPIRAL BALEEN)

**Function and Features:**

- Ensures support and stabilization of thumb joints.
- Made of elastic textile fabric. Ensures full harmony thanks to its special design.
- Spiral baleen is used for support.

**Indication:**

- In treatment of chronic, post-operative, and post-traumatic tissue sensitivity.
- For stabilization in ulnar cruciate ligaments.
- Carpal Tunnel Syndrome.
- Rheumatoid arthritis



Size (Right-Left)

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-454 MANUFIX PLUS EXTRA  
(WRIST SPLINT WITH THUMP)

**Function and Features:**

- Ensures wrist forearm and thumb support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

**Indication:**

- MCP joint disorders, CMC joint degenerations, De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-931 DEQUFIX LUX THUMB SUPPORT

**Function and Features:**

- Air permeable and moist regulator.
- Designed to cover wrist and thumb.
- Support and stabilization of thumb joints

**Indication:**

- Tendinitis
- Skier finger
- De Quervain Syndrome



Size (Right-Left)

S	M	L	XL
14 cm	20 cm		
20 cm	24 cm		



BA-455 **BOXER FRACTURE SPLINT**

**Functions and Features:**

- The Biofix Boxer Splint is an easy to adjust splint that positions and protects the 4th and 5th metacarpals in the treatment of "boxer fractures", and offers superior comfort and support to enhance compliance. The malleable palmar and dorsal aluminum stays are easily bent to immobilize the patient's hand. The products design protects and supports the fourth and fifth metacarpals without restricting flexion.

**Indications:**

- Boxer fracture



Size Universal (Right-Left)



BA-453 **THUMB SPLINT**

**Functions and Features:**

- Limits IBP / CMC joint movements.
- Made of soft and breathable fabric.
- Easy-to-use and adjustable thanks to elastic tapes.
- De Quervain tendinitis
- Gamekeeper thumb
- Collateral ligament immobilization



Size (Right-Left)



BA-118 **DYNAMIC FINGER SPLINT**

**Function and Features:**

- Made of stainless steel spring covered with eloxal.
- Points that touch the finger are made of special plastic with rounded edges.
- Produced in two sizes.
- The center point of the spring is attached using a screw, while the other edges are assembled with stainless rivet.
- Used in tendon cuts and nerve injuries.

**Indication:**

- Improvement of muscle immobility,
- Improvement of joint immobility,
- Ulnar nerve damages,
- Rehabilitation in post-operative period



Size  
S | M | L



BA-104 **MALLET FINGER SPLINT**

**Function and Features:**

- Mallet finger splint is a finger orthosis made of special material and is used to keep the DIP joint at extraction position.
- Finger is free from the point of PIP joint.
- Extensor tendon injuries.
- Finger Splint is used in fractures.
- Used in injuries and twists.

**Indication:**

- Extensor tendon treatment, mallet finger deformity. Finger bone fractures. Finger injuries. Finger nerve damages.



Size  
1 - 2 - 3 - 4 - 5 - 5,5 - 6



BA-101 FROG FINGER SPLINT

**Function and Features:**

- Used for fixing the fingers at required position.
- Internal surface is covered with soft plastazote.
- Can be formed if required.



Size  
S | M | L



BA-102 BASEBALL FINGER SPLINT

**Function and Features:**

- Used to keep the finger at extension position.
- Used in finger injuries, tendon splitting cases, and low joint cases.
- Size and numbers varies.



Size  
S | M | L



BA-183 PATIENTS HANDS AND FEET STRAP BELT

**Functions and Features:**

- Made of soft fabric with padding, breathable and comfortable to wear
- Size: adjustable, suitable for wrist or ankle circumference: 8.6"~11.4"/22~29cm
- Velcro close, easy and convenient to use
- Use this medical restraint strap to prevent patients from self-injury to ensure safety and smooth care or treatment
- Unconscious bedridden patients, mental disorder people, cognitive impairment people etc. applicable to.





BA-285 JEWETT / HYPEREXTENSION CORSET

**Functions and Features:**

- Limits the flexion of the body based on 3 point principle, and ensures hyperextension.
- Stabilizes the spine against lateral flexion and rotation movements.
- A mobile hip bar that ensure greater comfort while seated.
- Four different sizes at the side bars, which can be adjusted.
- Easy use with the B-fast lock.
- Adjustable chest bars.

**Indication:**

- Stable compression fractures in lumbar and inferior thoracic vertebra.
- Immobilization after intervertebral disc surgery.
- Axial stabilization in thoracic and osteoporotic compression fractures of lumbar spine.
- Temporary orthotic care after open surgeries on non-stable vertebral body fracture.
- Permanent orthotic care for spinal tumors and metastases which are not appropriate for surgery.
- Post-operative orthotic care after vertebral body tumor and metastases surgeries.
- Care after spinal decompression surgery with or without internal fixation.



**Size**

S	M	L
75 cm	90 cm	105 cm
90 cm	105 cm	120 cm



BA-287 TLSO CORSET FABRIC

**Functions and Features:**

- Made of rigid and air-permeable fabric.
- Limits flexion-extension of the body.
- Avoids lateral bending and rotation.
- Polyethylene rigid supports are available at the back and front sides.
- Limits extraction through thorax pad.

**Indication:**

- To mobilize the patient after surgery,
- In compression fractures, not requiring surgery.



**Size**

S	M	L
75 cm	90 cm	105 cm
90 cm	105 cm	120 cm



BA-280 DORSOFIX PLUS  
(DORSOLUMBAR CORSET)

**Functions and Features:**

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.

**Indication:**

- Can be used in compression fractures and herniated disks, for post-operative support, and in physiotherapy as a stable support.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-270 DORSOFIX  
(DORSOLUMBAR CORSET)

**Functions and Features:**

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support.

**Indication:**

- Compression fractures and herniated disks, post-operative support,
- In physiotherapy, as a stable support.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-310 POSTURFIX

**Functions and Features:**

- Ensures support to fix the posture.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Two plastic baleens are available at the back side.
- Bands on the shoulder ensure postural support, and avoid kyphoses.

**Indication:**

- Used for weak and lazy thoracic zone muscles, and ad stimulant and straightener in posture disorders.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-313 POSTURVEST

- It is very important to develop a good walking, standing and sitting posture. This posture correct belt is a good choice for you. Adopting high quality material, which is breathable and comfortable to wear. It can help you correct poor posture, stand straighter, prevention of myopia and keep a correct sitting, standing, walking posture. Such a good product, you are worth it!
- Dopting high quality material, which is breathable, durable and comfortable to wear.
- Fashionable Y-shape cross belt design on the back will get better effect,it can help to straight your back and stand upright. Ergonomic shoulder design will protect your shoulder from damaging caused by friction.
- The hook and loop on the waist belt is convenient to adjust the size and easy to wear and take off.
- Providing support for your shoulders and back, which will help to prevent humpback and curvature of the spine and prevention of myopia, correct poor postures to keep a correct sitting, standing, walking posture.
- Help to balance stress on the shoulder, back, waist, abdomen, relieve muscle fatigue and back pain, and protect the waist, back and spine to promote your image and grow healthily



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-20502 LUMBOSOFT  
(LUMBOSACRAL CORSET)

**Functions and Features:**

- Supports and stabilizes the lumbosacral region.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Easy-to-wear thanks to hand bands on the front right and left sides of the corset.
- Back height: 30 cm, Front height: 16 cm.
- Offers more rigid support and stabilization with the extra bands.

**Indication:**

- Lumbar spine instability, lumbar stress syndromes, post-operative treatment after disc surgery.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-265 LUMBOFIX

**Functions and Features:**

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbosacral.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 32 cm, Front height: 20 cm.

**Indication:**

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome, post-operative support and stabilization.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



**BA-260 LUMBOFIX PLUS**  
(LUMBOSACRAL CORSET)

**Functions and Features:**

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoid sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Produced in accordance with the anatomy of the body.
- Back height: 32 cm, Front height: 18 cm.

**Indication:**

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago - sciatic, spondylolisthesis, post-operative immobilization.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



**BA-250 LUMBOFIX**  
(LUMBOSACRAL CORSET)

**Functions and Features:**

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support, and ensure full harmonization with the body.
- Back height: 26 cm, Front height: 18 cm.

**Indication:**

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago - sciatic.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



**BA-20501 LUMBOSOFT**  
(WEAVED LUMBOSACRAL CORSET)

**Functions and Features:**

- Supports the lumbar region.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Compatible with the anatomy, and comfortable. Ensures anatomically appropriate support.
- Easy to wear and remove thanks to the velcros on the front side.

**Indication:**

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



**BA-245 LUMBOLIGHT EVO**

**Functions and Features:**

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbar.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 26 cm, Front height: 16 cm.

**Indication:**

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome,



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



**BA-246 LUMBOLIGHT EVO PLUS**  
(THE SPECIAL DESIGN EXTRA BELT ON THE WAIST STRENGTHENS THE SUPPORT)

**Functions and Features:**

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbar.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 26 cm, Front height: 16 cm.

**Indication:**

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome,



Size					
S	M	L	XL	XXL	
75 cm	83 cm	91 cm	99 cm	107 cm	
83 cm	91 cm	99 cm	107 cm	115 cm	



**BA-240 LUMBOLIGHT PLUS**  
(LUMBOSACRAL CORSET WITH AND EXTRA BELT)

**Functions and Features:**

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

**Indication:**

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



Size					
S	M	L	XL	XXL	
75 cm	83 cm	91 cm	99 cm	107 cm	
83 cm	91 cm	99 cm	107 cm	115 cm	



**BA-20503 LUMBO-X**

**Function and Features:**

- Supports the lumbar region.
- A triangle silicon pad-support appropriate for the anatomy of the body is placed on the back side of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Compatible with the anatomy, and comfortable. Silicon support appropriate for the anatomical structure is available.
- Easy to wear and remove thanks to the velcros on the front side.

**Indication:**

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome.



Size					
S	M	L	XL	XXL	
75 cm	83 cm	91 cm	99 cm	107 cm	
83 cm	91 cm	99 cm	107 cm	115 cm	



**BA-230 LUMBOUGHT**  
(LUMBOSACRAL CORSET 26 CM)

**Functions and Features:**

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two adjustable steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset offers full support.

**Indication:**

- Can be used in lumbago, sciatic, and herniated disk problems, for post-operative support, and in physiotherapy as a stable support.



Size					
S	M	L	XL	XXL	
75 cm	83 cm	91 cm	99 cm	107 cm	
83 cm	91 cm	99 cm	107 cm	115 cm	





BA-933 STRONG LUMBAR BELT WITH SUSPENDERS

**Functions and Features:**

- Lumbar belt with braces
- Stabilization and support of the lower back


**Indications:**

- During the handling of loads, labor protection in continuous overexertion or lumbar



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm




BA-302 STERNAFIX  
(STERNUM CORSET)

**Functions and Features:**

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

**Indication:**

- To support sternum zone in post-operative period.



**Size**

S	M	L	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm




BA-306 STERNUM BODY CORSET WITH STRAP

**Functions and Features:**

- Ensures immobilization after bypass operation.
- Made of highly air-permeable elastic fabric. Avoids allergy, transfers the moisture to outside.
- Can be used to increase the stabilization by tightening the corset with hand, provided that the patient is to cough after operation.
- Fixed with the shoulder bands to avoid sliding.

**Indication:**

- After bypass operation.



**Size**

S	M	L	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm




BA-307 STERNAFIX PLUS

**Functions and Features:**

- Ensures sternum support.
- Can be used to increase the stabilization by tightening the corset with hand, provided that the patient is to cough after operation.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

**Indication:**

- To support sternum zone in post-operative period.



**Size**

S	M	L	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm




BA-225 PREGNANT CORSET

**Functions and Features:**

- Helps to avoid waist pain during pregnancy.
- Designed to support the abdominal area, and to lift it upwards.
- Ensures maximum comforts thanks to the ergonomic design of the corset.
- Support can be increased or decreased with the elastic bands on the sides.

**Indication:**

- Lumbar and abdominal support during pregnancy.

Size		S	M	L	XL	XXL
	Waist	75 cm	83 cm	91 cm	99 cm	107 cm
	Hip	83 cm	91 cm	99 cm	107 cm	115 cm



BA-228 PREGNANT COMFORT

**Functions and Features:**

- This Maternity Girdle provides adequate support to pregnant women
- It supports both the back and bellies of women to prevent body aches
- The girdle is adjustable and comfortable to wear since it is breathable
- For a list of sizes and measurements, please scroll down to the description area to find the sizing chart
- Material: Cotton
- Abdominal Circumference

Size : M: 70-100cm L: 80-110cm XL: 90-120cm XXL: 100-130



BA-226 LUMBOSOFT SACRO  
(SACRO ILIAC SUPPORT BELT)

**Functions and Features:**

- SACRO-ILIAC alleviates your pain and releases tension around the sacroiliac joints. An elastic mesh material and two tensioning straps straighten your sacral bone, thereby relieving the ligaments. Additional relief is provided by two massage pads on the back. These pads target specific trigger points during movement and stimulate the ligaments, muscles and circulation.

**Indications:**

- sacroiliac joint disease (SIJ syndrome)
- sacroiliac joint arthritis
- sacroiliac joint instability
- myalgia and tendinopathies in the pelvic region
- pelvic girdle instability
- structural impairment after spondylodesis on the lumbar spine (LWS)
- recurrent sacroiliac joint blockage



BA-227 LUMBOLIGHT SACRO  
(SACRO ILIAC SUPPORT BELT)

**Functions and Features:**

- SACRO-ILIAC alleviates your pain and releases tension around the sacroiliac joints. An elastic BREATHBAL material and SPECIAL DESING tensioning straps straighten your sacral bone, thereby relieving the ligaments. Additional relief is provided by two massage pads on the back. The SILICON pads target specific trigger points during movement and stimulate the ligaments, muscles and circulation.

**Indications:**

- Sacroiliac joint disease (SIJ syndrome)
- sacroiliac joint arthritis
- sacroiliac joint instability
- myalgia and tendinopathies in the pelvic region
- pelvic girdle instability
- structural impairment after spondylodesis on the lumbar spine (LWS)
- recurrent sacroiliac joint blockage

Size		S	M	L	XL	XXL
	Waist	75 cm	83 cm	91 cm	99 cm	107 cm
	Hip	83 cm	91 cm	99 cm	107 cm	115 cm



BA-300 ABDOMINAFIX

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-3003-V ABDOMINAFIX ANATOMIC

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL
75 cm	83 cm	91 cm	99 cm
83 cm	91 cm	99 cm	107 cm



BA-3002 ABDOMINAFIX UNI.

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL
75 cm	83 cm	91 cm	99 cm
83 cm	91 cm	99 cm	107 cm



BA-3003 ABDOMINAFIX UNI.

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL
75 cm	83 cm	91 cm	99 cm
83 cm	91 cm	99 cm	107 cm



BA-3004 ABDOMINAFIX UNI.

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL
75 cm	83 cm	91 cm	99 cm
99 cm	107 cm	115 cm	



BA-301 ABDOMINAFIX UNI.

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL
75 cm	83 cm	91 cm	99 cm
99 cm	107 cm	115 cm	



BA-309 ECO ABDOMINAL

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.



**Size**

S	M	L	XL
75 cm	83 cm	91 cm	99 cm
99 cm	107 cm	115 cm	



BA-290 UMBILICAL HERNIA CORSET

**Function and Features:**

- A support that compresses the hernia area is available.
- Made of air-permeable, elastic, cotton fabric.

**Indication:**

- Used in umbilical hernias.



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-305 COLOSTOMY CORSET

**Functions and Features:**

- Manufactured of flexible and cotton fabric.
- 26 cm long corset.
- With the aid of the flexible fabric and 2 removeable plastic underwires at each sides of the stoma opening.
- Prevents any unwanted movements of the stoma bag and potential leakages which may occur during mundane and sporting activities.
- Offers perception of comfort and assurance to the patient.
- Suitable for wear underneath daily clothing.

**Indication**

- Protecting the ostomy collection bag
- Post-ostomy rehabilitation and support purposes
- Decrease risks of parastomal hernia



**Size**

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-500 SINGLE TRUSS

**Functions and Features:**

- Supports the hernia in inguinal area.
- The pad that supports the hernia is made of rubber. Covered with fabric.

**Indication:**

- Used for hernia, and in post-operative period.



Size Universal



BA-501 DOUBLE TRUSS

**Functions and Features:**

- Supports the hernia in inguinal area.
- The pad that supports the hernia is made of rubber. Covered with fabric.

**Indication:**

- Used for hernia, and in post-operative period.



Size Universal



BA-504 SUSPENSORY

**Functions and Features:**

- Avoids vesicas to lie over, and supports the ones that lied over.
- Made of cotton fabric. Has waist belt, and elastic bands that avoid sliding.

**Indication:**

- Used to support testicles in sports activities.
- Used in vesica varicoses and inflammation
- After varicosele operation.



**Size**

6	7	8	9	10
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA-414 **HIPFIX UNI.**  
(POST-OPERATIVE HIP ABDUCTION ORTHESIS)

**Functions and Features:**

- It is an anatomically designed hip abduction orthosis made of thermoplastic material, which stabilizes the hip joint.
- Internal surface that touches the skin is covered with towel cloth.
- Flexion can be adjusted up to 120 degree, using 15 degree steps.
- Abduction angle can be adjusted.

**Indication:**

- Ensures stabilization in total hip prosthesis or revision.
- Used in femur head resection,
- After femoral or acetabular osteotomy,
- Recurrent hip dislocation cases.

ONE SIZE



BA-415 **HIPFIX UNI.**  
(POST-OPERATIVE HIP ABDUCTION ORTHESIS)

**Functions and Features:**

- It is an anatomically designed hip abduction orthosis made of thermoplastic material, which stabilizes the hip joint.
- Internal surface that touches the skin is covered with towel cloth.
- Flexion can be adjusted up to 120 degree, using 15 degree steps.
- Abduction angle can be adjusted.

**Indication:**

- Ensures stabilization in total hip prosthesis or revision.
- Used in femur head resection,
- After femoral or acetabular osteotomy,
- Recurrent hip dislocation cases.



Size	
M	L
75 cm	95 cm
95 cm	125 cm





BA-524 **GENUROM ECO**  
(ROM ORTHOSIS WITH ADJUSTABLE ANGLE)

**Functions and Features:**

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

**Indication:**

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move within a certain range. In cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture



Size Universal



BA-525 **GENUROM COMFORT**  
(ROM KNEE ORTHOSIS WITH ADJUSTABLE ANGLE)

**Functions and Features:**

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg. side bars can be adjusted according to the height of the patient.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

**Indication:**

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move within a certain range, in cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture.



Size Universal



BA-400 **IMMOFIX**  
(KNEE MOBILIZER)

**Functions and Features:**

- Ensures immobilization at 0 degree extraction of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of three different panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

**Indication:**

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



**Size**

S	M	L
28 cm	31 cm	34 cm
31 cm	34 cm	37 cm



S - 48 cm  
M - 56 cm  
L - 62 cm



BA-401 **ECO KNEE IMMOBILIZER**

**Functions and Features:**

- Ensures immobilization at 0 degree extraction of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of one panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

**Indication:**

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



**Size**

S	M	L
28 cm	31 cm	34 cm
31 cm	34 cm	37 cm



S - 48 cm  
M - 56 cm  
L - 62 cm





**BA-20103 GENUFIX STABLE**  
(TEXTILE KNEEPAD WITH HINGE)

**Function and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect.
- Supports and stabilizes the knee.
- Medial and lateral joint supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.

**Indication:**

- Collateral ligament injuries. Knee instability.
- ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-30301 GENUAIRTEC STABLE**  
(KNEEPAD WITH HINGE)

**Functions and Features:**

- Ensures medial and lateral stabilization and support.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.

**Indication:**

- Collateral ligament injuries, knee instability,
- ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-30302 GENUAIRTEC STABLE OPEN**  
(OPEN KNEEPAD WITH HINGE)

**Functions and Features:**

- Supports and stabilizes the knee.
- Steel joint in medial and lateral supports the collateral ligaments, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.
- Open at the front side, easy-to-wear and remove.

**Indication:**

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-20101 GENUFIX PLUS**  
(LIGAMENT AND PATELLA SUPPORTED WEAVED KNEE PAD)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect.
- One spiral baleen is available in medial and lateral

**Indication:**

- Advanced twists on the knee, swelling, post-operative and post-traumatic irritations, patella subluxation Chondromalacia, Patellar tendinitis



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-30201 GENUAIRTEC PLUS**  
(LIGAMENT AND PATELLA SUPPORTED KNEE PAD)

**Function and Features:**

- Ensures patella and ligament stabilization and support.
- Made of air-permeable, skin-friendly airtec (spacer) fabric that avoids sweating.
- Silicon support reduces pain and edema thanks to its massage effect.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization.

**Indication:**

- Advanced twists on the knee, swelling,
- Post-operative and post-traumatic irritations,
- Patella subluxation, Chondromalacia, Patellar tendinitis



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-30203 GENUAIRTEC LIGA**  
(LATERAL J SUPPORTED KNEE PAD)

**Functions and Features:**

- Fixes the patella from lateral to medial using the "J" shaped silicon support.
- Spiral baleens in medial and lateral increase stability.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

**Indication:**

- Lateral patellar subluxation, dislocation or patella femoral disorder.



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-30204 GENUAIRTEC-X**  
(ANTERIOR CRUCIATE LIGAMENT SUPPORTED KNEE PAD)

**Function and Features:**

- Ensures stabilization in anterior cruciate ligaments, and patella.
- Offers extra support to anterior cruciate ligaments with 4 elastic bands around the patella.
- Avoids patella compression thanks to open patella.
- Made of air-permeable airtec fabric that avoids sweating.
- One spiral baleen is available in medial and lateral.

**Indication:**

- Patella subluxation, Chondromalacia, Patellar tendinitis, cruciate ligament injuries.



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



**BA-20102 GENUFIX**  
(PATELLA SILICON SUPPORTED WEAVED KNEE PAD)

**Functions and Features:**

- Supports and stabilizes the patella.
- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect

**Indication:**

- Advanced twists on the knee, swelling, patella subluxation, Chondromalacia, Patellar tendinitis.



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-30101 **GENUAIARTEC**  
(PATELLA SILICON SUPPOTED KNEE PAD)

**Function and Features:**

- Controls patilla with silicon support.
- Avoids patella compression thanks to open patella.
- Made of air-permeable, airX® fabric that avoids sweating.
- One spiral baleen is available in medial and lateral.

**Indication:**

- Patella subluxation, Chondromalacia, Patellar tendinitis



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-30104 **GENUAIARTEC TENDON**

**Functions and Features:**

- Patella tendon is supported with anatomically appropriate pad.
- Reduces the pain thanks to the compression.
- Ensures extra stabilization with the upper band.
- Made of air-permeable airtec fabric that avoids sweating.

**Indication:**

- Osgood-schlatter disease, Patellar tendinitis, Chondromalacia



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-940 **GENU X PATELLAR BAND**

**Function and Features:**

- Support is provided on the patella tendon with a silicone pad suitable for the anatomic shape.
- It is aimed to reduce pain with compression.
- Air permeable and moist regulator.

**Indication:**

- Patellar tendinitis
- Chondromalacia



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm





**BA-20201 MALLEOFIX**  
(MALLEOLAR SUPPORTED ELASTIC WEAVED ANKLE PAD)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that press around the malleolar reduce the pain through massage effect, and avoid edema.

**Indication:**

- Acute ankle twists, in case of any pain, Soft tissue traumas, Arthritis cases.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



**BA-20201-B MALLEOFIX PLUS**  
(MALLEOLAR SUPPORTED ELASTIC WEAVED ANKLE PAD AND BAND)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that press around the malleolar reduce the pain through massage effect, and avoid edema.

**Indication:**

- Acute ankle twists, in case of any pain, Soft tissue traumas, Arthritis cases.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



**BA-20203 ACCHILLFIX**  
(ACHILLES TENDON SUPPORTED ANKLE PAD)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- The silicon pads that are designed in accordance with the anatomical structure of Achilles tendon
- Offers support with massage and compression effect.
- The silicon heel pads that can be placed in the shoes reduces the Achilles tension, and relieves the pain.

**Indication:**

- Used in acute and chronic inflammation of Achilles tendon, post-operative and post-traumatic irritation, Achillodynia, and Achillobursitis.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



**BA-30501 MALLEOARTEC**  
(MALLEOLAR SUPPORTED ANKLE PAD)

**Function and Features:**

- Made of non-allergic, air-permeable airtex fabric that avoids sweating.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.

**Indication:**

- Acute ankle twists, in case of any pain,



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-30502 ELASTIC ANKLE PAD WITH 8 BANDAGES

**Function and Features:**

- Made of non-allergic, air-permeable airtec fabric that avoids sweating.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.

**Indication:**

- Used in acute ankle twist,
- In painful cases,
- And arthritis cases.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-30503 MALLEOARTEC  
(LIGAMENT SUPPORTED ANKLE PAD)

**Functions and Features:**

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate rigid plastic pieces.
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

**Indication:**

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-588 FOAMCAST PLUS  
(ANKLE STABILIZATION ORTHOSIS)

**Functions and Features:**

- Ensures medial and lateral stabilization in ankles. Easy-to-use as it allows movement in sagittal plane. External plastics are shaped anatomically.
- Fits the ankle. Latex foam covered with fabric takes the anatomical shape of the ankle.
- Ensures extra stabilization with the extra bandage that covers the ankle to metatarsus.

**Indication:**

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, and post-operative cases.

Size Universal



BA-589 FOAMCAST  
(ANKLE STABILIZATION ORTHOSIS)

**Functions and Features:**

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External plastics are shaped anatomically. Fits the ankle.
- Latex foam coveren with fabric takes the anatomical shape of the ankle.

**Indication:**

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, ad post-operative cases.

Size Universal



**BA-538 MALLEOAIR**  
(ANKLE STABILIZATION ORTHOSIS)

**Functions and Features:**

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External parts are made of plastics, and shaped anatomically.
- Internal air bladder is blown up as required, and ensures full support.

**Indication:**

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, and post-operative cases.



Universal



**BA-519 ECO ROM WALKER**

**Functions and Features:**

- Stabilizes the ankle with angle adjustable hinge, and steel bars.
- Dorsiflexion and plantar flexion can be adjusted through 10 degree ranges.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.

**Indication:**

- Used in Achilles tendon repairs, acute twist of ankles, Serious soft tissue injuries, Stable fractures of angles, and non-healing fractures in aistal tibia and fibula.



Size

	S	M	L
→	36 cm	39 cm	42 cm
←	39 cm	42 cm	45 cm



**BA-520 ACHILLES BOOT**

**Functions and Features:**

- Ankle plantar flexion height can be adjusted using the pads on the heel.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.

**Indication:**

- Used in Achilles tendon repairs, Acute twist of ankles, Serious soft tissue injuries, Stable fractures in ankles, And non-healing fractures in aistal tibia and fibula.



Size

	S	M	L
→	36 cm	39 cm	42 cm
←	39 cm	42 cm	45 cm



**BA-521 ROM WALKER**

**Functions and Features:**

- Stabilizes the ankle with angle adjustable hinge, and steel bars.
- Dorsiflexion and plantar flexion can be adjusted through 10 degree ranges.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.

**Indication:**

- Used in Achilles tendon repairs, acute twist of ankles, Serious soft tissue injuries, Stable fractures of angles, and non-healing fractures in aistal tibia and fibula.



Size

	S	M	L
→	36 cm	39 cm	42 cm
←	39 cm	42 cm	45 cm



**BA-420 FOOT GUARD  
(FOOT POSITION SPLINT)**

**Function and Features:**

- Made of polypropylene material.
- Keeps the ankle at 90°.
- Used in hemiplegia, paraplegia, peroneal nerve paralysis, and ankle instability.
- Avoids the foot to experience plantar flexion in patients in intensive-care unit.
- Furthermore, the portable T band under the foot avoids it to tend to right or left.
- The T band can be removed when required.
- Avoids eschar injuries.



Size			
S	M	L	
36 cm	39 cm	42 cm	
39 cm	42 cm	45 cm	



**BA-421 NIGHT SPLINT**

**Functions and Features:**

- Adjustable flexion bands ensures Achilles tendon and plantar fascia to stretch gently.
- Soft foam internal part can be washed, and the product is air-permeable to avoid warming up.
- Easy to wear, remove, and adjust thanks to the velcros and claps.

**Indication:**

- Plantar fasciitis
- Achilles tendinitis



Size			
S	M	L	
36 cm	39 cm	42 cm	
39 cm	42 cm	45 cm	



**BA-108 HALLUX VALGUS SHOES**

**Functions and Features:**

- Post-operative shoe designed in accordance with the anatomical form of the foot for the purpose of supporting mobilisation after Hallux Valgus surgery. It is manufactured from skin friendly, antibacterial materials. Its nonskid, platform sole speeds up the healing process preventing pressure on the operated region. The foot contact surface is soft. It does not put stress on metatarsus and phalanges, normalizes the extra stress on the heel as well. The upper surface is totally made of breathable fabric. Specially designed hook and loop fastener allows it to be used on both feet. It reduces possible post-operative edema on foot.



Size			
S	M	L	
36 cm	39 cm	42 cm	
39 cm	42 cm	45 cm	



**BA-536 AFO**

**Function and Features:**

- An orthosis made of polypropylene that can be worn with shoes and is used to keep the ankle angle at 90 degrees.
- Used in mild foot deformities caused by nerve injuries such as hemiplegia and disc hernia, and by peroneal nerve disorders.
- Has three different sizes.
- Right and left options are available.
- State the size of your calf.



Size			
S	M	L	
36 cm	39 cm	42 cm	
39 cm	42 cm	45 cm	





BA-583-2 **PLASTER SLIPPER**

**Functions and Features:**

- Makes it easy to walk on your plaster covered feet thanks to its special base
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.



Size			
S	M	L	
36 cm	39 cm	42 cm	
39 cm	42 cm	45 cm	



BA-583 **ROCKER PLASTER SLIPPER**

**Functions and Features:**

- Makes it easy to walk on your plaster covered feet thanks to its special base (Rocker).
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.



Size			
S	M	L	
36 cm	39 cm	42 cm	
39 cm	42 cm	45 cm	



BA-584 **SKIN TRACTION KIT**

**Functions and Features:**

- Ready to use kit provides non-invasive traction for non-operative fracture or post-operative limbs, including all items needed for the safe application of sustained traction weight support and suitable for use with a variety of traction methods.

Size; Adult and child



BA-105 **TRANSFER ARCH SUPPLEMENT**

**Functions and Features:**

- Used to support longitudinal arch in patients with pes planus or fallen arches



Size (1 - 2)



BA-106 **METATARTUS PAD**

**Functions and Features:**

- Used to support metatarsal arch in patients with pes planus or fallen arches



Size (1 - 2)



BA-537 SILICONE METATARSUS PILLOW

**Functions and Features:**

- Made of 10% gel silicone.
- Maximum elastic.
- The product is orthopedic and hypoallergenic.
- Relieves the pain in metatarsus areas, and avoid rash.
- Passed through the second toe and placed on the area where pain and rash is felt.
- Avoids callus on the same area.



Size Universal



BA-541 EPIN HEEL SUPPORTER

**Functions and Features:**

- Serve as a pad on the heel, and reduces the stress on the Achilles tendon;
- Reduces the irritation and pain (Epin Calcaneus) on the heel.
- Made of polyurethane, and has three different sizes.
- Bottom part is 100% rubber, not sponge.



Size ( 2 - 3 )



BA-542 HALLUX VALGUS NIGHT SPLINT

**Functions and Features:**

- Internal surface is plastazote, while the outer surface is rigid plastic.
- Avoids Hallux Valgus deformity to progress with three point principle.
- Used to bend toe outwards in post-operative process in the Hallux Valgus problems.
- Used in bed, not with shoes.



Size Universal



BA-544 ELASTIC HALLUX VALGUS NIGHT SPLINT

**Functions and Features:**

- Avoids Hallux Valgus deformity to progress with the pull applied on the toe.
- Can be used in shoes.

**Indication:**

- Used for Bunion, hallux valgus, during post-operative period.



Size Universal



BA-546 FABRIC COVERED SILICONE SUPPORT

**Function and Features:**

- In addition to the features of the product coded BA 701,
- Ensures comfort as it is fabric covered, avoids odor, and reduces sweating.

**Indication:**

- Flat feet
- Plantar fasciitis
- Pain on the heels
- Pain on metatarsus



**Size**

S	M	L	XL	XXL
35 cm	38 cm	41 cm	44 cm	47 cm
37 cm	40 cm	43 cm	46 cm	49 cm



BA-682 SILICONE HEEL PAD

**Function and Features:**

- Used for heel support.
- Raises the heel in achilles tendon operations, and reduces the stress on achilles ligament.
- Can be used in jogging as well.
- Can be easily used in training and daily shoes.
- The product is portable, and should be used in all shoes.
- It is appropriate for the patients receiving heel treatment to use the products in the house.
- Made of 100% silicone. The product is orthopedic and hypoallergenic.



**Size**

S	M	L
35 cm	38 cm	41 cm
37 cm	40 cm	43 cm



BA-685 BUNION SHIELD

**Function and Features:**

- Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

ONE SIZE



BA-683 SILICONE INTERDIGITAL APPARATUS

**Functions and Features:**

- Preserves the health gained (treatment process) using Hallux Valgus Night Splint in posture disorders of the toe.
- Avoids the irritation caused by overlapping toes.
- Made of silicon material.
- Anti-allergenic.
- Easily washed with water (cold water is recommended).
- Furthermore, ease of use is ensured if powdered with talcum powder.



**Size**

S	M	L
35 cm	38 cm	41 cm
37 cm	40 cm	43 cm



BA-701 SILICON SUPPORT

**Function and Features:**

- Used for painful plantars and heels, hyperkeratosis, plantar fasciitis, and as pain absorber.
- Avoids the pains caused by certain disorders such as pescavus.
- Avoids partial pressure and tension on the feet due to remaining standing position for a long time, and posture disorders.
- As the density is adjusted according to the foot anatomy (blue parts are less dense), the product protects sensitive and boned areas by increasing the adaptation to plantar.
- Absorbs the impacts on knee and ankle joints in dynamic activities such as jogging, running, and jumping.
- Avoids the pain likely to occur on the same joint due to imbalanced pressure.
- Can be used during sports, daily life, and all activities.
- Made of 100% silicone, and is hypoallergenic.
- Physical life can be prolonged by washing the product with cold water every other day, and drying.



**Size**

S	M	L	XL	XXL
35 cm	38 cm	41 cm	44 cm	47 cm
37 cm	40 cm	43 cm	46 cm	49 cm



BA-702-2 CORK INSOLE WITH EPIN CALCANIUS

**Function and Features:**

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



**Size:** (35 - 43)



BA-702 CORK INSOLE

**Function and Features:**

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



**Size:** (35 - 43)



BA-800 STEEL SUPPORT

**Functions and Features:**

- Used in the treatments of Pes Planus, Pes Cavus, Metatarsalgia, and Epin Calcainen deformities.
- Supported with stainless steel, and has apparent supports.
- Preserves elasticity and dynamism for a long time.
- The outer surface is covered with artificial leather.



**Size**

S	M	L
35 cm	38 cm	41 cm
37 cm	40 cm	43 cm



Varifix Classic Medical Compression Stockings are made from high quality microfilament yarn and elastan which has a soft fabric that breathes on the skin and offer perfect fit and excellent wearing comfort. Due to the firm knitted fabric this quality is optimal for patients with soft connective tissue and for the treatment of serious indications.

**Features :**

- CCL 1, 2, 3
- Elegant appearance
- Microfilament yarns give the stocking excellent breathability.
- High moisture transportation for pleasant wearing comfort
- Perfect fit because of soft stretch material. It has been manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis.
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- High wearing comfort and high efficiency
- Perfect sizing. Seven sizes and two length
- Soft and smooth touch
- Opaque appearance
- For patients with soft connective tissue
- For patient with more severe phlebological indications
- Skin friendly – Tested to Oeko-Tex Standard 100



### VARIFIX BELOW-KNEE STOCKINGS

- BA-716 Varicose stocking, below knee, open toe CCL1
- BA-717 Varicose stocking, below knee, closed toe CCL1
- BA-718 Varicose stocking, below knee, open toe CCL2
- BA-719 Varicose stocking, below knee, closed toe CCL2

COLOR  



### VARIFIX THIGH STOCKINGS

- BA-720 Varicose stocking, above knee, open toe with silicone CCL1
- BA-721 Varicose stocking, above knee, closed toe with silicone CCL1
- BA-722 Varicose stocking, above knee, open toe with silicone CCL2
- BA-723 Varicose stocking, above knee, closed toe with silicone CCL2

COLOR  



### VARIFIX PANTYHOSE


- BA-724 Varicose stocking, pantyhose, open toe CCL1
- BA-725 Varicose stocking, pantyhose, closed toe CCL1
- BA-726 Varicose stocking, pantyhose, open toe CCL2
- BA-727 Varicose stocking, pantyhose, closed toe CCL2

COLOR  



### VARIFIX MATERNITY COMPRESSION STOCKINGS


- BA-726 H
- BA-727 H



Leg Circumference	SIZES						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	22-25	23-28	26-32	28-35	30-37	32-38	34-39
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cE	32-38	33-41	37-44	39-47	41-50	43-52	45-54
cF	36-46	37-51	42-57	45-60	49-65	52-68	55-70
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79

Varifix Anti-Embolism Stockings Economical solution for patients and hospitals. Varifix Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients. Anti-embolism stocking which are produced as knee-length, thigh length and thigh length with waist belt types, has a measurement system that covers wide measurement range. Also, different colours are given to the inspection hole of each size of stockings to provide easy separation. Stocking has an inspection hole which controls the blood circulation under the foot fingers. Thigh length anti-embolism stocking has dot silicone topband which prevents the stocking from slipping and not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti-embolism stocking.

- 18-22 mmHG • Ecrú colour • Colour coded inspection hole • High wearing comfort and efficacy
- Easy put on and take off • Freedom of movements • Soft, elastic toe and heel
- Perfect sizing. Six sizes and two length • Skin friendly – Tested to Oeko-Tex Standard 100
- With upper inspection hole and open toe are produced upon request



Leg Circumference	SIZES					
	S	M	L	XL	XXL	
cB	20-22	22-24	24-26	26-28	28-30	
cB1	23-28	26-32	28-35	30-37	32-38	
cC	30-37	33-40	35-43	37-46	39-49	
cD	28-35	32-38	34-41	36-44	38-47	
cE	33-41	37-44	39-47	41-50	43-52	
cF	37-51	42-57	45-60	49-65	52-68	
cG	46-61	49-64	52-68	46-71	60-74	

BA-728  
ANTI-EMBOLISM STOCKING  
BELOW KNEE



BA-729  
ANTI-EMBOLISM STOCKING  
ABOVE KNEE



BA-729-S  
ANTI-EMBOLISM STOCKING  
ABOVE KNEE  
WITH SILICON BAND




BA-730  
ANTI-EMBOLISM STOCKING, INGUINAL



BA-750  
VARIFIX TRAVEL STOCKING



BA-748  
VARIFIX TRAVEL STOCKING



CCL- 1 (18-21 mmHg)

**BIOFIX**  
**ONE SIZE**

**PRODUCT**





BD-40501-U ANKLE SUPPORT WITH ELASTIC BANDAGE

**Function and Features:**

- Made of non-allergic, air-permeable airtec fabric that avoids sweating.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.

**Indication:**

- Acute ankle twists, in case of any pain,



Size Universal



BD-40503-U ANKLE SUPPORT LIGAMENT SUPPORTED WITH PAD

**Functions and Features:**

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate rigid plastic pieces.
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

**Indication:**

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size Universal



BD-589-U ANKLE SUPPORT LIGAMENT SUPPORTED WITH SPIRAL BALENS

**Functions and Features:**

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate spiral baleens
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtec fabric that avoids sweating.

**Indication:**

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size Universal



BD-40402-U UNIFIX CALF (LOWER CALF)

**Functions and Features:**

- Made of air-permeable, skin-friendly airtec (spacer) fabric support the lower calf.
- (standard size)
- Ensures support with compression feature, and reduces the pain.

**Indication:**

- In myotonia,
- For protective purposes during and after sportive activities,
- To warm and compress the area.



Size Universal



BD-40401-U **UNIFIX THIGH**  
(UPPER CALF)

**Functions and Features:**

- Made of air-permeable, skin-friendly airtec (spacer) fabric. to support the upper calf.
- (standard size)
- Ensures support with compression feature, and reduces the pain.

**Indication:**

- In myotonia,
- For protective purposes during and after sportive activities,
- To warm and compress the area.



Size Universal



BD-40201-U **UNIFIX LIGA**  
(LIGAMENT AND PATELLA SUPPORTED KNEE PAD)

**Function and Features:**

- Made of 3.5 mm, air-permeable, skin-friendly airtec (spacer) fabric.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization, and ensure additional support.

**Indication:**

- Advanced twists on the knee, swelling, post-operative and post-traumatic irritations, patella subluxation, Chondromalacia, Patellar tendinitis, sports activities.



Size Universal



BD-40101-U **UNIFIX PATELLA**  
(PATELLA SUPPORTED KNEE PAD)

**Functions and Features:**

- Controls patilla with silicon support.
- Avoids patella compression thanks to open patella.
- Made of 3.5 mm, air-permeable, skin-friendly airtec (spacer) fabric.
- Standard size. Appropriate for all patients.

**Indication:**

- Patella subluxation, Chondromalacia, Patellar tendinitis



Size Universal



BD-40301-U **UNIFIX STABIL**  
(UNIVERSAL KNEE PAD WITH HINGE)

**Functions and Features:**

- Supports and stabilizes the knee.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of 3.5 mm skin-friendly airtec (spacer) fabric.

**Indication:**

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



Size Universal



BD-30105-U GENUAIRTEC BAND

**Functions and Features:**

- Patella tendon is supported with an anatomically appropriate silicon pad.
- Reduces the pain thanks to the compression.
- Made of air-permeable airtec fabric that avoids sweating.

**Indication:**

- Osgood-schlatter disease, Patellar tendinitis, Chondromalacia



Size Universal




BD-240-U LUMBOLIGHT UNI.  
(LUMBOSACRAL CORSET WITH AN EXTRA BELT)

**Functions and Features:**

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

**Indication:**

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



Size

S	M	L	XL
75 cm	99 cm	99 cm	115 cm




BD-260-U LUMBOFIX UNI.  
(LUMBOSACRAL CORSET)

**Functions and Features:**

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Produced in accordance with the anatomy of the body.
- Back height: 30 cm, Front height: 18 cm.

**Indication:**

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago – sciatic, spondylolisthesis, post-operative immobilization.



Size

S	M	L	XL
75 cm	99 cm	99 cm	115 cm



BD-310-U POSTURFIX

**Functions and Features:**

- Ensures support to fix the posture.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Two plastic baleens are available at the back side.
- Bands on the shoulder ensure postural support, and avoid kyphoses.

**Indication:**

- Used for weak and lazy thoracic zone muscles, and ad stimulant and straightener in posture disorders.



Size

S	M	L	XL
75 cm	99 cm	99 cm	115 cm



BD-300-U ABDOMINAFIX

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.

Size	
S	M   L   XL
75 cm	99 cm
99 cm	115 cm



BD-302-U STERNAFIX  
(STERNUM CORSET)

**Functions and Features:**

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

**Indication:**

- To support sternum zone in post-operative period.

Size	
S	M   L   XL
75 cm	99 cm
99 cm	115 cm



BD-280-U DORSOFIX UNI.  
(DORSOLUMBAR CORSET)

**Functions and Features:**

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.

**Indication:**

- Can be used in compression fractures and herniated disks, for post-operative support, and in physiotherapy as a stable support.

Size	
S	M   L   XL
75 cm	99 cm
99 cm	115 cm



BD-170-U CLAVICULAFIX UNI.  
(CLAVICLE BANDAGE)

**Function and Features:**

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

**Indication:**

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.

ONE SIZE



BD-300-N ABDOMINAL BELT NEOPREN

**Functions and Features:**

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

**Indication:**

- Used to support thorax or abdomen area in post-operative period.

Size	
S	M   L   XL
75 cm	89 cm
99 cm	115 cm



BD-160-U VELPOFIX UNI.

**Functions and Features:**

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

**Indication:**

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.

ONE SIZE



BD-110-U (WITH WAIST BELT)

**Function and Features:**

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

**Indication:**

- Used when immobilization of arm and shoulder is required.

ONE SIZE



BD-115-U ARM SLING ADJUSTBALE

**Function and Features:**

- Made of air-permeable airtec spacer fabric.
- Ensures immobilization of arm
- Easy to wear and remove.

**Indication:**

- Used when immobilization of arm and shoulder is required.

ONE SIZE



BD-40600-U UNIFIX THUMB

**Function and Features:**

- Supports the thumb.
- Produced as a single body from material.
- Keeps the thumb at resting position.

**Indications:**

- Thumb collateral ligament injuries.
- Root arthrosis.
- Rheumatoid arthritis.



Size (Right-Left)



BD-40601-U UNIFIX MANU  
(UNIVERSAL WRIST SPLINT)

**Function and Features:**

- Ensures wrist support and stabilization.
- Made of 3 mm spacer material. Reduces the pain in the arthritis thanks to the heating feature of spacer
- An adjustable aluminum support on the volar surface of the wrist, and an extra bandage around the wrist are available.

**Indication:**

- For immobilization in wrist pains, for carpal tunnel syndrome, tendinitis, tenosynovitis, tenovaginitis, injuries, for protection in post-operative period, and after plaster treatment in fractures.



Size (Right-Left)



BD-180-U EPICONDYLITIS TENNIS ELBOW

**Functions and Features:**

- Produced as a standard size elbow band from soft fabric, and offers epicondylitis support.
- Epicondylitis support is ensured with a pad. This support
- Reduces the stress.

**Indication:**

- Medial and lateral epicondylitis, epicondylitis humeri (tennis elbow)



Size (Universal)



BD-460-U UNIVERSAL WRIST SPLINT

**Function and Features:**

- Ensures wrist support and stabilization.
- An adjustable aluminum support on the volar surface of the wrist, and an extra bandage around the wrist are available.

**Indication:**

- For immobilization in wrist pains, for carpal tunnel syndrome, tendinitis, tenosynovitis, tenovaginitis, injuries, for protection in post-operative period, and after plaster treatment in fractures.



Size (Right-Left)

S	M	L	XL
14 cm	20 cm	20 cm	24 cm



BD-40602-U UNIFIXMANU PLUS  
(WRIST SPLINT WITH THUMB SUPPORT)

**Functions and Features:**

- Ensures wrist and thumb support and stabilization.
- An adjustable aluminum baleen is available for volar support.
- Easy to apply on the patient thanks to its velcro tapes.
- Aluminum baleen is used in stabilization at neutral position. Produced as a single body.

**Indication:**

- Used in injuries, twists, rheumatic and traumatic arthritis, tendinitis, carpal tunnel syndrome pains, and for support in all plaster treatment cases requiring wrist rest.



Size (Right-Left)



BD-30603-U EPICONDYLITIS  
ELBOW BANDAGE WITH PAD

**Functions and Features:**

- The bandage that ensures epicondylitis support is made of air-permeable, three dimensional airtec spacer fabric.
- Reduces edema with compression.
- Reduces sweating, and avoids allergy.
- Epicondylitis support is ensured with an extra band and pad. This support reduces the stress on tendon, and ensures relief of pain.

**Indication:**

- Medial and lateral Epicondylitis (Tennis Elbow), Golfer's Elbow



PEDIATRICS



BA-511-P FOAM COLLAR CERVICAL ORTHOSIS

**Functions and Features:**

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

**INDICATION:**

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.

ONE SIZE



BA-110-P DELUX ARM SLING  
(WITH WAIST BELT)

**Function and Features:**

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

**Indication:**

- Used when immobilization of arm and shoulder is required.



Size (Right-Left)



BA-160-P VELPOFIX  
(VELPEAU BANDAGE)

**Functions and Features:**

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

**Indication:**

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.



Size (Right-Left)



BA-170-P CLAVICULAFIX  
(CLAVICLE BANDAGE)

**Function and Features:**

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

**Indication:**

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

XS	XXS	XXXS
50 cm	40 cm	40 cm
70 cm	50 cm	ve aşağı



BA-30601-P **MANUAIRTEC**  
(WRIST SPLINT)

**Functions and Features:**

- Ensures wrist support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Aluminum support that keeps the wrist at neutral position is available. Ensures wrist support and stabilization.
- Stability is increased with the extra tape on the wrist.

**Indication:**

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



**Size (Right-Left)**

XS	XXS	XXXS
13 cm	11 cm	11 cm
15 cm	13 cm	ve Aşağı



BA-30602-P **MANUAIRTEC PLUS**  
(WRIST SPLINT WITH THUMB)

**Function and Features:**

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Ensures wrist and thumb support and stabilization.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

**Indication:**

- MCP joint disorders, CMC joint degenerations
- De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



**Size (Right-Left)**

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-502-P **TRUSS FOR CHILDREN**

**Functions and Features:**

- Used for inguinal hernia, and in post-operative period.



**Size Universal**  
Right-Left



BA-503-P **TRUSS FOR CHILDREN**

**Functions and Features:**

- Used for inguinal hernia, and in post-operative period.



**Size Universal**



BA-411-P **HIP-X CRUSER HIP ABDUCTION SPLINT**

**Function and Features:**

- Positions the hip joint at flexion and abduction.
- Ensures femur head to stay in acetabulum.
- External part is made of polyethylene while the internal part is made of air-permeable spongy fabric.
- All materials can be washed. Easy-to-apply, light.
- Single orthoses would be enough during treatment, as the adjustment room is wide.

**Indication:**

- Coksit, hip dislocation.



Size			
S	M	L	
32 cm	37 cm	42 cm	
37 cm	42 cm	46 cm	



BA-416-P **PAVLIK BANDAGE**

**Functions and Features:**

- Positions both hip joints at flexion and abduction.
- Avoids adduction and extraction.
- This position and movement ensure spontaneous reduction of the hip, and acetabular development in displaced hip.
- Ensures femur head to stay in acetabulum during free movements of the child.
- Easy-to-wear and remove; made of velcro bands and cotton column.

**Indication:**

- Hip displacement, subluxation.



Size				
S	M	L	XL	
24 cm	34 cm	43 cm	52 cm	
34 cm	43 cm	52 cm	61 cm	



BA-419-P **FREJKA**

**Functions and Features:**

- Made of 12 mm plastozote material.
- Allows compulsory hip and leg movements of the baby, as it is made of elastic material.

**Indication:**

- Coksit, hip dislocation.



Size					
S	M	L	XL	XXL	
27 cm	32 cm	37 cm	42 cm	46 cm	
32 cm	37 cm	42 cm	46 cm	53 cm	



BA-409-P **PAVLIK HARNES TUBINGER TYPE**

**Functions and Features:**

- Apply for developmental dysplasia (or dislocation) of the hip (DDH) which is an abnormal development of the hip joint. The ball at the top of the thighbone is not stable within the hip socket and the ligaments of the hip joint may also be stretched and loose. This used to be called Congenital Dislocation of the Hip (CDH).

**Features:**

- 1 It has straps that are fastened around the baby's legs and held up by shoulder and chest straps. This holds the hips and knees up with the legs apart. This is the best position for the hip joint to be in. It allows contact between the thigh and hip bones and helps strengthen the muscles and ligaments of the hip while it is developing.
- 2 Soft TPR material to make baby feel comfortable.
- 3 Adjustable pearl strap fits most of the babies

Size: Universal



BA-553-P DENNIS BROWN ORTHESIS

**Functions and Features:**

- Controls dorsiflexion and abduction angle of the foot by assembling special PEV boot on an aluminum bar.
- Abduction angle is stable, while the external rotation is adjusted.
- Ankle dorsiflexion angle should be adjusted by the technician according to the instructions of the doctor.

**Indication:**

- Pes ekinovarus (PEV), Hip Dysplasia

**Size:**

- The product is combined according to shoe size and shoulder width of the patient.  
To order the product separately: BA 553 - Dennis Brown Splint \*BA 553-1 Dennis Brown Boot.



BA-400-P IMMOFIX  
(KNEE MOBILIZER)

**Functions and Features:**

- Ensures immobilization at 0 degree extraction of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of three different panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

**Indication:**

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.

Size	XS	XXS
	40 cm	30 cm



BA-30301-P GENUAIRTEC STABLE  
(KNEE PAD WITH HINGE)

**Functions and Features:**

- Ensures medial and lateral stabilization and support.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.

**Indication:**

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros

Size	XS	XXS	XXXS
	24 cm	21 cm	19 cm
	26 cm	23 cm	21 cm



BA-30201-P GENUFIX PLUS  
(LIGAMENT AND PATELLA SUPPORTED KNEE PAD)

**Function and Features:**

- Ensures patella and ligament stabilization and support.
- Made of air-permeable, skin-friendly airtec (spacer) fabric that avoids sweating.
- Silicon support reduces pain and edema thanks to its massage effect.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization.

**Indication:**

- Advanced twists on the knee, swelling,
- Post-operative and post-traumatic irritations,
- Patella subluxation, Chondromalacia, Patellar tendinitis

Size	XS	XXS	XXXS
	24 cm	21 cm	19 cm
	26 cm	23 cm	21 cm



BA-940-P GENU X PATELLAR BAND

**Function and Features:**

- Support is provided on the patella tendon with a silicone pad suitable for the anatomic shape.
- It is aimed to reduce pain with compression.
- Air permeable and moist regulator.

**Indication:**

- Patellar tendinitis
- Chondromalacia

ONE SIZE



BA-702-P CORK SUPPORT

**Function and Features:**

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



Size: (23-34)



BA-589-P FOAMCAST  
(ANKLE STABILIZATION ORTHOSIS)

**Functions and Features:**

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External plastics are shaped anatomically. Fits the ankle.
- Latex foam coveren with fabric takes the anatomical shape of the ankle.

**Indication:**

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, ad post-operative cases.

Size Universal



BA-583-P ROCKER PLASTER SLIPPER

**Functions and Features:**

- Makes it easy to walk on your plaster covered feet thanks to its special base (Rocker).
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.

ONE SIZE



BA-549-P **ARM PIT CRUTCH**

**Product Info:**

- Arm pit support is available. No burr likely to harm your arm pit is found on the plastic support.
- Handle support is available.
- Handle height can be adjusted.
- Height can be adjusted.



BA-548-P **MOBILE WALKER**

**Product Info:**

- An appropriate walker for elderly, patients having troubles with walking, patients experienced foot paralysis, and it can be used to rehabilitate the patients receiving physiotherapy treatment.
- It can be folded thanks to the mobile joints of the angled pipes, and ensures easy and safe walking by lifting two legs and pushing forward.
- Easy-to-use.
- Offers an easy experience thanks to the high quality, soft sponge knobs appropriate for the anatomical structure of the hand.
- Height can be adjusted according to the demands of the patients, thanks to the click-fit legs.
- Does not slide on wet surface thanks to the rubber bottom supports.
- Does not scratch or damage the gentle surfaces such as parquet.
- Metal parts can be wiped with a moist duster.
- Ensures easy and safe walking by lifting two legs and pushing forward.



BA-550-P **FOREARM CRUTCH**

**Product Info:**

- Arm support is made of rigid plastic.
- Hence, no burr likely to harm your arm is found on the plastic support.
- Made of light aluminum alloyed profile.
- Height can be adjusted.



SPORTS

BA-916 ELBOW SLEEVE SPORT WITH SILICON PAD

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- Silicone pads that compresses condyles reduce oedema and oedema related pain.

**Indications:**

- Lateral or Medial Epicondylitis
- Elbow Dislocation
- Post-operative or Post-traumatic inflammation



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-905 ELBOW SLEEVE LUX

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

**Indications:**

- Elbow Dislocation
- Post-operative or Posttraumatic inflammation
- Osteoarthritis, Arthritis
- Tenosynovitis



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-915 ELBOW SLEEVE SPORT

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.

**Indications:**

- Elbow Dislocation
- Post-operative or Post-traumatic inflammation



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-20303 EPITEX SPORT  
(ELASTIC KNITTED ELBOW SUPPORT)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.

**Indication:**

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm





BA-20602 WRIST SLEEVE BASIC

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- High wearing comfort.

**Indications:**

- Wrist soft tissue injuries
- Oedema and pain around the wrist
- Wrist support



**Size (Right-Left)**

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-907 WRIST SLEEVE LUX

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression. The inner surface is made of cotton outer surface polyamide thread.
- Air permeable and moist regulator.
- High wearing comfort.

**Indications:**

- Wrist soft tissue injuries
- Edema and pain around the wrist
- Wrist support



**Size (Right-Left)**

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-917 WRIST SLEEVE SPORT

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- High wearing comfort.

**Indications:**

- Wrist soft tissue injuries
- Oedema and pain around the wrist
- Wrist support



**Size (Right-Left)**

S	M	L	XL
14 cm	16 cm	18 cm	20 cm
16 cm	18 cm	20 cm	22 cm



BA-911 KNEE SLEEVE SPORT WITH SILICON PAD

**Functions and Features:**

- Support and stabilization of patella
- It is produced with circular knitting technique that provides balanced compression.
- Air permeable and moist regulator.
- Supports the circulatory system and stimulates muscle regeneration.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

**Indications:**

- Patella subluxation
- Chondromalacia
- Osteoarthritis, arthritis
- Post-traumatic and post-operative
- Patellar tendinitis



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-20105 **GENUFIX SPORT**  
(ELASTIC KNITTED KNEE PAD)

**Functions and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

**Indication:**

- To relief the pain on knees, and to reduce edema.



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-900 **KNEE SLEEVE LUX**

**Functions and Features:**

- It is produced with circular knitting technique that provides balanced compression and can stretch in both directions.
- Thanks to equal compression, it reduces oedema and oedema related pain.
- Air permeable and moist regulator.
- Supports the circulatory system and stimulates muscle regeneration.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

**Indications:**

- Oedema and pain in and around Patella



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-910 **KNEE SLEEVE SPORT**

**Functions and Features:**

- It is produced with circular knitting technique that provides balanced compression and can stretch in both directions.
- Thanks to equal compression, it reduces oedema and oedema related pain.
- Air permeable and moist regulator.
- Supports the circulatory system and stimulates muscle regeneration.

**Indications:**

- Oedema and pain in and around Patella



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-920 **THIGTS SUPPORT**

**Functions and Features:**

- Helps relieve pain from an injured or weak shin.
- Designed to conform naturally to the contours of the shin.
- Effective support for strained muscles and torn hamstrings.
- Hold the top of support using both hands, one on each side.
- Insert foot into support and pull up to the shin.

**Indications:**

- Muscle weakness



**Size**

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BA-921 CALF SUPPORT

**Functions and Features:**

- Aids pain relief due to varicose veins, post surgery, pregnancy and restricted circulation.
- Provides moderate compression and support.
- Recommended for weak muscles or leg veins by surgery or lack of exercise.
- Hold the top of support using both hands, one on each side.
- Insert foot into support and pull up to the calf.

**Indications:**

- Muscle weakness
- Relief muscle fatigue



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	



BA-922 ELBOW COMPRESSION SLEEVE

**Functions and Features:**

- Protects your elbows from shocks and bumps during contact sports, while providing constant compression and soothing warmth. Lightweight and comfortable fabric is designed to stretch and contour for use and multi-purpose protection.

**Size:**

- Small :20.3-27.9 cm.
- Medium 25.4-33.0 cm
- Large 30.4-38.1 cm.
- Xlarge 35.5-43.2 cm



BA-923 ELBOW COMPRESSION SLEEVES WITH WRIST

**Functions and Features:**

- Protects your elbows from shocks and bumps during contact sports, while providing constant compression and soothing warmth. Lightweight and comfortable fabric is designed to stretch and contour for use and multi-purpose protection.

**Size:**

- Small :20.3-27.9 cm.
- Medium 25.4-33.0 cm
- Large 30.4-38.1 cm.
- Xlarge 35.5-43.2 cm



BA-914 ANKLE SLEEVE SPORT WITH SILICON PAD

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- It reduces oedema and oedema related pain with compression.
- It can be used as a protection against sprains.
- The edema and hematoma around the malleol are reduced with the malleol pads that apply compression around the malleol.

**Indications:**

- Acute ankle sprains
- Soft tissue traumas
- Arthritis



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	



BA-20205 **MALLEOFIX SPORT**  
(ELASTIC WEAVED ANKLE PAD)

**Function and Features:**

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Reduces the pain and removes the edema thanks to the compression feature and massage effect.

**Indication:**

- Used in acute ankle twists, in case of any pain, soft tissue traumas, and arthritis cases.



**Size**

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-903 **ANKLE SLEEVE LUX**

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression. The inner surface is made of cotton outer surface polyamide thread.
- Air permeable and moist regulator.
- It reduces oedema and oedema related pain with compression.
- It supports healthy and correct walking.

**Indications:**

- Acute ankle sprains
- Soft tissue traumas
- Osteoarthritis, Arthritis



S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



BA-913 **ANKLE SLEEVE SPORT**

**Functions and Features:**

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- It reduces oedema and oedema related pain with compression.
- As a protection against sprains
- High wearing comfort.

**Indications:**

- Acute ankle sprains
- Soft tissue traumas
- Arthritis



S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



THE  
H  
E  
R  
E  
M  
A  
T

BA-316 ELBOWHEAT

**Functions and Features:**

- Keeps the elbow warm .
- Made of acrylic ,elastic fabric

**Indications:**

- Used in rheumatic sciatic,and elbow pains



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	



BA-317 ELBOWHEAT LONG

**Functions and Features:**

- Keeps the elbow warm .
- Made of acrylic ,elastic fabric

**Indications:**

- Used in rheumatic sciatic,and elbow pains



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	



BA-303 HEAT CORSET

**Functions and Features:**

- Ensures abdominal and waist support.
- Made of elastic fabric.
- Closed on the front part.
- Keeps warm due to heating feature.

**Indication:**

- Used to keep the affected areas hot in rheumatic and sciatic pains.



Size					
S	M	L	XL	XXL	
75 cm	83 cm	91 cm	99 cm	107 cm	
83 cm	91 cm	99 cm	107 cm	115 cm	



BA-303-V WITH VELCRO BAND

**Functions and Features:**

- Made of elastic fabric.
- Used to keep waist and abdominal area, and kidneys warm, in muscle spasms, and for rheumatic and sciatic pains.
- Can be opened and closed with the velcro on the front side.

**Indication:**

- Used to keep the affected areas hot in rheumatic and sciatic pains.



Size					
S	M	L	XL	XXL	
75 cm	83 cm	91 cm	99 cm	107 cm	
83 cm	91 cm	99 cm	107 cm	115 cm	



**BA-304 GENUWOOL**  
(WOLL KNEE PAD)

**Functions and Features:**

- Keeps the knees warm.
- Made of wool, elastic fabric.

**Indication:**

- Used in rheumatic, sciatic, and knee pains.



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	



**BA-312 GENUHEAT LONG**

**Functions and Features:**

- Keeps the leg warm .
- Made of acrylic ,elastic fabric

**Indications:**

- Used in rheumatic sciatic, and leg pains



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	



**BA-315 CALFHEAT**

**Functions and Features:**

- Keeps the calf warm .
- Made of acrylic ,elastic fabric

**Indications:**

- Used in rheumatic sciatic, and calf pains



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	



**BA-314 TIGHTSHEAT**

**Functions and Features:**

- Keeps the thigh warm .
- Made of acrylic ,elastic fabric

**Indications:**

- Used in rheumatic sciatic, and thighs pains



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	



PILOT W



**Product Info:**

- Made of visco elastic material.
- Pillow case is available. Case can be removed and washed.

BA-609  
VISCO NECK PILLOW



BA-610  
VISCO CONTOUR NECK PILLOW



BA-615  
VISCO ANATOMIC NECK PILLOW



BA-611  
VISCO WAIST PILLOW



BA-616  
VISCO ANATOMIC WAIST PILLOW



BA-612  
VISCO SEATING WHEEL



BA-613  
VISCO SEATING WHEEL U FORM



BA-619  
VISCO SEATING CUSHION COCCYX



BA-620  
VISCO SEATING CUSHION



BA-617  
VISCO CYLINDER PILLOW



BA-618  
VISCO CYLINDER PILLOW HALF



BA-621  
VISCO SUPPORT PILLOW



BA-614  
VISCO TRAVEL NECK PILLOW



BA-622  
VISCO VEHICLE NECK PILLOW



BA-612 - P  
POLIURETAN SEATING WHEEL



BA-555 TOILET RISER

**Functions and Features:**

- Dual function stability frame and raised toilet seat. A toilet safety frame is a great way to ensure that safety is maintained for vulnerable individuals in the bathroom. Elderly, disabled and mobility restricted individuals may need an extra level of support when using the bathroom in order to maintain independence and self-esteem

**Endication:**

- Suitable for assisted living and aged care
- 19.5 in. - 21 in. L x 21.5 in. W x 25 in. - 30 in. H
- Lightweight and portable, easy to relocate
- White powder-coated steel frame
- Height adjustable legs
- Non-slip feet
- Supports up to 250 lbs. 114 kg.



BA-852 TOILET SEAT BOOSTER FOR CLOSET

**Product Info:**

- BA-853 Toilet Seat Booster (With Legs)
- Toiler Seat Booster with Cover
- Toiler seat booster is first class, white, and plastic.
- Used for the patients having trouble in standing up and sitting down after hip, knee, and waist surgery.
- Can be used in all toilets.
- Boosts the sitting height of the toilet for 11 cm.
- The toilet seat booster has key apparatus, anti-bacterial, and does not hold any bacteria.
- The product has a cover, and this cover can open and close the seat.



BA-548 **MOBILE WALKER**

**Product Info:**

- An appropriate walker for elderly, patients having troubles with walking, patients experienced foot paralysis, and it can be used to rehabilitate the patients receiving physiotherapy treatment.
- It can be folded thanks to the mobile joints of the angled pipes, and ensures easy and safe walking by lifting two legs and pushing forward.
- Easy-to-use.
- Offers an easy experience thanks to the high quality, soft sponge knobs appropriate for the anatomical structure of the hand.
- Height can be adjusted according to the demands of the patients, thanks to the dick-fit legs.
- Does not slide on wet surface thanks to the rubber bottom supports.
- Does not scratch or damage the gentle surfaces such as parquet.
- Metal parts can be wiped with a moist duster.
- Ensures easy and safe walking by lifting two legs and pushing forward.



BA-549 **ARM PIT CRUTCH**

**Product Info:**

- Arm pit support is available. No burr likely to harm your arm pit is found on the plastic support.
- Handle support is available.
- Handle height can be adjusted.
- Height can be adjusted.



BA-550 **FOREARM CRUTCH**

**Product Info:**

- Arm support is made of rigid plastic.
- Hence, no burr likely to harm your arm is found on the plastic support.
- Made of light aluminum alloyed profile.
- Height can be adjusted.



BA-551 **TRIPOD WALKING STICK**

**Product Info:**

- Balance walking stick, tripod aluminum.
- Made of aluminum.
- Height can be adjusted.
- Handle is appropriate for the shape of the hand.
- Has three legs.
- Can be used in both sides (right and left).
- Ensures stable walking support.



**Biofix**<sup>®</sup>  
Orthopedic Soft Products

**Varifix**<sup>®</sup>  
Medical compression stockings

**Bodyflex**<sup>®</sup>  
Liposuction Garments

**FootFlex**<sup>®</sup>  
Foot Orthotic and orthopedic Footwear

 **Bioaktif**  
Support For Motion

Bioaktif Tıbbi Ürünler ve Teknik Tekstil San. Tic. Ltd. Şti.  
Haraççı Mah. Bayrak Cad. No: 26 Arnavutköy - İstanbul  
Tel: +90 212 567 90 34 Pbx [www.bioaktif.com](http://www.bioaktif.com)