



2024

CATALOG



Biofix[®]

Orthopedic Soft Products



Hello everyone,


Bioaktif was established in 1995, originally operated as prosthesis orthosis workshop. Our company began to produce orthopedic soft products in 1999 in order to reduce import-dependency in ready-to-use orthopedic products and to produce new ones.

In the meantime, in 2005 Bioteks narrow woven and textile limited company was established and began to produce orthopedic elastic fabrics as raw material is the most important factor in production of quality products and product development. In 2007, machinery for flatted knitted bangades was established, and the company became a pioneer in the sector.

Orthopedic soft products and orthopedic memory pillows under the brand of Biofix, orthopedic slippers under the brand of Footflex, after-surgery garments and after surgery body supports under the Bodyflex, medical compression stockings under the Varifix brand hit stores not only in Turkish market but also in the markets all around the world.

Our service is innovative and of high quality. We stick to our comitments in order to develop our company constantly and to produce world class products. Our always prior objective is to go beyond the expectations of our customers, to maximize the customer satisfaction and find the most appropriate solutions to patient problems.

Yours faithfully,
Mehmet BÜYÜKÇANGA





Contents

Cervical Zone Orthosis

3-6

Shoulder and Elbow Orthosis

7-17

Wrist Orthosis

18-26

Waist and Upper Body Orthosis

27-41

Hip Orthosis

42-43

Knee Orthosis

44-53

Ankle Orthosis

54-60

Foot Health

61-65

Orthopedic Pillows

66-68

Walking Aids

69-70

Biofix[®]

Orthopedic Soft Products

Biofix[®]
Orthopedic Soft Products



CERVICAL ZONE ORTHOSIS



BA 506

PHILADELPHIA COLLAR

Functions and Features:

- Ensures cervical zone support and stabilization.
- Anatomically appropriate for thorax, mandibula and occiput.
- Stabilizes cervical spines in accordance with the anatomical position.
- Consists of two parts (front and back).
- Has plastic enhancers in front and back.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcros on the sides.
- Reduced perspiration thanks to the air holes.

Indication:

- Cervical trauma
- Disk hernia
- Cervical stabilization
- Support and stabilization of cervical spine in pre and post-operative periods.



Size

XS/0	S/1	M/2	L/3	XL/4
28 cm	33 cm	38 cm	42 cm	46 cm
32 cm	37 cm	42 cm	46 cm	50 cm

BA 507

PHILADELPHIA COLLAR (Open Tracheostomy)

Functions and Features:

- Ensures cervical zone support and stabilization.
- Anatomically appropriate for thorax, mandibula and occiput.
- Stabilizes cervical spines in accordance with the anatomical position.
- Consists of two parts (front and back).
- Has plastic enhancers in front and back.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcros on the sides.
- Reduced perspiration thanks to the air holes.

Indication:

- Cervical trauma
- Disk hernia
- Cervical stabilization
- Support and stabilization of cervical spine in pre and post-operative periods.



Size

XS/0	S/1	M/2	L/3	XL/4
28 cm	33 cm	38 cm	42 cm	46 cm
32 cm	37 cm	42 cm	46 cm	50 cm



BA 508

NELSON PLASTOZOTE COLLAR

Functions and Features:

- Supports cervical zone.
- Limits the neck in the direction of flexion.
- Made of hypoallergenic plastazote material.
- Easy to wear and remove thanks to the velcro

Indication:

- Disk hernia
- Cervical muscle spasm.
- Cervical spine support.
- For cervical disc hernias, injuries, and protection in post-trauma period.



Size

XS/0	S/1	M/2	L/3	XL/4
28 cm	33 cm	38 cm	42 cm	46 cm
32 cm	37 cm	42 cm	46 cm	50 cm



BA 509

CERVICAL COLLAR WITH CHIN

Functions and Features:

- Ensures cervical zone support and immobilization.
- Covers the chin to limit flexion.
- Made of plastic.
- Covered with artificial leather.
- Has a foam pad appropriate for the anatomical structure of the chin.

Indication:

- In cervical disc hernia,
- Disk hernia
- Cervical muscle spasm, Cervical spine support
- Injuries, for recovery and protection in post-surgery period.



Size

XS/0	S/1	M/2	L/3	XL/4	XXL/5
28 cm	33 cm	38 cm	42 cm	46 cm	50 cm
32 cm	37 cm	42 cm	46 cm	50 cm	54 cm



BA 510

CERVICAL COLLAR WITHOUT CHIN

Functions and Features:

- Ensures cervical zone support and immobilization.
- Made of plastic, Covered with artificial leather.
- Adjustable neck height.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm,
- Cervical spine support.
- Adjusted according to the size of neck circle.



Size

XS/0	S/1	M/2	L/3	XL/4	XXL/5
28 cm	33 cm	38 cm	42 cm	46 cm	50 cm
32 cm	37 cm	42 cm	46 cm	50 cm	54 cm



BA 511

FOAM COLLAR

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size

XS/0	S/1	M/2	L/3	XL/4
28 cm	33 cm	38 cm	42 cm	46 cm
32 cm	37 cm	42 cm	46 cm	50 cm



BA 511 U

FOAM CERVICAL COLLAR ORTHOSIS

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size

6 cm
8 cm
10 cm
12 cm





BA 515

NELSON COLLAR (FABRIC COATED)

Functions and Features:

- Ensures cervical zone support and stabilization. Made of hypoallergenic plastazote material.
- Covered with cotton, non-sweating fabric.

Indication:

- Traumatic cases, Osteoarthritis, Cervical tendinitis.
- Muscle spasm, disk hernia, for protection in post-trauma period.



Size

XS/0	S/1	M/2	L/3	XL/4
28 cm	33 cm	38 cm	42 cm	46 cm
32 cm	37 cm	42 cm	46 cm	50 cm

FOAM COLLAR UNI.

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size (Universal)

BA 516 U



BA 511 P

FOAM COLLAR PEDIATRIC

Functions and Features:

- Supports cervical zone.
- Made of foam.
- Covered with cotton stockinette.
- Has extra plastic support.

INDICATION:

- Neck pains, disk hernia, cervical muscle spasm, cervical spine support.



Size (Universal)

MINERVA CTO

Functions and Features:

- Ensures cervical zone support and stabilization.
- Cervical supports that support chin and occipital zone are made of plastic, and the product has body supports that support chest and back. A connection was ensured between body and cervical supports using metal pieces.
- The product gets support from the shoulder through chin and shoulder parts, and ensures extraction and fraction to the cervical zone, if required.

Indication:

- Used for protective purposes after Osteochondrosis, Spondylose, Spondylarthrose, Tortikolis and other operations.



Size

S	M	L
34 cm	38 cm	42 cm
38 cm	42 cm	46 cm

BA 513





Biofix®
Orthopedic Soft Products

SHOULDER & ELBOW ORTHOSIS



BA 413 UNI

EPIROM COMFORT ROM ELBOW ORTHOSIS UNI.

Functions and Features:

- An elbow orthosis with adjustable angle, designed in accordance with the anatomical elbow shape.
- Inner surface that touches the skin is made of foam-backed fabric that reduces sweating.
- The mechanism used for angle adjustments is an easy-to-use pushing system.
- Orthosis controlled elbow ensures flexion and extension.
- Flexion (0-15-30-45-60-75-90-105) and extension (-15-0-15-30-45-60-75-90)
- Angle can be adjusted. Orthosis length can be adjusted according to the length of arm.
- Aluminum support that avoids the wrist to bend is available and length can be adjusted.
- Shoulder belt offers comfort in use.

Indication:

- Support for rehabilitation in elbow contracture.
- Protection and support for rehabilitation in elbow fracture.
- Limited and controlled movement in the elbows.
- Conservative treatment after elbow dislocation, in post-op period for ligament instability.
- Post-op rehabilitation after elbow endoprosthesis.



Size (Universal Right and Left)

BA 412

HUMERUS SERRIANTO ORTHOSIS UNI.

Functions and Features:

- Made of 2 mm polyethylene, and covered with foam fabric on the inside.
- This fabric can be removed and washed.
- Ensures full harmony by surrounding the humerus with integrated band system.

Indication:

- With the aim of the protection of Humerus proximal fractures after operation, plaster treatment and stabilization.



Size (Universal)

S/1	M/2
19 cm	26 cm
25 cm	30 cm



BA 20505

SHOULDERFIX PLUS

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Elastic tapes on the shoulder.
- Strengthens the muscles by compressing the tissues with its elastic structure, avoids swelling.
- Limits rotation and spreading movements of the shoulders.
- Avoids humerus head to distance from shoulder joint, and limits the injuring movements.

Indication:

- Shoulder bone fractures • Shoulder arthritics.
- Rehabilitation after shoulder implant prosthesis treatment
- Post-operative and post-traumatic irritations.



Size

S	M	L	XL	XXL
22 cm	24 cm	27 cm	30 cm	33 cm
24 cm	27 cm	30 cm	33 cm	37 cm



BA 110

DELUX ARM SLING (With Waist Belt)

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.



Size	XS	S	M	L	XL
	25 cm	30 cm	35 cm	40 cm	45 cm
	30 cm	35 cm	40 cm	45 cm	50 cm



BD 110 U

ARM SLING (With Waist Belt) UNI.

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.



Size (Universal)



BA 110 P UNI

DELUX ARM SLING (With Waist Belt) PEDIATRIC

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.



Size (Universal)





BA 110 A


DELUX ARM SLING (With Waist Belt)

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Spongy fabric is used in the inner surface to reduce sweating.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder and waist slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Can be used in both arms.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.

Size	XS	S	M	L	XL
	25 cm 30 cm	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm	45 cm 50 cm


LUX ARM SLING (With Waist Belt)

Function and Features:

- Made of air-permeable airtex spacer fabric.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Easy to wear and remove.
- Shoulder and waist belts do not stretch, and do not cause discomfort for the patient.

Indication:

- Used when immobilization of arm and shoulder is required.

Size	XS	S	M	L	XL
	25 cm 30 cm	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm	45 cm 50 cm

BA 115



BD 115 U

ARM SLING ADJUSTABLE

Function and Features:

- Made of 4 mm, spongy, air-permeable, foam-backed fabric.
- Designed in accordance with the anatomical structure of the arm.
- Shoulder slings are made of soft fabric that does not stretch.
- Ensures immobilization of arm and shoulder.
- Easy to wear and remove.

Indication:

- Used when immobilization of arm and shoulder is required.

	Size (Universal)
---	------------------

BA 116


ECO ARM SLING (With Waist Belt)

Function and Features:

- Made of air-permeable airtec spacer fabric.
- Ensures immobilization of arm and shoulder.
- Shoulder movements are limited with the waist belt.
- Designed to be used in both arms.

Indication:

- Used when immobilization of arm and shoulder is required.

Size	XS	S	M	L	XL
	25 cm 30 cm	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm	45 cm 50 cm



BA 140


ACTIVE SLING 30° (Arm Sling with Pad)

Function and Features:

- Used in immobilization and positioning of shoulder.
- Made of 6 mm, spongy, air-permeable, foam-backed fabric.
- The pad that abducts the shoulder at 30° is made of 25x12 foam.
- Designed to ensure comfort in long-term treatment of arm, and to offer armpit air circulation.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

Indication:

- Rotator cuff, capsular rupture, and bankart repair. Used in shoulder joint dislocation and soft tissue injury.

Size	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm



BA 140 A


ACTIVE SLING 30° (Arm Sling with Pad)

Function and Features:

- Used in immobilization and positioning of shoulder.
- Made of 6 mm, spongy, air-permeable, foam-backed fabric.
- The pad that abducts the shoulder at 30° is made of 25x12 foam.
- Designed to ensure comfort in long-term treatment of arm, and to offer armpit air circulation.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

Indication:

- Rotator cuff, capsular rupture, and bankart repair. Used in shoulder joint dislocation and soft tissue injury.

Size	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm





BA 150


ACTIVE SLING 45°-60° (Arm Sling with Pad)

Functions and Features:

- Used in immobilization and positioning of shoulder.
- Six mm inner surface is made of foam-backed fabric that reduces sweating.
- Ensured 45° or 60° abduction angle on shoulder.
- Designed to ensure comfort in long-term treatment of arm.
- Ensures effective healing, enabling positioning of the arm in neutral position.
- Helps in avoiding internal rotation contractions in post-operative period.

Indication:

- In shoulder prosthesis operations, rotator cuff, capsular rupture, and bankart repairs. Used in shoulder joint dislocation and soft tissue injury.

Size	S	M	L
	30 cm 35 cm	35 cm 40 cm	40 cm 45 cm


VELPOFIX (Velpeau Bandage)

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.

Size	XS	S	M	L	XL
	60 cm 70 cm	70 cm 80 cm	80 cm 90 cm	90 cm 100 cm	100 cm 110 cm

BA 160



BD 160 U



VELPOFIX (Velpeau Bandage) UNI.

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.

Size (Universal)


BD 160 E

VELPOFIX (Velpeau Bandage) ECO

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.



Size (Universal)



BA 160 P UNI

VELPOFIX (Velpeau Bandage) PEDIATRIC

Functions and Features:

- Ensures immobilization of shoulder and clavicle.
- Made of 4 mm, air-permeable, foam-backed fabric.
- Reduces sweating with the air-permeable spongy fabric.
- A tape that fixes humerus to the body is available.
- Offers limited exercise opportunity when the tape is removed at the advanced phase of the treatment.

Indication:

- Clavicle fractures. Shoulder joint dislocation, soft-tissue injuries, stabilization of shoulder and arms.



Size (Universal)



BA 170

CLAVICULAFIX (Clavicle Bandage)

Function and Features:

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

XS	S	M	L	XL	XXL
60 cm	70 cm	80 cm	90 cm	100 cm	110 cm
70 cm	80 cm	90 cm	100 cm	110 cm	120 cm





BD 170 U

CLAVICULAFIX (Clavicle Bandage) UNI.

Function and Features:

- Made of 8 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size (Universal)

CLAVICULAFIX (CLAVICLE BANDAGE) PEDIATRIC

Function and Features:

- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

XS	XXS	XXXS
50 cm	40 cm	< 40 cm
70 cm	50 cm	

BA 170 P



BA 171

CLAVICULAFIX PLUS (Clavicle Bandage)

Functions and Features:

- Designed in accordance with anatomic structure.
- Made of 12 mm, air permeable, foam-backed fabric.
- Cross bands under the armpit gives hypertension on the shoulders, and ensures the bone to heal accurately.
- Non-slip fabric used and easy-to-wear feature ensure great comfort even while laying on back.

Indication:

- Used in clavicle fractures, clavicle and sternoclavicular dislocations.



Size

XS	S	M	L	XL	XXL
60 cm	70 cm	80 cm	90 cm	100 cm	110 cm
70 cm	80 cm	90 cm	100 cm	110 cm	120 cm

EPIFIX ACTIVE (Epicondylitis Elbow Band)

BA 185

Functions and Features:

- Produced as a standard size elbow band from velcro-compatible fabric, and offers epicondylitis support.
- Epicondylitis support is anatomically ensured with a silicon pad.
- The support reduces the stress on tendon, and ensures relief of pain.

Indication:

- Medial and lateral epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



BA 186

EPITEC COMFORT (Epicondylitis Elbow Band)

Functions and Features:

- Produced as a standard size elbow band from velcro-compatible fabric, and offers epicondylitis support.
- Epicondylitis support is anatomically ensured with a silicon pad.
- The support reduces the stress on tendon, and ensures relief of pain.

Indication:

- Medial and lateral epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



BD 180 U

EPICONDYLITIS ELBOW BAND UNI.

Function and Features:

- Produced as a standard size elbow band from soft fabric, and offers epicondylitis support.
- Epicondylitis support is ensured with a pad.
- This support
- Reduces the stress.

Indication:

- Medial and lateral epicondylitis, epicondylitis humeri (tennis elbow)



Size (Universal)



BA 20301



EPIFIX PLUS (Silicon Backed Elbow Support With Extra Band)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that apply pressure on condyls relieve the pain.

Indication:

- Lateral and Medial Epicondylitis (Tennis Elbow) Elbow dislocation. Chronic post-operative and post-traumatic inflammation.



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	

BA 20302-3

EPIFIX BASIC (Elastic Knitted Elbow Support)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Appropriate for use in sport.

Indication:

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	



BA 20302-4



EPIFIX BASIC (Elastic Knitted Elbow Support)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Appropriate for use in sport.

Indication:

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	

BA 20302-5

EPIFIX BASIC (Elastic Knitted Elbow Support)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Appropriate for use in sport.

Indication:

- Elbow dislocation, chronic post-operative and post-traumatic inflammation.



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	



BA 905

EPITEX LUX ELBOW SLEEVE

Functions and Features:

- It is produced with circular knitting technique providing balanced compression.
- Air permeable and moist regulator.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

Indications:

- Elbow Dislocation
- Post-operative or Posttraumatic inflammation
- Osteoarthritis, Arthritis
- Tenosynovitis



Size					
S	M	L	XL	XXL	
19 cm	21 cm	23 cm	25 cm	27 cm	
21 cm	23 cm	25 cm	27 cm	29 cm	



BD 30603 U

EPICONDYLITIS ELBOW BANDAGE WITH PAD UNI.

Functions and Features:

- The bandage that ensures epicondylitis support is made of air-permeable, three dimensional airtex spacer fabric.
- Reduces edema with compression.
- Reduces sweating, and avoids allergy.
- Epicondylitis support is ensured with an extra band and pad. This support reduces the stress on tendon, and ensures relief of pain.

Indication:

- Medial and lateral Epicondylitis (Tennis Elbow), Golfer's Elbow



Size (Universal)



WRIST ORTHOSIS

Biofix®
Orthopedic Soft Products

BA 20601

MANU COMFORT (Knitted Wrist Splint)

Functions and Features:

- Produced with an advanced technology using smart knitting-waving comfort polyamide threads.
- Air-permeable, breathable. Does not hold moisture, transfers the moisture to outside.
- Hygienic, non-allergic, anti-bacterial.
- Shaped in accordance with the anatomical structure, fits the hand and wrist.
- Offers high comfort and harmony. A formable thermoplastic splint is available.
- The silicon pad on the volar side of the hand ensures support, and relieves the pain.

Indication:

- Tenosynovitis
- Carpal instability
- CTS
- After plaster treatment of wrist fractures
- Rheumatoid arthritis
- Ganglion pain
- Soft tissue traumas on the wrist



Size (Right-Left)

S	M	L	XL	XXL
14 cm 16 cm	16 cm 18 cm	18 cm 20 cm	20 cm 22 cm	22 cm 24 cm



BA 907

WRIST SLEEVE LUX

Functions and Features:

- It is produced with circular knitting technique providing balanced compression. The inner surface is made of cotton outer surface polyamide thread.
- Air permeable and moist regulator.
- High wearing comfort.

Indication:

- Wrist soft tissue injuries
- Edema and pain around the wrist
- Wrist support



Size

S	M	L	XL	XXL
14 cm 16 cm	16 cm 18 cm	18 cm 20 cm	20 cm 22 cm	22 cm 24 cm



BA 460

MANU LACE (Wrist Splint With Lace)

Functions and Features:

- Supports sensorimotor movements, and ensures muscle stability of the wrist.
- Made of air-permeable, good-for-skin, elastic textile fabric.
- Ensures comfortable stabilization thanks to the adjustable laces that wrap the wrist.
- Limits the hand-wrist movement. Fully adjusts with the anatomical structure of the wrist.
- Keeps your wrist at neutral position, relieves you pain, helps through rehabilitation process

Indication:

- Acute and chronic irritation (e.g. tendovaginitis, tendomyopathy, minor joint effusion, arthralgia).
- Post-operative and post-trauma irritation.
- Wrist Injuries.
- Rheumatoid Arthritis.
- Carpal Tunnel Syndrome
- Degenerative Wrist Disorders (e.g. arthritis)
- Sense of weakness.



Size (Right-Left)

S-M / 1	L-XL / 2
14 cm 18 cm	18 cm 22 cm





BD 460 U

MANU LACE (Wrist Splint With Lace) UNI.

Functions and Features:

- Supports sensorimotor movements, and ensures muscle stability of the wrist.
- Made of air-permeable, good-for-skin, elastic textile fabric.
- Ensures comfortable stabilization thanks to the adjustable laces that wrap the wrist.
- Limits the hand-wrist movement. Fully adjusts with the anatomical structure of the wrist.
- Keeps your wrist at neutral position, relieves you pain, helps through rehabilitation process

Indication:

- Acute and chronic irritation (e.g. tendovaginitis, tendomyopathy, minor joint effusion, arthralgia).
- Post-operative and post-trauma irritation.
- Wrist Injuries.
- Rheumatoid Arthritis.
- Carpal Tunnel Syndrome
- Degenerative Wrist Disorders (e.g. arthritis)
- Sense of weakness.



Size (Universal)

S/1	L/2
14 cm	18 cm
18 cm	22 cm

MANUFIX (Elastic Wrist Splint)

Function and Features:

- Ensures wrist support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Has anatomically shaped aluminum support on the volar surface of the wrist.
- Easy to wear and remove thanks to the velcros.

Indication:

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

S	M	L	XL	XXL
14 cm	16 cm	18 cm	20 cm	22 cm
16 cm	18 cm	20 cm	22 cm	24 cm

BA 440



BA 30601



MANUAIRTEC (Wrist Splint)

Functions and Features:

- Ensures wrist support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Aluminum support that keeps the wrist at neutral position is available. Ensures wrist support and stabilization.
- Stability is increased with the extra tape on the wrist.

Indication:

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

S	M	L	XL	XXL
14 cm	16 cm	18 cm	20 cm	22 cm
16 cm	18 cm	20 cm	22 cm	24 cm

BA 30601 P

MANUAIRTEC (Wrist Splint) PEDIATRIC

Functions and Features:

- Ensures wrist support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Aluminum support that keeps the wrist at neutral position is available. Ensures wrist support and stabilization.
- Stability is increased with the extra tape on the wrist.

Indication:

- Wrist soft tissue traumas, carpal tunnel syndrome, instabilities, wrist fractures, after plaster treatment, rheumatoid arthritis, ganglion pain.



Size (Right-Left)

XS	XXS	XXXS
13 cm	11 cm	<11 cm
15 cm	13 cm	



BD 40601 U

UNIFIX MANU (Universal Wrist Splint) UNI.

Function and Features:

- Ensures wrist support and stabilization.
- Made of 3 mm airtec spacer fabric. Reduces the pain in the arthritis thanks to the heating feature of neoprene.
- An adjustable aluminum support on the volar surface of the wrist, and an extra bandage around the wrist are available.

Indication:

- For immobilization in wrist pains, for carpal tunnel syndrome, tendinitis, tenosynovitis, tenovaginitis, injuries, for protection in post-operative period, and after plaster treatment in fractures.



Size (Universal Right-Left)



BA 932

MANUFIX LUX WRIST SUPPORT

Function and Features:

- Air permeable and moist regulator.
- It covers the wrist, stabilizes and provides thumb movement.
- There is immobilization with ligament supports on the inner and outer sides of the palm.

Indication:

- Carpal Tunnel Syndrome
- Painful wrist conditions
- Wrist strain



Size (Universal)

S-M / 1	L-XL / 2
14 cm	18 cm
18 cm	22 cm



BA 450

MANUFIX PLUS (Wrist Splint With Thumb)

Function and Features:

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic and air-permeable textile fabric that does not harm the skin.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

Indication:

- MCP joint disorders, CMC joint degenerations, De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

S	M	L	XL	XXL
14 cm	16 cm	18 cm	20 cm	22 cm
16 cm	18 cm	20 cm	22 cm	24 cm



BA 30602

MANUAIRTEC PLUS (Wrist Splint With Thumb)

Function and Features:

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Ensures wrist and thumb support and stabilization.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

Indication:

- MCP joint disorders, CMC joint degenerations
- De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

S	M	L	XL	XXL
14 cm	16 cm	18 cm	20 cm	22 cm
16 cm	18 cm	20 cm	22 cm	24 cm



BA 30602 P

MANUAIRTEC PLUS (Wrist Splint With Thumb)

Function and Features:

- Ensures wrist and thumb support and stabilization.
- Made of non-allergic, air-permeable airtec spacer fabric that avoids sweating.
- Ensures wrist and thumb support and stabilization.
- Aluminum support that keeps the wrist at neutral position is available.
- Required abduction angle can be adjusted using the aluminum support in the thumb.

Indication:

- MCP joint disorders, CMC joint degenerations
- De Quervain, and other cases where carpal tunnel and De Quervain are observed together.



Size (Right-Left)

XS	XXS	XXXS
13 cm	11 cm	<11 cm
15 cm	13 cm	



BD 40602 U

UNIFIXMANU PLUS (Wrist Splint With Thumb Support) UNI.

Functions and Features:

- Ensures wrist and thumb support and stabilization.
 - Made of neoprene fabric - Reduces the pain in the arthritis thanks to the heating feature of neoprene.
 - An adjustable aluminum baleen is available for volar support.
 - Easy to apply on the patient thanks to its velcro tapes.
- Aluminum baleen is used in stabilization at neutral position. Produced as a single body.

Indication:

- Used in injuries, twists, rheumatic and traumatic arthritis, tendinitis, carpal tunnel syndrome pains, and for support in all plaster treatment cases requiring wrist rest.



Size (Universal Right-Left)



BA 931

DEQUFIX LUX THUMB SUPPORT

Function and Features:

- Air permeable and moist regulator.
- Designed to cover wrist and thumb.
- Support and stabilization of thumb joints

Indication:

- Tendinitis
- Skier finger
- De Quervain Syndrome



Size (Universal)

S-M / 1	L-XL / 2
14 cm	18 cm
18 cm	22 cm



BA 451 U

DEQUFIX (Thumb Fixing Splint) UNI.

Function and Features:

- Ensures support and stabilization of thumb joints.
- Made of elastic textile fabric.
- Has a formable aluminum support.
- Fixes CMC and MC joints.

Indication:

- In treatment of chronic, post-operative, and post-traumatic tissue sensitivity.
- For stabilization in ulnar cruciate ligaments.
- Carpal Tunnel Syndrome.
- Rheumatoid arthritis.
- Root arthrosis.



Size (Universal Right-Left)



BA 453

THUMB SPLINT

Functions and Features:

- Limits IBP / CMC joint movements.
- Made of soft and breathable fabric.
- Easy-to-use and adjustable thanks to elastic tapes.
- De Quervain tendinitis
- Gamekeeper thumb
- Collateral ligament immobilization



Size (Universal)



BD 40600 U



UNIFIX

Function and Features:

- Produced as a single body from neoprene material.
- Keeps the thumb at resting position.

Indications:

- Thumb collateral ligament injuries.
- Root arthrosis.
- Rheumatoid arthritis.



Size (Universal)

BA 455

BOXER FRACTURE SPLINT

Functions and Features:

- The Biofix Boxer Splint is an easy to adjust splint that positions and protects the 4th and 5th metacarpals in the treatment of "boxer fractures", and offers superior comfort and support to enhance compliance.
- The malleable palmar and dorsal aluminum stays are easily bent to immobilize the patient's hand.
- The products design protects and supports the fourth and fifth metacarpals without restricting flexion.

Indications:

- Boxer fracture



Size (Universal Right-Left)



BA 183



PATIENTS HANDS AND FEET STRAP BELT

Functions and Features:

- Made of soft fabric with padding, breathable and comfortable to wear
- Size: adjustable, suitable for wrist or ankle circumference:
8.6" ~ 11.4"/22 ~ 29cm
- Velcro close, easy and convenient to use
- Use this medical restraint strap to prevent patients from self-injury to ensure safety and smooth care or treatment
- Unconscious bedridden patients, mental disorder people, cognitive impairment people etc. applicable to.

Size (Universal)

BA 430

MANUSTABIL SPASTIC WRIST SPLINT

Function and Features:

- Ensures neutral position in hands and wrists.
- Made of 3 mm polyethylene using vacuum method.
- Internal surface is covered with 3 mm plastozote.
- Durable, can be used in long-term treatments, easy-to-wear and remove.

Indication:

- For protections after fractures. For paralyzed patients, in wrist problems, and in cases where hand is required to be rested.



Size (Right-Left)

S	M	L
14 cm	16 cm	18 cm
16 cm	18 cm	20 cm



BA 431

MANU ANTI-SPASTISITE (Anti-Spasticity Wrist Splint Separator)

Functions and Features:

- Keeps the hand, wrist, and fingers at neutral position.
- Made of 3 mm polyethylene using vacuum method.
- Internal surface is covered with 3 mm plastazote.
- Durable, can be used in long-term treatments, easy-to-wear and remove.
- Ideal for patients with spasticity due to its channeled finger design and resistance to stress.
- Ensures full harmonization and control with the special-designed tape system.

Indication:

- For patients with spasticity, protective after fractures. For paralyzed patients, in wrist problems, and in cases where hand is required to be rested.



Size (Right-Left)

S	M	L
14 cm 16 cm	16 cm 18 cm	18 cm 20 cm



BA 101



FROG FINGER SPLINT

Function and Features:

- Used for fixing the fingers at required position.
- Internal surface is covered with soft plastazote.
- Can be formed if required.



Size

S	M	L
---	---	---

BASEBALL FINGER SPLINT

Function and Features:

- Used to keep the finger at extension position.
- Used in finger injuries, tendon splitting cases, and low joint cases.
- Size and numbers varies.



Size

S	M	L
---	---	---

BA 102





BA 104

MALLET FINGER SPLINT

Function and Features:

- Mallet finger splint is a finger orthosis made of special material and is used to keep the DIP joint at extraction position.
- Finger is free from the point of PIP joint.
- Extensor tendon injuries.
- Finger Splint is used in fractures.
- Used in injuries and twists.

Indication:

- Extensor tendon treatment, mallet finger deformity. Finger bone fractures. Finger injuries. Finger nerve damages.



Size

1 2 3 4 5 5.5 6

DYNAMIC FINGER SPLINT

Function and Features:

- Made of stainless steel spring covered with eloxal.
- Points that touch the finger are made of special plastic with rounded edges.
- Produced in two sizes.
- The center point of the spring is attached using a screw, while the other edges are assembled with stainless rivet.
- Used in tendon cuts and nerve injuries.

Indication:

- Improvement of muscle immobility,
- Improvement of joint immobility,
- Ulnar nerve damages,
- Rehabilitation in post-operative period



Size

S M L

BA 118





WAIST & UPPER BODY ORTHOSIS

Biofix®

Orthopedic Soft Products



BA 20501

LUMBOSOFT (Knitted Lumbosacral Corset)

Functions and Features:

- Supports the lumbar region.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Compatible with the anatomy, and comfortable. Ensures anatomically appropriate support.
- Easy to wear and remove thanks to the velcros on the front side.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 240

LUMBOLIGHT PLUS (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA 245

LUMBOLIGHT EVO

Functions and Features:

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbar.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 26 cm, Front height: 16 cm.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome,



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BD 240 U

LUMBOLIGHT UNI. (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm



BD 240 US

**LUMBOLIGHT UNI.
(Lumbosacral Corset With An Extra Belt)**

Function and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove. Back height is 26 cm.

Indication:

- Acute back pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm



BA 20503

LUMBOSACRAL CORSET

Function and Features:

- Supports the lumbar region.
- A triangle silicon pad-support appropriate for the anatomy of the body is placed on the back side of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Compatible with the anatomy, and comfortable. Silicon support appropriate for the anatomical structure is available.
- Easy to wear and remove thanks to the velcros on the front side.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm





BD 243 U

WAIST SUPPORT BLACK (17 cm)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, nylon fabric that avoids sweating.
- Four plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove.

Indication:

- Acute waist pain, degenerative disc, paravertebral muscle spasm.

Size	S-M / 1		L-XL / 2	
	75 cm	99 cm	99 cm	115 cm

LUMBOLIGHT UNI ECO. (Lumbosacral Corset With An Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, polyester fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Easy to wear and remove.

Indication:

- Acute waist pain, degenerative disc, paravertebral muscle spasm, lumbago - sciatic.

Size	S	M	L	XL	XXL
	75 cm 83 cm	83 cm 91 cm	91 cm 99 cm	99 cm 107 cm	107 cm 115 cm

BA 240 E



BA 20502

LUMBOSOFT (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable, three-dimensional textile fabric that avoids moisture.
- Easy-to-wear thanks to hand bands on the front right and left sides of the corset.
- Back height: 30 cm, Front height: 16 cm.
- Offers more rigid support and stabilization with the extra bands.

Indication:

- Lumbar spine instability, lumbar stress syndromes, post-operative treatment after disc surgery.

Size	S	M	L	XL	XXL
	75 cm 83 cm	83 cm 91 cm	91 cm 99 cm	99 cm 107 cm	107 cm 115 cm

BA 246

LUMBOLIGHT EVO PLUS

Functions and Features:

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbar.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 26 cm, Front height: 16 cm.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, disorders, post-discectomy syndrome, ligamentous



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA 230

LUMBOLIGHT (Lumbosacral Corset 26 cm)

Functions and Features:

- Supports and stabilizes the lumbar region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two adjustable steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset offers full support.

Indication:

- Can be used in lumbago, sciatic, and herniated disk problems, for post-operative support, and in physiotherapy as a stable support.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA 260

LUMBOFIX PLUS (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support, and ensure full harmonization with the body.

Indication:

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago - sciatic, spondylolisthesis, post-operative immobilization.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm





BD 260 U

LUMBOFIX UNI. (Lumbosacral Corset with an Extra Belt)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Produced in accordance with the anatomy of the body.
- Back height: 30 cm, Front height: 18 cm.

Indication:

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago sciatic, spondylolisthesis, post-operative immobilization.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm

LUMBOFIX UNI. (Lumbosacral Corset)

Functions and Features:

- Supports and stabilizes the lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.
- Produced in accordance with the anatomy of the body.
- Back height: 30 cm, Front height: 18 cm.

Indication:

- Acute and chronic waist pain, degenerative disc, hypermobile and degenerative changes, lumbago sciatic, spondylolisthesis, post-operative immobilization.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm

BD 260 US



BA 265

LUMBOFIX EVO

Functions and Features:

- Ensures strong lumbosacral support and stabilization.
- Four flexible steel baleens appropriate for the anatomy are available at the back of the corset.
- Can be used in both male and female patients thanks to its greater contour differences between hip and waist.
- Offers effective support in especially elder and overweight patients, as it compresses the belly, and supports the lumbosacral.
- Relieves the pain by reducing the muscular stress.
- Made of air-permeable fabric that avoids moisture.
- Offers extra comfort as it is produced with a bonding technique, avoiding any stitch.
- Easy to wear and remove thanks to hand bands on the front panel of the corset.
- Back height: 32 cm, Front height: 20 cm.

Indication:

- Lumbago, disc hernia, sciatic, risk prevention in sports and work activities, ligamentous disorders, post-discectomy syndrome, post-operative support and stabilization.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 933

STRONG LUMBAR CORSET

Functions and Features:

- Ensures stabilization and support of the lower back.

Indications:

- During the handling of loads, labor protection in continuous overexertion or lumbar



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BA 280

DORSOFIX PLUS (Dorsolumbar Corset)

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.

Indication:

- Can be used in compression fractures and herniated disks, for post-operative support, and in physiotherapy as a stable support.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



BD 280 U

DORSOFIX UNI. (Dorsolumbar Corset)

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic, inner surface cotton fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- The extra belt on the waist strengthens the support.

Indication:

- Can be used in compression fractures and herniated disks, for post-operative support, and in physiotherapy as a stable support



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm





BA 270

DORSOFIX (Dorsolumbar Corset)

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support.

Indication:

- Compression fractures and herniated disks, post-operative support,
- In physiotherapy, as a stable support.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

FABRIC TISO

Functions and Features:

- Made of rigid and air-permeable fabric.
- Limits flexion-extension of the body.
- Avoids lateral bending and rotation.
- Polyethylene rigid supports are available at the back and front sides.
- Limits extension through thorax pad.

Indication:

- To mobilize the patient after surgery,
- In compression fractures, not requiring surgery.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 287



BA 286

JEWEEET / HYPEREXTENSION CORSET

Functions and Features:

- Limits the flexion of the body based on 3 point principle, and ensures hyperextension.
- Stabilizes the spine against lateral flexion and rotation movements.
- A mobile hip bar that ensure greater comfort while seated.
- Four different sizes at the side bars, which can be adjusted.
- Easy use with the B-fast lock.
- Adjustable chest bars.

Indication:

- Stable compression fractures in lumbar and inferior thoracic vertebra.
- Immobilization after intervertebral disc surgery.
- Axial stabilization in thoracic and osteoporotic compression fractures of lumbar spine.
- Temporary orthotic care after open surgeries on non-stable vertebral body fracture.
- Permanent orthotic care for spinal tumors and metastases which are not appropriate for surgery.
- Post-operative orthotic care after vertebral body tumor and metastases surgeries.
- Care after spinal decompression surgery with or without internal fixation.



Size

S	M	L
75 cm	90 cm	105 cm
90 cm	105 cm	120 cm

POSTURFIX

Functions and Features:

- Ensures support to fix the posture.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Two plastic baleens are available at the back side.
- Bands on the shoulder ensure postural support, and avoid kyphoses.

Indication:

- Used for weak and lazy thoracic zone muscles, and ad stimulant and straightener in posture disorders.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 310



BD 310 U



POSTURFIX UNI.

Functions and Features:

- Ensures support to fix the posture.
- Made of air-permeable, elastic, cotton fabric that avoids sweating.
- Two plastic baleens are available at the back side.
- Bands on the shoulder ensure postural support, and avoid kyphoses.

Indication:

- Used for weak and lazy thoracic zone muscles, and ad stimulant and straightener in posture disorders.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm

POSTURVEST

- Made from high quality material, which is breathable, durable and comfortable to wear.
- Fashionable Y-shape cross belt design on the back will get better effect, it can help to straight your back and stand upright. Ergonomic shoulder design will protect your shoulder from damaging caused by friction.
- The hook and loop on the waist belt is convenient to adjust the size and easy to wear and take off.
- Help to balance stress on the shoulder, back, waist, abdomen, relieve muscle fatigue and back pain, and protect the waist, back and spine to promote your image and grow healthily.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 313 S



BA 313 T



BA 314 T

POSTURVEST ECO

Functions and Features:

- Supports and stabilizes the dorsal and lumbosacral region.
- Made of air-permeable, elastic fabric that avoids sweating.
- Shoulder bands avoid rotation and kyphoses.
- A soft pad is placed at the back of the corset to support lordosis.
- Two steel baleens and two plastic baleens are available at the back side.
- Cross bands at the back of the corset increase the support.

Indication:

- Compression fractures and herniated disks, post-operative support,
- In physiotherapy, as a stable support.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 314 S



PREGNANT CORSET

Functions and Features:

- Helps to avoid waist pain during pregnancy.
- Designed to support the abdominal area, and to lift it upwards.
- Ensures maximum comforts thanks to the ergonomic design of the corset.
- Support can be increased or decreased with the elastic bands on the sides.

Indication:

- Lumbar and abdominal support during pregnancy.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 225



BA 300



ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm



ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm

BD 300 U



BD 300 U S



ABDOMINAFIX BLACK UNI.

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm

ABDOMINAFIX ECO

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic polyester fabric that avoids sweating.

Indications:

- Used to support thorax or abdomen area in post-operative period.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 300 E





BA 3002

ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.


Indication:

- Used to support thorax or abdomen area in post-operative period.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm



BA 3003

ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm



BA 3004

ABDOMINAFIX

Functions and Features:

- Ensures thorax or abdomen support, according to the area on which it is used.
- Made of air-permeable, elastic fabric that avoids sweating.

Indication:

- Used to support thorax or abdomen area in post-operative period.



Size

S-M / 1	L-XL / 2
75 cm	99 cm
99 cm	115 cm

BA 302


STERNAFIX (Sternum Corset)

Functions and Features:

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

Indication:

- To support sternum zone in post-operative period.



Size					
S	M	L	XL	XXL	
75 cm	85 cm	95 cm	105 cm	115 cm	
85 cm	95 cm	105 cm	115 cm	125 cm	



BD 302 U


STERNAFIX (Sternum Corset)

Functions and Features:

- Ensures sternum support.
- Made of air-permeable, elastic fabric that avoids sweating.
- The bands on the shoulders avoid the corset to slide downwards.
- Chest pad ensures extra support.

Indication:

- To support sternum zone in post-operative period.



Size		
S-M / 1	L-XL / 2	
75 cm	99 cm	
99 cm	115 cm	



BA 306


STERNUM CORSET WITH STRAP

Functions and Features:

- Ensures immobilization after bypass operation.
- Made of highly air-permeable elastic fabric. Avoids allergy, transfers the moisture to outside.
- Can be used to increase the stabilization by tightening the corset with hand, provided that the patient is to cough after operation.
- Fixed with the shoulder bands to avoid sliding.

Indication:

- After bypass operation.



Size					
S	M	L	XL	XXL	
75 cm	85 cm	95 cm	105 cm	115 cm	
85 cm	95 cm	105 cm	115 cm	125 cm	





BA 303 V

WOOL ELASTIC CORSET WITH VELCRO BAND

Functions and Features:

- Made of wool, elastic fabric.
- Used to keep waist and abdominal area, and kidneys warm, in muscle spasms, and for rheumatic and sciatic pains.
- Can be opened and closed with the velcro on the front side.

Indication:

- Used to keep the affected areas hot in rheumatic and sciatic pains.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

WOOL CORSET

Functions and Features:

- Ensures abdominal and waist support.
- Made of wool, elastic fabric.
- Closed on the front part.
- Keeps warm due to heating feature.

Indication:

- Used to keep the affected areas hot in rheumatic and sciatic pains.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 303



BA 305



COLOSTOMY CORSET

Functions and Features:

- Manufactured of flexible and cotton fabric.
- 26 cm long corset.
- With the aid of the flexible fabric and 2 removeable plastic underwires at each sides of the stoma opening,
- Prevents any unwanted movements of the stoma bag and potential leakages which may occur during mundane and sporting activities.
- Offers perception of comfort and assurance to the patient.
- Suitable for wear underneath daily clothing.

Indication

- Protecting the ostomy collection bag
- Post-ostomy rehabilitation and support purposes
- Decrease risks of parastomal hernia



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

UMBILICAL HERNIA CORSET

Function and Features:

- A support that compresses the hernia area is available.
- Made of air-permeable, elastic, cotton fabric.

Indication:

- Used in umbilical hernias.



Size

S	M	L	XL	XXL
75 cm	83 cm	91 cm	99 cm	107 cm
83 cm	91 cm	99 cm	107 cm	115 cm

BA 290



BA 500

SINGLE TRUSS

Functions and Features:

- Supports the hernia in inguinal area.
- The pad that supports the hernia is made of rubber. Covered with fabric.

Indication:

- Used for hernia, and in post-operative period.



Size Universal - Right / Left

DOUBLE TRUSS

Functions and Features:

- Supports the hernia in inguinal area.
- The pad that supports the hernia is made of rubber. Covered with fabric.

Indication:

- Used for hernia, and in post-operative period.



Size Universal

BA 501



BA 504

SUSPENSORY

Functions and Features:

- Avoids testicles to lie over, and supports the ones that lie over.
- Made of cotton fabric. Has waist belt, and elastic bands that avoid sliding.

Indication:

- Used to support testicles in sports activities.
- Used in vesica varicoses and inflammation
- After varicocele operation.



Size

5	6	7	8	9	10
75 cm	83 cm	91 cm	99 cm	107 cm	115 cm
83 cm	91 cm	99 cm	107 cm	115 cm	123 cm





BA 419 P

FREJKA PILLOW

Functions and Features:

- Made of 12 mm plastozote material.
- Allows compulsory hip and leg movements of the baby, as it is made of elastic material.

Indication:

- Coksit, hip dislocation.



Size

S	M	L	XL	XXL
27 cm 32 cm	32 cm 37 cm	37 cm 42 cm	42 cm 46 cm	46 cm 53 cm


PAVLIC BANDAGE

Functions and Features:

- Positions both hip joints at flexion and abduction.
- Avoids adduction and extension.
- This position and movement ensure spontaneous reduction of the hip, and acetabular development in displaced hip.
- Ensures femur head to stay in acetabulum during free movements of the child.
- Easy-to-wear and remove; made of velcro bands and cotton column.

Indication:

- Hip displacement, subluxation.



Size

S	M	L	XL
24 cm 34 cm	34 cm 43 cm	43 cm 52 cm	52 cm 61 cm

BA 416 P



BA 411 P



HIP-X CRUSER HIP ABDUCTION BRACE

Function and Features:

- Positions the hip joint at flexion and abduction.
- Ensures femur head to stay in acetabulum.
- External part is made of polyethylene while the internal part is made of air-permeable spongy fabric.
- All materials can be washed. Easy-to-apply, light.
- Single orthoses would be enough during treatment, as the adjustment room is wide.

Indication:

- Coxitis, hipdislocation.



Size


S	M	L
32 cm 37 cm	37 cm 42 cm	42 cm 46 cm

BA 409 P

PAVLIK HARNESS TUBINGER TYPE

Functions and Features:

- It has straps that are fastened around the baby's legs and held up by shoulder and chest straps.
- This holds the hips and knees up with the legs apart. This is the best position for the hip joint to be in. It allows contact between the thigh and hip bones and helps strengthen the muscles and ligaments of the hip while it is developing.
- Soft TPR material to make baby feel comfortable.
- Adjustable pearl strap fits most of the babies.



Size			
1	2	3	
24 cm	43 cm	61 cm	80 cm
43 cm	61 cm		



BA 415 U

HIPFIX (Post-Operative Hip Abduction Orthosis)

Functions and Features:

- It is an anatomically designed hip abduction orthosis made of thermoplastic material, which stabilizes the hip joint.
- Internal surface that touches the skin is covered with towel cloth.
- Flexion can be adjusted up to 120 degree, using 15 degree steps
- Abduction angle can be adjusted.

Indication:

- Ensures stabilization in total hip prosthesis or revision.
- Used in femur head resection,
- After femoral or acetabular osteotomy,
- Recurrent hip dislocation cases.



Size Universal -Right / Left





KNEE ORTHOSIS

BA 20103

GENUFIX STABLE (Knitted Kneepad With Hinge)

Function and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect.
- Supports and stabilizes the knee.
- Medial and lateral joint supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.

Indication:

- Collateral ligament injuries. Knee instability.
- ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	43 cm
31 cm	34 cm	37 cm	40 cm	43 cm	



BA 30301

GENUAI RTEC STABLE (Kneepad with Hinge)

Functions and Features:

- Ensures medial and lateral stabilization and support.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.

Indication:

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	43 cm
31 cm	34 cm	37 cm	40 cm	43 cm	



BA 30301 P

GENUAI RTEC STABLE (Kneepad with Hinge) PEDIATRIC

Functions and Features:

- Ensures medial and lateral stabilization and support.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of air-permeable, anti-allergic, skin-friendly airtec (spacer) fabric that avoids sweating.
- Airtec fabric that can be stretched to four directions ensures balanced compression.

Indication:

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



Size		
XS	XXS	XXXS
24 cm	21 cm	19 cm
26 cm	23 cm	21 cm





BD 40301 U

UNIFIX STABIL (Universal Knee Pad with Hinge)

Functions and Features:

- Supports and stabilizes the knee.
- Steel joint in medial and lateral supports the collateral connective tissues, and allows controlled movement.
- Upper and lower bands increase stabilization.
- Made of 3.5 mm spongy Neoprene fabric as one size.

Indication:

- Collateral ligament injuries, knee instability, ACL injuries, osteoarthritis, meniscus degenerations, gonarthros



Size Universal

IMMOFIX (Knee Mobilizer)

Functions and Features:

- Ensures immobilization at 0 degree extension of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of three different panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

Indication:

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



Size

S	M	L
28 cm	31 cm	34 cm
31 cm	34 cm	37 cm



S - 48 cm
M - 56 cm
L - 62 cm

BA 400



BA 400 P

IMMOFIX (Knee Mobilizer) PEDIATRIC

Functions and Features:

- Ensures immobilization at 0 degree extension of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of three different panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

Indication:

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



Size

XS	XXS	XXXS
24 cm	21 cm	19 cm
26 cm	23 cm	21 cm



XS - 40 cm
XXS - 30 cm
XXXS - 20 cm

ECO KNEE IMMOBILIZER

Functions and Features:

- Ensures immobilization at 0 degree extension of the knee.
- Supported with aluminum bars appropriate for the anatomy of the knee on medial-lateral and posterior.
- Formed of one panels, offers wide adjustment opportunity and can be used with any leg
- Easy-to-wear and remove.

Indication:

- Patella dislocation, pre-operative pains, ligament tissue injuries, after PCL and ACL surgeries, for protective purposes after tibia fractures, for knee immobilization when required.



Size

S	M	L
28 cm 31 cm	31 cm 34 cm	34 cm 37 cm



S - 48 cm
M - 56 cm
L - 62 cm

BA 401



BA 524

GENUROM ECO LONG (Rom Ortois with Adjustable Angle)

Functions and Features:

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

Indication:

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move within a certain range. In cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture



Size Universal



BA 524 K

GENUROM ECO SHORT (Rom Ortois with Adjustable Angle)

Functions and Features:

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg. side bars can be adjusted according to the height of the patient.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

Indication:

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move within a certain range, in cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture.



Size Universal





BA 525

GENUROM COMFORT (Rom Knee Ortois with Adjustable Angle)

Functions and Features:

- Ensures knee immobilization in post-operative period.
- Extension and flexion can be adjusted through 10 degree ranges.
- Offers the opportunity to do exercises in required ranges.
- Aluminum bar can be adjusted in accordance with the anatomic structure of the leg. Side bars can be adjusted according to the height of the patient.
- Highly comfortable thanks to 4 foam-backed fabric pads and these pads can be removed and washed.

Indication:

- After ACL, PCL, MCL and LCL operations, to control varus-valgus, in cases where the knee is allowed to move with in a certain range, in cases requiring physiotherapy, to avoid or treat contracture, anterior cruciate ligament tear, and in patella fracture.



Size Universal

GENUFIX PLUS (Ligament and Patella Supported Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect.
- One spiral baleen is available in medial and lateral

Indication:

- Advanced twists on the knee, swelling, post-operative and post-traumatic irritations, patella subluxation Chondromalacia, Patellar tendinitis



Size

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm

BA 20101



BA 30201

GENUAIIRTEC PLUS (Ligament and Patella Supported Knee Pad)

Function and Features:

- Ensures patella and ligament stabilization and support.
- Made of air-permeable, skin-friendly airtex (spacer) fabric that avoids sweating.
- Silicon support reduces pain and edema thanks to its massage effect.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization.

Indication:

- Advanced twists on the knee, swelling,
- Post-operative and post-traumatic irritations,
- Patella subluxation, Chondromalacia, Patellar tendinitis



Size

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm

BA 30201 P

GENUAIRTEC PLUS (Ligament and Patella Supported Knee Pad)

Function and Features:

- Ensures patella and ligament stabilization and support.
- Made of air-permeable, skin-friendly airtex (spacer) fabric that avoids sweating.
- Silicon support reduces pain and edema thanks to its massage effect.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization.

Indication:

- Advanced twists on the knee, swelling,
- Post-operative and post-traumatic irritations,
- Patella subluxation, Chondromalacia, Patellar tendinitis



Size

XS	XXS	XXXS
24 cm	21 cm	19 cm
26 cm	23 cm	21 cm



BA 30203

GENUAIRTEC LIGA (Lateral J Supported Knee Pad)

Functions and Features:

- Fixes the patella from lateral to medial using the "J" shaped silicon support.
- Spiral baleens in medial and lateral increase stability.
- Made of air-permeable, skin-friendly airtex fabric that avoids sweating.

Indication:

- Lateral patellar subluxation, dislocation or patella femoral disorder.



Size (Right / Left)

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BD 40201 U

UNIFIX LIGA (Ligament and Patella Supported Knee Pad)

Function and Features:

- Made of 3.5 mm, air-permeable, spongy Neoprene fabric.
- Avoids patella compression thanks to open patella.
- Two spiral baleens are available in medial and lateral.
- Upper and lower bands increase stabilization, and ensure additional support.

Indication:

- Advanced twists on the knee, swelling, post-operative and post-traumatic irritations, patella subluxation, Chondromalacia, Patellar tendinitis, sports activities.



Size Universal





BA 20102

**GENUFIX
(Patella Silicon Supported Weaved Knee Pad)**

Functions and Features:

- Supports and stabilizes the patella.
- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.
- Silicon support reduces pain and edema thanks to its massage effect

Indication:

- Advanced twists on the knee, swelling, patella subluxation, Chondromalacia, Patellar tendinitis.



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	

GENUAIRTEC (Patella Silicon Supported Knee Pad)

Function and Features:

- Controls patella with silicon support.
- Avoids patella compression thanks to open patella.
- Made of air-permeable, airX® fabric that avoids sweating.
- One spiral baleen is available in medial and lateral.

Indication:

- Patella subluxation, Chondromalacia, Patellar tendinitis



Size					
S	M	L	XL	XXL	
28 cm	31 cm	34 cm	37 cm	40 cm	
31 cm	34 cm	37 cm	40 cm	43 cm	

BA 30101



BD 40101 U

UNIFIX PATELLA (Patella Supported Knee Pad)

Functions and Features:

- Controls patella with silicon support.
- Avoids patella compression thanks to open patella.
- Made of 3.5 mm, air-permeable, spongy Neoprene fabric.
- Standard size. Appropriate for all patients.

Indication:

- Patella subluxation, Chondromalacia, Patellar tendinitis



Size Universal

BA 30104

GENUAI RTEC TENDON

Functions and Features:

- Patella tendon is supported with anatomically appropriate pad.
- Reduces the pain thanks to the compression.
- Ensures extra stabilization with the upper band.
- Made of air-permeable airtec fabric that avoids sweating.

Indication:

- Osgood-schlatter disease, Patellar tendinitis, Chondromalacia



Size

S	M	L	XL	XXL
28 cm	31 cm	34 cm	37 cm	40 cm
31 cm	34 cm	37 cm	40 cm	43 cm



BD 30105 U

GENUAI RTEC BAND UNI.

Functions and Features:

- Patella tendon is supported with an anatomically appropriate silicon pad.
- Reduces the pain thanks to the compression.
- Made of air-permeable airtec fabric that avoids sweating.

Indication:

- Osgood-schlatter disease, Patellar tendinitis, Chondromalacia



Size Universal



BA 940

GENU X PATELLAR BAND

Functions and Features:

- Support is provided on the patella tendon with a silicone pad suitable for the anatomic shape.
- It is aimed to reduce pain with compression.
- Air permeable and moist regulator.

Indication:

- Patellar tendinitis
- Chondromalacia



Size

S-M / 1	L-XL / 2
28 cm	34 cm
31 cm	37 cm





BA 20104-3

GENUFIX BASIC (Elastic Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

Indication:

- To relief the pain on knees, and to reduce edema.



Size

S	M	L	XL	XXL
28 cm 31 cm	31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 43 cm

GENUFIX BASIC (Elastic Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

Indication:

- To relief the pain on knees, and to reduce edema.



Size

S	M	L	XL	XXL
28 cm 31 cm	31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 43 cm

BA 20104-4



BA 20104-5



GENUFIX BASIC (Elastic Knitted Knee Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Reduces pain and edema through compression.
- Decompression is ensured with a special knitting technique in popliteal area.
- Air-permeable and moist regulator.

Indication:

- To relief the pain on knees, and to reduce edema.



Size

S	M	L	XL	XXL
28 cm 31 cm	31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 43 cm

KNEE SLEEVE LUX

Functions and Features:

- It is produced with circular knitting technique that provides balanced compression and can stretch in both directions.
- Thanks to equal compression, it reduces oedema and oedema related pain.
- Air permeable and moist regulator.
- Supports the circulatory system and stimulates muscle regeneration.
- Stimulation of proprioceptors has a positive effect on the sensorimotor system and hence muscle control.

Indication:

- Oedema and pain in and around Patella



Size	S	M	L	XL	XXL
	28 cm	31 cm	34 cm	37 cm	40 cm
	31 cm	34 cm	37 cm	40 cm	43 cm

BA 900



BA 304

GENUWOOL (Wool Knee Pad)

Functions and Features:

- Keeps the knees warm.
- Made of wool, elastic fabric.

Indication:

- Used in rheumatic, sciatic, and knee pains.



Size	S	M	L	XL	XXL
	28 cm	31 cm	34 cm	37 cm	40 cm
	31 cm	34 cm	37 cm	40 cm	43 cm

UNIFIX NEOPRENE CALF SUPPORT

Functions and Features:

- Made of air-permeable, spongy neoprene fabric to support the calf.
- Ensures support with compression feature, and reduces the pain.

Indication:

- In myotonia,
- For protective purposes during and after sportive activities,
- To warm and compress the area.



Size Universal

BD 40402 U



BD 40401 U

UNIFIX NEOPRENE THIGH SUPPORT

Functions and Features:

- Made of air-permeable, spongy neoprene fabric to support the calf.
- Ensures support with compression feature, and reduces the pain.

Indication:

- In myotonia,
- For protective purposes during and after sportive activities,
- To warm and compress the area.



Size Universal



ANKLE ORTHOSIS

BA 20201

MALLEOFIX (Malleolar Supported Elastic Weaved Ankle Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Silicon pads that press around the malleolar reduce the pain through massage effect, and avoid edema.

Indication:

- Acute ankle twists, in case of any pain, Soft tissue traumas, Arthritis cases.



Size

S	M	L	XL	XXL
19 cm 21 cm	21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm



BD 40501 U

MALLEOLAR SUPPORTED ANKLE SUPPORT UNI.

Function and Features:

- Made of non-allergic, air-permeable airtex fabric that avoids sweating.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.

Indication:

- Acute ankle twists, in case of any pain,



Size Universal



BD 589 U

ANKLE SUPPORT LIGAMENT SUPPORTED WITH SPIRAL BALEENS

Functions and Features:

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate spiral baleens
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtex fabric that avoids sweating.

Indication:

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size Universal





BD 40503 U

ANKLE SUPPORT LIGAMENT SUPPORTED WITH PAD UNI.

Functions and Features:

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate rigid plastic pieces.
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtex fabric that avoids sweating.

Indication:

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size Universal

MALLEOAIIRTEC (Ligament Supported ankle Pad)

Functions and Features:

- Ensures ankle support and stabilization.
- Internal and external malleolar avoids twists by being supported with anatomically appropriate rigid plastic pieces.
- Eight bandages on the ankle increase stabilization.
- Reduces edema and pain thanks to the compression.
- Made of air-permeable, skin-friendly airtex fabric that avoids sweating.

Indication:

- Used in stable malleolar fractures, medial and lateral ligament injuries,
- Acute ankle twist, chronic ankle instability,
- Ankle distortion indication.



Size

S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm

BA 30503



BA 20203

ACCHILLFIX (Achilles Tendon Supported Ankle Pad)

Functions and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- The silicones that are designed in accordance with the anatomical structure of Achilles tendon
- Offers support with massage and compression effect.
- The silicon heel pads that can be placed in the shoes reduces the Achilles tension, and relieves the pain.

Indication:

- Used in acute and chronic inflammation of Achilles tendon, post-operative and post-traumatic irritation, Achillodynia, and Achillobursitis.

Size



S	M	L	XL	XXL
19 cm	21 cm	23 cm	25 cm	27 cm
21 cm	23 cm	25 cm	27 cm	29 cm



UNIVERSAL (Neoprene Ankle Pad)

Functions and Features:

- Made of air-permeable, Neoprene fabric.
- Supports the ankle in slight ankle twists, and soft tissue injuries.
- Reduces the pain and edema thanks to the compression.
- (standard size)

Indication:

- Used in acute ankle twist,
- In case of pain, and in arthritis cases.



Size Universal

BD 40502 U



BA 20204-5

MALLEOFIX BASIC (Elastic Ankle Support)

Function and Features:

- Knitted with three-dimensional knitting-weaving technique that ensures balanced compression and stretches to both sides, which is appropriate for the anatomical structure.
- Air-permeable and moist regulator.
- Reduces the pain and removes the edema thanks to the compression feature and massage effect.

Indication:

- Used in acute ankle twists, in case of any pain, soft tissue traumas, and arthritis cases.



Size

S	M	L	XL	XXL
19 cm 21 cm	21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm

BA 20204-3

BA 20204-4

ANKLE SLEEVE LUX

Function and Features:

- It is produced with circular knitting technique providing balanced compression.
- The inner surface is made of cotton outer surface polyamide thread.
- Air permeable and moist regulator.
- It reduces oedema and oedema related pain with compression.
- It supports healthy and correct walking.

Indications:

- Acute ankle sprains
- Soft tissue traumas
- Osteoarthritis, Arthritis



Size

S	M	L	XL	XXL
19 cm 21 cm	21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm

BA 903





BA 538 UNI

MALLEOAIR (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External parts are made of plastics, and shaped anatomically.
- Internal air bladder is blown up as required, and ensures full support.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, and post-operative cases.



Size Universal



BA 589 UNI

FOAMCAST (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External plastics are shaped anatomically. Fits the ankle.
- Latex foam coveren with fabric takes the anatomical shape of the ankle.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, ad post-operative cases.

Size Universal



BA 589 P UNI

FOAMCAST (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization of ankles.
- Easy-to-use as it allows movement in sagittal plane.
- External plastics are shaped anatomically. Fits the ankle.
- Latex foam coveren with fabric takes the anatomical shape of the ankle.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, ad post-operative cases.

Size Universal

BA 588 UNI

FOAMCAST PLUS (Ankle Stabilization Orthosis)

Functions and Features:

- Ensures medial and lateral stabilization in ankles. Easy-to-use as it allows movement in sagittal plane. External plastics are shaped anatomically.
- Fits the ankle. Latex foam covered with fabric takes the anatomical shape of the ankle.
- Ensures extra stabilization with the extra bandage that covers the ankle to metatarsus.

Indication:

- Used in chronic ankle twists, treatments of stable malleolar fractures, medial and lateral ligament injuries, ankle distortion, and post-operative cases.

Size Universal



BA 519

ECO ROM WALKER

Functions and Features:

Stabilizes the ankle with angle adjustable hinge, and steel bars.

- Dorsiflexion and plantar flexion can be adjusted through 10 degree ranges.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.

Indication:

- Used in Achilles tendon repairs, acute twist of ankles, Serious soft tissue injuries, Stable fractures of angles, and non-healing fractures in aistal tibia and fibula.

Size Universal



S	M	L
36	39	42
39	42	45



BA 522

AIR ACCILLIES BOOT LONG

Functions and Features:

- Ankle plantar flexum height can be adjusted using the pads on the heel.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.
- Air Inflatable

Indication:

- Used in Achilles tendon repairs, Acute twist of ankles, Serious soft tissue injuries, Stable fractures in ankles, And non-healing fractures in aistal tibia and fibula.



Size	S	M	L
	36	39	42
	39	42	45





BA 522 K

AIR ACCHILLES BOOT SHORT

Functions and Features:

- Ankle plantar flexion height can be adjusted using the pads on the heel.
- The product is light and easy-to-use.
- Pads can be washed.
- The plantar is rocker, and does not slide.
- Offers the opportunity to do exercises as it can be removable.
- Air inflatable

Indication:

- Used in Achilles tendon repairs, Acute twist of ankles, Serious soft tissue injuries, Stable fractures in ankles, And non-healing fractures in aistal tibia and fibula.



Size

S	M	L
36	39	42
39	42	45

FOOT GUARD (Foot Position Splint)

Function and Features:

- Made of polypropylene material.
- Keeps the ankle at 90°.
- Used in hemiplegia, paraplegia, peroneal nerve paralysis, and ankle instability.
- Avoids the foot to experience plantar flexion in patients in intensive-care unit.
- Furthermore, the portable T band under the foot avoids it to tend to right or left.
- The T band can be removed when required.
- Avoids eschar injuries.



Size

S	M	L
36	39	42
39	42	45

BA 420



BA 536

AFO

Function and Features:

- An orthosis made of polypropylene that can be worn with shoes and is used to keep the ankle angle at 90 degrees.
- Used in mild foot deformities caused by nerve injuries such as hemiplegia and disc hernia, and by peroneal nerve disorders.
- Has three different sizes.
- Right and left options are available.
- State the size of your calf.



Size (Right - Left)

S	M	L
36	39	42
39	42	45

BA 583

ROCKER PLASTER SLIPPER

Functions and Features:

- Makes it easy to walk on your plaster covered feet thanks to its special base (Rocker).
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.



Size			
S	M	L	
36	39	42	
39	42	45	



BA 583-2

PLASTER SLIPPER

Functions and Features:

- Makes it easy to walk on your plaster covered feet thanks to its special base (Rocker).
- Avoids the plaster to be broken, base to contaminate and slide.
- Easy-to-use with the adjustable velcro fastener system.



Size			
S	M	L	
36	39	42	
39	42	45	



BA 108

HALLUX VALGUS SHOES

Functions and Features:

- Post-operative shoe designed in accordance with the anatomical form of the foot for the purpose of supporting mobilisation after Hallux Valgus surgery.
- It is manufactured from skin friendly, antibacterial materials. Its nonskid, platform sole speeds up the healing process preventing pressure on the operated region.
- The foot contact surface is soft. It does not put stress on metatarsus and phalanges, normalizes the extra stress on the heel as well. The upper surface is totally made of breathable fabric.
- Specially designed hook and loop fastener allows it to be used on both feet. It reduces possible post-operative edema on foot.



Size			
S	M	L	
36	39	42	
39	42	45	





BA 701

SILICONE INSOLE

Function and Features:

- Used for painful plantars and heels, hyperkeratosis, plantar fasciitis, and as pain absorber.
- Avoids the pains caused by certain disorders such as pescausus.
- Avoids partial pressure and tension on the feet due to remaining standing position for a long time, and posture disorders.
- As the density is adjusted according to the foot anatomy (blue parts are less dense), the product protects sensitive and boned areas by increasing the adaptation to plantar.
- Absorbs the impacts on knee and ankle joints in dynamic activities such as jogging, running, and jumping.
- Avoids the pain likely to occur on the same joint due to imbalanced pressure.
- Can be used during sports, daily life, and all activities.
- Made of 100% silicone, and is hypoallergenic.
- Physical life can be prolonged by washing the product with cold water every other day, and drying.



Size

XS	S	M	L	XL	XXL
32	35	38	41	44	47
34	37	40	43	46	49

SILICONE INSOLE WITH ARCH SUPPORTED

Function and Features:

- Used for painful plantars and heels, hyperkeratosis, plantar fasciitis, and as pain absorber.
- Avoids the pains caused by certain disorders such as pescausus.
- Avoids partial pressure and tension on the feet due to remaining standing position for a long time, and posture disorders.
- As the density is adjusted according to the foot anatomy (blue parts are less dense), the product protects sensitive and boned areas by increasing the adaptation to plantar.
- Absorbs the impacts on knee and ankle joints in dynamic activities such as jogging, running, and jumping.
- Avoids the pain likely to occur on the same joint due to imbalanced pressure.
- Can be used during sports, daily life, and all activities.
- Made of 100% silicone, and is hypoallergenic.
- Physical life can be prolonged by washing the product with cold water every other day, and drying.



Size

XS	S	M	L	XL	XXL
32	35	38	41	44	47
34	37	40	43	46	49

BA 701-1



BA 546

FABRIC COVERED SILICONE SUPPORT

Function and Features:

- In addition to the features of the product coded BA 701,
- Ensures comfort as it is fabric covered, avoids odor, and reduces sweating.

Indication:

- Flat feet
- Plantar fasciitis
- Pain on the heels
- Pain on metatarsus



Size

S	M	L	XL	XXL
35	38	41	44	47
37	40	43	46	49

SILICONE HEEL PAD

Function and Features:

- Used for heel support.
- Raises the heel in achilles tendon operations, and reduces the stress on achilles ligament.
- Can be used in jogging as well.
- Can be easily used in training and daily shoes.
- The product is portable, and should be used in all shoes.
- It is appropriate for the patients receiving heel treatment to use the products in the house.
- Made of 100% silicone. The product is orthopedic and hypoallergenic.



Size

S	M	L
35	38	41
37	40	43

BA 682



BA 683

SILICONE INTERDIGITAL APPARATUS

Functions and Features:

- Preserves the health gained (treatment process) using Hallux Valgus Night Splint in posture disorders of the toe.
- Avoids the irritation caused by overlapping toes.
- Made of silicon material.
- Anti-allergenic.
- Easily washed with water (cold water is recommended).
- Furthermore, ease of use is ensured if powdered with talcum powder.



Size

S	M	L
35 37	38 40	41 43



BA 537

SILICONE METATARSUS PILLOW

Functions and Features:

- Made of 10% gel silicone.
- Maximum elastic.
- The product is orthopedic and hypoallergenic.
- Relieves the pain in metatarsus areas, and avoid rash.
- Passed through the second toe and placed on the area where pain and rash is felt.
- Avoids callus on the same area.



Size Universal



BA 105

TRANSVERS ARCH SUPPLEMENT

Functions and Features:

- Used to support longitudinal arch in patients with pes planus or fallen arches



Size (1 - 2)



BA 106

METATARTUS PAD

Functions and Features:

- Used to support metatarsal arch in patients with pes planus or fallen arches



Size (1 - 2)



BA 685



BUNION SHIELD

Function and Features:

- Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.


Size Universal

BA 541

EPIN HEEL SUPPORT

Functions and Features:

- Serve as a pad on the heel, and reduces the stress on the Achilles tendon;
- Reduces the irritation and pain (Epin Calcanien) on the heel.
- Made of polyurethane, and has three different sizes.
- Bottom part is 100% rubber, not sponge.



BA 703

CORK INSOLE WITH EPIN CALCANIUS

Function and Features:

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



Size

33 - 34 | 35 - 36 | 37 - 38 | 39 - 40 | 41 - 42 | 43 - 44



BA 702

CORK SUPPORT

Function and Features:

- Longitudinal arch support and metatarsus support are made of cork.
- The outer surface is covered with artificial leather.



Size

33 - 34 | 35 - 36 | 37 - 38 | 39 - 40 | 40 - 41 | 42 - 43 | 44 - 45



BA 800

INSOLE WITH STEEL SUPPORT

Functions and Features:

- Used in the treatments of Pes Planus, Pes Cavus, Metatarsalgia, and Epin Calcainen deformities.
- Supported with stainless steel, and has apparent supports.
- Preserves elasticity and dynamism for a long time.
- The outer surface is covered with artificial leather.



Size

33 - 34 | 35 - 36 | 37 - 38 | 39 - 40 | 41 - 42 | 43 - 44



BA 544 U

ELASTIC HALLUX VALGUS NIGHT SPLINT

Functions and Features:

- Avoids Hallux Valgus deformity to progress with the pull applied on the toe.
- Can be used in shoes.

Indication:

- Used for Bunion, hallux valgus, din post-operative period.



Size Universal (Right – Left)



BA 542 UNI

HALLUX VALGUS NIGHT SPLINT

Functions and Features:

- Internal surface is plastazote, while the outer surface is rigid plastic.
- Avoids Hallux Valgus deformity to progress with three point principle.
- Used to bend toe outwards in post-operative process in the Halluks Valgus problems.
- Used in bed, not with shoes.



Size Universal



BD 542 U

HALLUX VALGUS NIGHT SPLINT UNI.

Functions and Features:

- Internal surface is plastazote, while the outer surface is rigid plastic.
- Avoids Hallux Valgus deformity to progress with three point principle.
- Used to bend toe outwards in post-operative process in the Halluks Valgus problems.
- Used in bed, not with shoes.



Size Universal





The advertisement features a light blue background with a subtle pattern of bamboo leaves and several pillows. A large, light-colored, rectangular pillow is the central focus. A line points from the text 'Memory foam Visco' to the bottom left corner of this pillow. Another line points from the text 'Bamboo Pillowcase' to the top right corner of the same pillow. At the bottom of the image, there is a cross-section diagram of a human head and neck resting on a pillow, showing the internal structure of the pillow and the spine. Below this diagram are three icons: a hand icon labeled 'MEMORY', a wavy line icon labeled 'ERGONOMIC', and a leaf icon labeled 'ECOLOGICAL'.

Biofix[®]

ORTHOPEDIC VISCO PILLOWS

Bamboo
Pillowcase

Memory foam
Visco



MEMORY



ERGONOMIC



ECOLOGICAL

Product Info:

- Pillow case is available. Case can be removed and washed

**BA 609
VISCO NECK PILLOW**



Dimensions: 59 x 40 x 12-10 cm

**BA 610
VISCO CONTOUAR NECK PILLOW**



Dimensions: 58 x 39 x 11-9 cm

**BA 615
VISCO ANATOMIC NECK PILLOW**



Dimensions: 53 x 33 x 13 cm

**BA 623
VISCO NECK PILLOW SMALL**



Dimensions: 54 x 35 x 10-8 cm

**BA 614
VISCO TRAVEL NECK PILLOW**



Dimensions: 31 x 25 cm

**BA 622
VISCO VEHICLE NECK PILLOW**



**BA 611
VISCO WAIST PILLOW**



Dimensions: 37 x 34 cm

**BA 616
VISCO ANATOMIC WAIST PILLOW**



Dimensions: 53 x 42 cm

BA 612
VISCO SEATING CUSHION



Dimensions: 44 x 44 cm

BA 613
VISCO U FORM SEATING CUSHION



Dimensions: 44 x 40 cm

BA 619
VISCO COCCYX SEATING CUSHION



Dimensions: 43 x 37 x 9 cm

BA 620
VISCO SEATING CUSHION



Dimensions: 45 x 42 x 8 cm

BA 612 P
POLYURETHANE SEATING CUSHION



Dimensions: 44 x 44 cm

BA 617
VISCO CYLINDER PILLOW



Dimensions: 43 x 17 cm

BA 618
VISCO CYLINDER PILLOW HALF



Dimensions: 43x18 cm

BA 621
VISCO SUPPORT PILLOW



BA 548

MOBILE WALKER

Product Info:

- An appropriate walker for elderly, patients having troubles with walking, patients experienced foot paralysis, and it can be used to rehabilitate the patients receiving physiotherapy treatment.
- It can be folded thanks to the mobile joints of the angled pipes, and ensures easy and safe walking by lifting two legs and pushing forward.
- Easy-to-use.
- Offers an easy experience thanks to the high quality, soft sponge knobs appropriate for the anatomical structure of the hand.
- Height can be adjusted according to the demands of the patients, thanks to the click-fit legs.
- Does not slide on wet surface thanks to the rubber bottom supports.
- Does not scratch or damage the gentle surfaces such as parquet.
- Metal parts can be wiped with a moist duster.
- Ensures easy and safe walking by lifting two legs and pushing forward.



BA 549

ARM PIT CRUTCH

Product Info:

- Arm pit support is available. No burr likely to harm your arm pit is found on the plastic support.
- Handle support is available.
- Handle height can be adjusted.
- Height can be adjusted.



BA 550

FOREARM CRUTCH

Product Info:

- Arm support is made of rigid plastic.
- Hence, no burr likely to harm your arm is found on the plastic support.
- Made of light aluminum alloyed profile.
- Height can be adjusted.





BA 551

TRIPOD WALKING STICK

Product Info:

- Balance walking stick, tripod aluminum.
- Made of aluminum.
- Height can be adjusted.
- Handle is appropriate for the shape of the hand.
- Has three legs.
- Can be used in both sides (right and left).
- Ensures stable walking support.

SINGLE BALANCE STICK

Product Info:

- Single Aluminum Walking Stick
- Made of aluminum.
- Height can be adjusted.
- A handle appropriate for the shape of the hand is available.

BA 552



BA 553

DENNIS BROWN ORTHESIS

Functions and Features:

- Controls dorsiflexion and abduction angle of the foot by assembling special PEV boot on an aluminum bar.
- Abduction angle is stable, while the external rotation is adjusted.
- Ankle dorsiflexion angle should be adjusted by the technician according to the instructions of the doctor.

Indication:

- Pes ekinovarus (PEV), Hip Dysplasia

Size:

- The product is combined according to shoe size and shoulder width of the patient.

To order the product separately: BA 553 - Dennis Brown Splint *BA 553-1 Dennis Brown Boot.



BA 852

TOILET SEAT BOOSTER FOR CLOSET

Product Info:

- BA-853 Toilet Seat Booster (With Legs)
- Toiler Seat Booster with Cover
- Toiler seat booster is first class, white, and plastic.
- Used for the patients having trouble in standing up and sitting down after hip, knee, and waist surgery.
- Can be used in all toilets.
- Boosts the sitting height of the toilet for 11 cm.
- The toilet seat booster has key apparatus, anti-bacterial, and does not hold any bacteria.
- The product has a cover, and this cover can open and close the seat.



Privately produced for you...



Our company has been producing for orthopaedic market since 1999.

We continue to insert for secure product quality and variety.

We use appropriate material to produce healthy products.

Starting with supply of gain, we are completing all phases of production line, such as narrow fabric weaving. Flat knitting and sewing.

So maintaining full products quality, we are capable of producing semi-finished and finished good which are privately required by our costumers.



We have advantages of ;

- Height quality,
- Rich variety of products,
- Being closed to the market,
- Quick services capability,
- Reasonable prices.

NOTES

A large, light gray rectangular area with a torn paper edge effect, containing horizontal dashed lines for writing notes. The edges are irregular and jagged, giving it the appearance of a piece of paper that has been torn from a notebook. The interior of the rectangle is filled with horizontal dashed lines, providing a guide for writing.



bioaktif



BİOAKTİF Tıbbi Ürünler ve Teknik Tekstil Sanayi Ticaret Limited Şirketi
Haraççı Mah. Bayrak Cad. No: 26 34381 Arnavutköy - İstanbul - Turkey
Phone: +90 212 567 92 34 / +90 533 215 41 84
Fax: +90 212 674 88 92 E-mail: info@bioaktif.com
www.bioaktif.com